
dr. aktive **XL** CRYO PRO™

- EN** User manual
- FR** Manuel d'utilisation
- DE** Bedienungsanleitung
- IT** Istruzioni per l'uso
- ES** Manual de empleo
- NL** Gebruikershandleiding
- PT** Manual do usuário

IFU - Cryo Pro XL
4648001820 - 02/2025

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Thank you for purchasing our product. To ensure proper operation and use, please read this manual carefully prior to use.

Please keep this manual properly for future reference.

Contents

Composition_	/P04
Warning and cautions_	/P05
Intended use_	/P06
Contraindications_	/P06
Operating the system_	/P07 - 09
Cleaning and Disinfection _	/P09
Warranty_	/P10
Disposal and recycling_	/P10
Cuffs _	/P11 - 12
Accessories_	/P12 - 13
Troubleshooting_	/P13 - 14
Specifications_	/P14
Symbols used_	/P15 - 16

Composition

① Main Unit	
② Connector Hose DA03020024	
③ AC Adaptor	
④ Drain Gadget DA03020009	

Warning and Cautions

Important Information

- Please read the instruction manual carefully before use. Use according to the instruction manual and shall not be used for other purposes.
- The applied part is the cuff of different body parts, when using, a medical isolation pad should be placed between the cuff and the skin.
- Use only the cuff and power adapter supplied by the manufacturer.
- This product is under certain pressure when used. It is forbidden to use on any wound.
- This product is portable, take care when moving. Do not drop or press this product. Do not disassemble or destroy this product. Do not change the moving parts of product. Do not use if the electrical parts are damaged/leaked.
- If the button symbols are not clear , please find the manufacturer.
- To avoid damage, do not operate the system if water is not filled in the cooler.
- The product is not refrigeration equipment, need to add pure ice water or ice cubes to adjust the temperature. Do not add hot water, hot compress is prohibited.
- To avoid risk of electrical shock, always turn off the system and disconnect the power line cord from its electrical outlet when not in use or before adding or emptying ice and water.
- Follow the treatment recommendations of your health care practitioner for the use of this product; improper or prolonged use of the product could result in tissue damage. For patients after surgery, it is not recommended to use 4 - 9 mode.
- Monitor the skin condition of affected area for any burning, itching, increased swelling, or pain. If any of these happens, or any changes in skin appearance occur such as blisters, increased redness or discoloration, users should stop using and consult a physician.
- Children, the elderly and those who have cognitive disabilities or communication barriers cannot use the product alone without super vision.
- When using this product, the cooler and the cuff should be placed at the same level, and the height difference between the two should not be more than 40CM.
- The main unit should not be tilted more than 10 degrees when in use. Do not use this product when water leaks.
- If device malfunctions or you feel uncomfortable, take off the cuff immediately and disconnect from power source.
- The equipment operates with noise, electromagnetic radiation, and electromagnetic interference to the environment.

Indented use

It is intended to be used for local treatment in cases where cold treatment is required and to temporarily reduce swelling and pain after surgery or injury. It is intended to be used by or on the order of licensed healthcare professionals in hospitals, outpatient clinics, athletic training or home settings.

Contradictions





Important Information

DO NOT USE THE PRODUCT BEFORE CONSULTING A DOCTOR IF ANY OF THE FOLLOWING APPLIES:

- Patients with Raynaud disease.
- Acute, unstable (untreated) fracture.
- Who have rheumatoid arthritis in the affected region.
- Who have open wounds, or localized unstable skin condition.
- Who are intolerant to cold compress or have cold hypersensitivity.
- The pregnant, patients with hyper tension, extreme low blood pressure or circulation disorders.
- Who have vascular diseases such as prior frostbite, diabetes, arteriosclerosis or ischemia.
- Who have hematological dyscrasias which affect thrombosis such as paroxysmal cold hemoglobinuria, cryoglobulinemia, sickle-cell disease or serum cold agglutinins.
- Who are in the acute stages of inflammatory phlebitis in the affected region.


Operating the system

1. Press the lid release button to open the cooler lid (1), add water to fill-line indicated on label within the cooler (2). ENSURE WATER IS ABOVE THE FILL-LINE. Add as much ice as you want but do not exceed the ice max. mark (3). After ice is added close the lid 4.


2. Apply the selected cuff (knee cuff as example). DO NOT APPLY DIRECTLY TO SKIN.


Add an isolation layer between the wrap and the skin (pyjamas, jogging bottoms etc.).

3. Connect cuff to cooler using connector hose. For best results, place the cooler at the same level as the cuff. The height difference between the two should not be more than 40cm.


4. Connect the AC adaptor to the cooler (12), then plug the AC adaptor (13).

5. Press the ON/OFF button for about 3-5 seconds, the screen should light up and show Welcome.

6. Press the Plus or Minus button to choose the mode.

Mode	Water circulationcycles	Time on (min)	Time off (min)
1	Low "ON:8s/10s/12s; OFF:10s"	20	20
2		40	40
3		60	60
4	Mid "ON:14s/16s/18s; OFF:12s"	20	20
5		40	40
6		60	60
7	High "ON:20s/22s/26s; OFF:14s"	20	20
8		40	40
9		60	60

Mode: There are 9 modes. Shows which program you have selected.

Press: There are 3 levels: low, medium, high. Represents the time of cycles of water circulation in the cuff during operation.

Time on: Time of active cycles of water circulation.

Time off: Pause time (no cycles of water circulation).

7. Press the Start/Pause button to start or pause the cooler.
When cooling you can press the button Start/Pause to pause.
8. Press the On/Off button for about 3-5 seconds to turn off the cooler.
9. After use, to drain the cuff, disconnect the cuff from the cooler. Connect to the cuff the drain gadget. Once the drain gadget is connected squeeze the cuff to push the water out through the drain gadget (14).
To drain the cooler, remove the lid and pour out the ice and water. Alternatively, you can use the drain gadget (15).

Cleaning and Disinfection

- The following surfaces must be cleaned and disinfected using wipe disinfection: Surfaces of the cooling device, surfaces of the tubes and surfaces of the cuff. Please use a foam disinfectant (alcohol-free or a solution with <5% alcohol).
- Optimal hygiene is ensured if the cuff is used exclusively for one patient at a time.
- Store the product in a cool, dry, well-ventilated, and clean environment, away from corrosive substances.
- To clean the inside of the tubes, fill the reservoir with warm water, Lysol, 70% isopropyl alcohol or white vinegar, then activate the Cryo Pro™ XL. Allow a few pump cycles to flush the system, then drain it thoroughly. ***Do not use soap or detergent, as this may create suds and damage the device and/or cuffs.***

Warranty

The Kinetec SAS guarantee is strictly limited to the free replacement or factory repair of the part or parts that are recognised to be defective.

Kinetec SAS guarantees Dr Aktive® Cryo Pro™ XL for 2 years against any manufacturing defect, from the date of purchase by the consumer.

Kinetec SAS is solely accredited to judge the application of the guarantee for its devices.

The guarantee cannot be applied if the device has been used abnormally or was used in conditions of use other than those contained in our user manual.

The guarantee shall also not apply in the event of deterioration or accident resulting from negligence, compliance or maintenance fault resulting from a transformation of the equipment or an attempt to repair equipment.

Disposable and Recycling

- **Packaging:** The packaging must be separated from plastic and paper/cardboard components and placed in specific recycling locations.
- **Cooler:** It contains electronic components, cables, plastic, steel and aluminium parts. When the device is no longer functional, dismantle and separate into groups of materials and place them in approved recycling units or return the device to Kinetec o MedFac® for destruction. Or contact local authorities to establish the appropriate method of disposing of parts and accessories that are potentially dangerous to the environment.
- **Cuff and hose:** Put it in specific recycling locations.

Cuffs

- 360 Knee Cuff

DA0302X201

- Shoulder Cuff

DA0302X202

- Ankle Cuff

DA0302X203

• **Universal Cuff**

DA0302X204

• **Back & Hip Cuff**

DA0302X207




• **Full Leg**

• **Full Long Leg Cuff**


DA0302X205

DA0302X206

Accessories

Portable Bag DA0302X301	
Power Bank 2.0 DA03020031	

Cuffs • 360 knee cuff - DA0302X201 • Shoulder cuff - DA0302X202 • Ankle cuff - DA0302X203 • Universal cuff - DA0302X204 • Back & Hip Cuff - DA0302X207 • Full Leg Cuff - DA0302X205 • Long Leg Cuff - DA0302X206	
Split Hose DA0302X303	

Troubleshooting

Missing accessory	Contact MedFac to ship the missing item.
Damaged item(s)	Replace the damaged structural component or return the product for repair.
Damage or Stains to product	1.Clean the bag, use damp towel to wipe device cover. 2. Replace the stained part by contacting MedFac
Flashing screen or white/blank screen	1.Disconnect from power supply, restart the machine; or try another adaptor 2.Contact MedFac for repair or exchange.
Does not turn it on	1. Check if the power supply of socket or outlet is normal. 2. Check if the adaptor is damaged (use any port to match the adaptor of 12V-2A DC output to test). 3. Check if there is any damage on display or button area. Replace broken parts accordingly. 4.Contact MedFac for repair or exchange.

No water pressure,	1. Machine is on, there is no water pressure nor circulation, but water pump is working and roaring. Then air exists inside pump, water can not be pumped due to atmospheric pressure. Shake the cooler, or press hose joint of cooler several times after hose disconnected from it , it will thus back to normal. 2.Machine is working, no water pressure or the pressure is too low, and circulation speed is slow; then check if the water inlet and outlet of water pump is blocked, get it cleaned. Check if cuff is overly folded to cause waterway clogged. 3.Machine is on, no water pressure nor circulation,and no sound is generated from equipment; then water pump is damaged and needs replacement. Contact for repair or exchange.
Water leakage	1. Check if hose is properly connected with cooler and cuff; replace the leaking part. 2. Contact for repair or exchange.
Buzzer alarm	1.Check if water line is too low and add more water accordingly; 2.Check if the water inlet is blocked for normal water pumping, clear it. 3.Inspection and replacement of water pump assembly and circuit board may be needed; Contact for customer service.

Specifications

Main unit size: 375x235x267mm

Main unit weight: around 3600g

Water capacity: 6.0L

DC Output: 12V/2A

AC power: 100-240V 50/60Hz

Pressure level: Low/Medium/High, adjustable

Time settings: 20/40/60 minutes, adjustable

Noise: below 55 decibels

Environment Requirements

Operating:

Temperature: 5°C - 40°C

Relative humidity: 10% - 90%

Atmospheric pressure: 700hPa - 1060hPa

Transport & Storage:

Temperature: -25°C - 55°C

Relative humidity: 10% - 90%

Atmospheric pressure: 700 - 1060hPa

Kinetec SAS


Zone Industrielle de Tournes

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08090 - Tournes

France

Symbols used

	Refer to instruction manual		Class II device
	Refer to instruction manual		Medical number
	Serial number		Date of manufacture
	Manufacturer		Medical Device
	BF TYPE device (protection from electric shocks)		For indoor use only
	Manufacture batch number		Keep Dry
	Direct current		Keep away from sunlight
IP22	Prevent solid foreign objects with a diameter greater than 12.5mm from entering; still prevent water droplets from entering when inclined at 15 degrees.		
	Disposal of product according to local regulations.		

	Caution Symbol Indicates that important cautionary information, such as warnings and precautions, cannot be presented on the medical device itself and are presented in the instructions for use.
	Temperature limits to which the device can be safely exposed.
	Humidity limitations that the device can be safely exposed.
	Unique Device Identifier, Indicates a carrier that contains Unique Device Identifier information.

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Composition

① Unité de contrôle	 A black, rectangular control unit with a handle on top and a small display or sensor on the front.
② Tuyau de connexion DA03020024	 A black, flexible connection hose with a metal fitting on one end and a multi-pin connector on the other.
③ Transformateur	 A black, rectangular transformer with a power cord and a multi-pin connector.
④ Accessoire de vidange DA03020009	 A grey, rectangular drainage accessory with a handle and a multi-pin connector.

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