
17SF300=L-22, 17SF300=R-22

DE	Gebrauchsanweisung	7
EN	Instructions for use	15
FR	Instructions d'utilisation	23
IT	Istruzioni per l'uso	31
ES	Instrucciones de uso	39
PT	Instruções de utilização	48
NL	Gebruiksaanwijzing	56
SV	Bruksanvisning	64
DA	Brugsanvisning	72
NO	Bruksanvisning	80
FI	Käyttöohje	87
PL	Instrukcja użytkowania	95
HU	Használati útmutató	103
CS	Návod k použití	111
RO	Instrucţiuni de utilizare	119
HR	Upute za uporabu	128
SK	Návod na použitie	135
BG	Инструкция за употреба	144
TR	Kullanım kılavuzu	152
EL	Οδηγίες χρήσης	160
RU	Руководство по применению	169
JA	取扱説明書	178
ZH	使用说明书	186
KO	사용 설명서	192


1

17SF300=L-22


17SF300=R-22

2


17AD316=22


3







4


5



12


13


14

15

16

17

INFORMATION

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INFORMATION

Date of last update: 2026-09-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide information on the processing of the 17SF300=L-22 and 17SF300=R-22 foot stirrups.

2 Safety-related information

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION**Excessive strain due to use on more than one user**

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on one user only.
- ▶ Observe the maintenance recommendations.

NOTICE**Thermal overloading of the product**

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE**Thermal overloading of product components**

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE**Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand).
- ▶ Do not expose the product to operating temperatures below **-10 °C (14 °F)** or above **+45 °C (113 °F)** (e.g. excessive sunlight, drying on a radiator).

3 Product description

Scope of delivery for 17SF300=L-22, 17SF300=R-22 (see fig. 1)

Item	Qty.	Designation
1	1	Bearing bushing
2	2	Shim plate
3	4	Spiral pin for shim plate (2 of which as service part)
4	2	Fixing screw for shim plate
5	1	Base body of the foot stirrup

Service set

Scope of delivery for 17AD316=22 (see fig. 2)		
Item	Qty.	Designation
1	2	Shim plate
2	4	Spiral pin for shim plate (2 of which as service part)
3	2	Fixing screw for shim plate
4	1	Bearing bushing

3.1 Available sizes

The product size is selected depending on the planned area of application (unilateral or bilateral, medial or lateral) and the user's body weight.

- 17SF300=L-22
- 17SF300=R-22

3.2 Maximum body weight

The product is approved for up to **125 kg**.

4 Intended use

4.1 Indications for use

The product is intended **exclusively** for orthotic treatment of the lower limbs of a user.

The product is suitable for unilateral or bilateral use.

Only use the product for the lamination or prepreg technique.

4.2 Target group

- Qualified personnel
- User

4.2.1 Qualification of qualified personnel

Users may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4.3 Indications

- Partial or total paralysis of the leg muscles
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

4.4 Product combinations

The 17SF300=L-22 and 17SF300=R-22 foot stirrups are intended **exclusively** for use with the orthotic ankle joint 17AD300=22-T.

4.5 Lifetime

The product is designed for a maximum lifetime of **6 years** when used as intended and assembled professionally.

5 Accessories

Lamination dummy set

Scope of delivery for 17AD313=22 (see fig. 3)		
Item	Qty.	Designation
1	1	Lamination dummy – lower section
3	1	Slotted oval head screw
4	1	Insertion aid for joint nut

Not included in 17AD313=22 scope of delivery (see fig. 3)		
Item	Qty.	Designation
2	1	Lamination dummy upper section *)
*) The upper section of the lamination dummy is not part of the lamination dummy set, but is in the scope of delivery of the 17AD300=22-T orthotic ankle joint.		

Assembly tool

Scope of delivery for 17AD312=22 (see fig. 4)			
Item	Qty.	Designation	Reference number
1	1	Assembly/disassembly aid for spiral pins	–
2	1	Assembly/disassembly aid for bearing bushing	30Y671=22

6 Preparing the product for use

The section is intended for qualified personnel.

CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- Observe the alignment, assembly and adjustment instructions.

The foot stirrup is disassembled on delivery.

6.1 Bending the foot stirrup

Required materials and tools: 17AD313=22 lamination dummy set, bending irons

> **Prerequisite:** The orthosis model has been prepared.

- 1) Position the foot stirrup **without a bearing bushing** in the upper section of the lamination dummy.
- 2) Set the upper section of the lamination dummy onto the lower section of the lamination dummy and secure it with the slotted oval head screw.
INFORMATION: Due to the design, the foot stirrup is brought to the neutral position in the upper section. The neutral alignment of the foot stirrup permits the maximum dorsal or plantar range of motion in the orthotic joint.
- 3) **Optionally**, use the 17AD310=22 shaping dummy of the 17AD300=22-T orthotic joint. For details, please refer to the instructions for use 647G1875=*. **INFORMATION: The shaping dummy guarantees the optimal distance of the modules to the orthosis model, taking into account the modules of the 17AD300=22-T orthotic joint that are used individually for the user.**
- 4) **Optionally**, use the bending irons to bend the foot stirrup; remove the foot stirrup from the lamination dummy for processing.
- 5) **NOTICE! Do not bend the insertion zone of the foot stirrup. There is a risk of breakage.**
- 6) **INFORMATION: Do not use the 17AD312=22 assembly tool to bend the foot stirrup.**
- 7) Respect the allowable bending zone (see fig. 5).
- 8) Grind to remove grooves and burrs.
- 9) Screw the lamination dummy onto the alignment insert.
- 10) Fabricate the orthosis.

6.2 Assembling the single components of the foot stirrup

Required materials and tools: 17AD312=22 assembly tool

The assembly/disassembly aid for the spiral pins consists of 2 identical parts and is clamped in the vice. Small magnets are integrated that ensure that the parts of the assembly/disassembly aid are held together and stay in position.

- The short cylinder pins are used to press the spiral pins into the foot stirrup, slide them in and position them precisely. They are marked with **IN** on the assembly/disassembly aid. This mark indicates that the short pins are used to insert the spiral pins.
- The long cylinder pins allow the spiral pins to be pushed out and removed during maintenance. They are marked with **OUT** on the assembly/disassembly aid.

This allows the assembly/disassembly aid to be used both for inserting and removing the spiral pins.

> **Prerequisite:** The foot stirrup is laminated into the orthosis.

- 1) Attach the shim plates to the foot stirrup (see fig. 6).
- 2) Use the fixing screw of the shim plate and the **T8** Torx bit to attach the shim plate to the foot stirrup so that the screw head of the fixing screw contacts the shim plate but the fixing screw is not fully tightened (see fig. 7).
- 3) **NOTICE! The shim plates should still be easy to move. Over-tightening them can cause damage to the roll pins and the assembly aid (see fig. 8).**
- 4) Repeat this step on the opposite side.
- 5) To install the spiral pins, use the assembly disassembly aid for spiral pins from the 17AD312=22 assembly tool set (see fig. 4, item 1).
- 6) **INFORMATION: The illustrations show the insertion of the spiral pins and the use of the short cylinder pins. For a better overview, the foot component has been omitted in some illustrations.**
- 7) Insert the foot stirrup into the contour of the assembly/disassembly aid for spiral pins (see fig. 9).
- 8) Put on the counter piece (see fig. 10).
- 9) Remove the short cylinder pins (see fig. 11).
- 10) Insert the spiral pins into the counter piece from the outside (see fig. 12).
- 11) Then insert the short cylinder pins into the assembly/disassembly aid from the outside (see fig. 13).
- 12) Clamp the assembly/disassembly aid in the vice (see fig. 14).
- 13) Turn the crank of the vice and press in the spiral pins to the stop.
- 14) Remove the assembly/disassembly aid with the foot stirrup and spiral pins from the vice.

- 15) Before handing the orthosis over to the user, clean the fixing screws of the shim plates using a degreasing cleaner, secure them with Loctite 241 and use the **T8** Torx bit to tighten them to a tightening torque of **1 Nm** (see fig. 15).
- 16) Attach the bearing bushing in the foot stirrup, using the bearing bushing assembly/disassembly aid to do so (see fig. 4, item 2).
- 17) Set the bushing onto the assembly/disassembly aid and position it in the foot stirrup (see fig. 16, see fig. 17).

INFORMATION: Make sure to insert the bearing bushing into the foot stirrup without canting (see fig. 18).

7 Cleaning

Always avoid contact with water containing salt, chlorine or soap.

Should contact nevertheless occur or the product be dirty, clean it promptly.

- 1) Clean the product with a damp cloth and mild soap when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

8 Maintenance

The section is intended for qualified personnel.

INFORMATION

The product may be exposed to increased loads on a user-specific basis.

- Shorten the maintenance intervals according to the expected loads.

The manufacturer requires an inspection of the product at least every **6 months** to verify functionality and check for wear.

- **Replace the single components of the foot stirrup every 6 months using the 17AD316=22 service set.**

The product is designed to be free of lubricants.

Do not use lubricants.

8.1 Removing the foot stirrup

Required materials: 17AD312=22 assembly tool set, 17AD316=22 service set

- 1) Loosen the joint screw on the orthotic ankle joint and separate the foot component from the orthotic ankle joint.
- 2) Remove the bearing bushing from the foot stirrup using the assembly/disassembly aid (see fig. 4, item 2). Insert the assembly/disassembly aid into the bearing bushing and press it out of the foot stirrup.

- 3) To disassemble the spiral pins, use the assembly/disassembly aid for spiral pins from the 17AD312=22 assembly tool set (see fig. 4, item 1).
- 4) Insert the foot stirrup into the contour of the assembly/disassembly aid for spiral pins (see fig. 9)
- 5) Put on the counter piece (see fig. 10).
- 6) Remove the long cylinder pins (see fig. 11).
- 7) Then insert the long cylinder pins into the assembly/disassembly aid from the outside (see fig. 13).
- 8) Clamp the assembly/disassembly aid in the vice, turn the crank of the vice and press the spiral pins out of the foot stirrup (see fig. 14).
- 9) Take the assembly/disassembly aid out of the vice, then remove the foot stirrup, cylinder pins and spiral pins.
- 10) Loosen the fixing screws of the shim plates using the **T8** Torx bit and take the plates off the foot stirrup.

For assembly of the foot stirrup, see the section "Assembling the foot stirrup".

Please contact the manufacturer in case of questions.

9 Disposal

Dispose of the product in accordance with national regulations.

10 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

10.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

10.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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