
17S100=LB-22-T

DE	Gebrauchsanweisung	4
EN	Instructions for use	8
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	17
ES	Instrucciones de uso	21
PT	Instruções de utilização	26
NL	Gebruiksaanwijzing	30
SV	Bruksanvisning	35
DA	Brugsanvisning	39
NO	Bruksanvisning	43
FI	Käyttöohje	47
PL	Instrukcja użytkowania	51
HU	Használati útmutató	56
CS	Návod k použití	60
RO	Instrucţiuni de utilizare	64
HR	Upute za uporabu	69
SK	Návod na použitie	73
BG	Инструкция за употреба	77
TR	Kullanım kılavuzu	82
EL	Οδηγίες χρήσης	86
RU	Руководство по применению	91
JA	取扱説明書	95
ZH	使用说明书	100
KO	사용 설명서	103

17S100=LB-22-T

501S136=M6x12-1

(Lieferumfang von/
Scope of delivery for
17AD300=22-T)

INFORMATION

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INFORMATION

Date of last update: 2026-09-07

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide information on the processing of the 17S100=LB-22-T lamination bar.

2 Safety-related information

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

CAUTION

Excessive strain due to use on more than one user

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on one user only.
- ▶ Observe the maintenance recommendations.

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand).
- ▶ Do not expose the product to operating temperatures below **-10 °C** (14 °F) or above **+45 °C** (113 °F) (e.g. excessive sunlight, drying on a radiator).

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C** (570 °F).
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

3 Product description

3.1 Available sizes

- 17S100=LB-22-T

3.2 Maximum body weight

The product is approved for up to **125 kg**.

4 Intended use

4.1 Indications for use

The product is intended **exclusively** for orthotic treatment of the lower limbs of a user.

The product is suitable for unilateral or bilateral use.

Only use the product for the lamination or prepreg technique.

4.2 Target group

- Qualified personnel
- User

4.3 Qualification

Users may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4.4 Indications

- Partial or total paralysis of the leg muscles
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

4.5 Product combinations

The 17S100=LB-22-T lamination bar can be combined with the 17AD300=22-T orthotic ankle joint.

4.6 Lifetime

The product is designed for a maximum lifetime of **6 years** when used as intended and assembled professionally.

5 Preparing the product for use

The section is intended for qualified personnel.

5.1 Processing

CAUTION

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

Bending the lamination bar

Prerequisite: The orthosis model has been fabricated specifically for the user.

To avoid a risk of breakage, do not bend or otherwise modify the insertion zone on the bar. The allowable bending zone begins at the slotted hole (see fig. 1).

Required tools and materials: bending iron

- 1) Clamp a bending iron in a vice.
- 2) Bend the bar with a bending iron, noting the allowable bending zone.

5.2 Assembly

Mounting the bar

- 1) Insert the lamination bar into the base body of the joint (see fig. 2).
- 2) Secure it with the two flat head screws and a **T 30** wrench, securing the screw connections with Loctite 241 and a tightening torque of **8 Nm**.

6 Cleaning

Always avoid contact with water containing salt, chlorine or soap.

Should contact nevertheless occur or the product be dirty, clean it promptly.

- 1) Clean the product with a damp cloth and mild soap when needed.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

7 Maintenance

The section is intended for qualified personnel.

INFORMATION

The product may be exposed to increased loads on a user-specific basis.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires an inspection of the product at least every **6 months** to verify functionality and check for wear.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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