beurer

MultiSense

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ENGLISH

Read these instructions for use carefully. Follow the warnings and safety notes. Keep these instructions for use for future reference. Make the instructions for use accessible to other users. If the device is passed on, provide the instructions for use to the next user as well.

Contents

1. Included in delivery 20
2. Signs and symbols 20
3. Intended Purpose 22
4. Warnings and safety notes 22
5. Information on measuring body temperature 24
 5.1 Influences on body temperature 24
 5.2 Influences from handling 24
 5.3 Influences on ear measurement 25
 5.4 Influences on forehead measurement 25
6. Device description 25
7. Initial use 26
8. Settings 26
9. Usage 26
 9.1 Measuring body temperature in the ear 27
 9.2 Measuring body temperature on the forehead 27
 9.3 Measuring the surface or room temperature 27
 9.4 Evaluating the results 28
 9.5 Stored measurements 28
10. Cleaning and maintenance 28
11. Accessories and/or replacement parts 29
12. Troubleshooting 29
13. Disposal 30
14. Technical specifications 30
15. Warranty/service 31

1. INCLUDED IN DELIVERY

Check that the exterior of the cardboard delivery packaging is intact and make sure that all contents are included in the delivery. Before use, ensure that there is no visible damage to the device and accessories and/or replacement parts, and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

- Thermometer
- Batteries, see chapter “Technical specifications”
- Storage box
- Forehead measurement cap
- Instructions for use
- Quick start guide

2. SIGNS AND SYMBOLS

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

⚠ DANGER

Indicates an imminent danger. If it is not avoided, it will result in death or serious injury.

⚠ WARNING

Indicates a potentially impending danger. If it is not avoided, there is a risk of death or serious injury.

⚠ CAUTION

Indicates a potentially impending danger. If it is not avoided, slight or minor injuries may result.

NOTICE

Indicates a potentially harmful situation. If it is not avoided, the device or something in its vicinity may be damaged.

Observe the instructions
Read the instructions before starting work and/or operating devices or machines

The electronic device must not be disposed of with household waste

Do not dispose of batteries containing hazardous substances with household waste

Manufacturer

Date of manufacture

CE labelling
This product satisfies the requirements of the applicable European and national directives.

Marking to identify the packaging material.
A = material abbreviation, B = material number:
1–7 = plastics, 20–22 = paper and cardboard

Separate the product and packaging elements and dispose of them in accordance with local regulations.

Device protected against foreign objects ≥12.5 mm and against water dripping at an angle

Direct current

Battery

Unique device identifier (UDI)
Identifier for unique product identification

Article number

Serial number

Swiss authorised representative

Medical device

Type number

Type BF applied part

Temperature limit

	Humidity limitation
	Importer symbol
	Protect from sunlight and heat
	Manufacturer in accordance with the Packaging Regulation (PPWR)

3. INTENDED PURPOSE

Intended Use

The thermometer is intended exclusively for measuring human body temperature in the ear canal (ear temperature mode) or on the forehead (forehead temperature mode). The thermometer is only intended for private use at home, not for commercial use.

Target group

The thermometer is intended for medical care in a home environment, not in professional healthcare facilities. It is suitable for use on people of all ages. Self-administration is possible from the age of 12.

Clinical benefits

The thermometer can be used in a home environment to measure the body temperature in the ear canal or on the forehead. The temperature is monitored in °C (or °F) to one decimal place.

Indications

The thermometer allows the user to quickly and easily measure human body temperature in the ear canal or on the forehead. This enables the body temperature to be monitored.

Contraindications

The measurement must not be taken

- in an ear affected by inflammatory diseases (e.g. discharging pus or secretion)
- after possible ear injuries (e.g. eardrum damage)
- in the healing period after surgery
- in an ear with pathological narrowing of the ear canal

4. WARNINGS AND SAFETY NOTES

General notes

⚠ WARNING

- Keep packaging material away from children. There is a risk of suffocation.
- Keep the device away from children, animals, and vermin.
- The measurements you take are for your information only – they are not a substitute for a medical examination! Discuss your measured values with your doctor and never make your own medical decisions based on them (e.g. regarding medicine doses).
- If you obtain unusual temperature results, please contact your GP straight away.
- The device must not be used on ears where the ear canal is too small.
- Children must not be allowed to use the device. Medical devices are not toys.
- Do not modify the device.
- The surface temperature measured may be different from the actual temperature inside the object being measured, e.g. in the case of baby bottles.

⚠ CAUTION

- If the device has been stored at temperatures below 15 °C (59 °F) or above 40 °C (104 °F), allow it to acclimatise for about 60 minutes at an ambient temperature of 15 °C - 40 °C (59 °F - 104 °F) before use.
- This thermometer is a sensitive electronic device. Do not drop the device or subject it to any mechanical impacts.
- Do not expose the device to direct sunlight, extreme temperatures or high humidity.
- Do not expose the device to liquids. The device is not waterproof. Avoid all direct contact with water or other liquids.
- Portable and mobile HF communication systems may interfere with this device
- Do not touch the measuring sensor with your fingers.

NOTICE

- Measure your temperature regularly to obtain your normal average temperature. If you suspect you have a fever, simply compare your normal average temperature with the temperature measured currently.
- The thermometer is only designed for the measuring area on the human body stated in the instructions for use.
- The thermometer is calibrated during production. If the thermometer is used according to the instructions for use, further calibration is not necessary.
- Each time before using the device, check that the measuring sensor is intact. If it is damaged, please contact your retailer or the service address.
- Should you have any questions about using the device, please contact either your retailer or Customer Services.
- The manufacturer is not liable for damage resulting from improper or incorrect use.
- Repairs must only be carried out by authorised service centres, otherwise the warranty is invalidated.

Notes on handling batteries

⚠ DANGER

- **Choking hazard!** Keep batteries out of the reach of children. Seek medical attention immediately if swallowed. Swallowing them may cause burns, severe internal injuries, and death.
- Never allow children to replace batteries without adult supervision.

⚠ WARNING

- **Risk of explosion! Risk of fire!** Failure to comply with the following points can result in personal injury or cause overheating, leakage, venting, breakage, explosion, or fire on the battery.
- This device contains non-rechargeable batteries which must not be charged.
- Do not throw batteries into a fire.
- Never charge, forcibly discharge, heat, disassemble, open, crush, deform, encapsulate, or modify batteries.
- Never short-circuit batteries or battery compartment contacts.
- Protect the batteries from direct sunlight, rain, heat, and water.
- Exposure of batteries to an environment with extremely high temperatures or an extremely low air pressure may result in explosion or leakage of flammable liquids and gases.
- Dispose of defective and discharged batteries immediately and properly (see chapter on disposal).
- Do not use modified or damaged batteries.

- Always select the correct battery type.
- Always insert the batteries correctly, taking into account the polarity (+ / -).
- Never mix batteries of different manufacturers, capacities (new and used), size, or type within a device.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- If fluid from a battery comes into contact with your skin or eyes, wash the affected areas with water and seek medical assistance.

⚠ CAUTION

- Store batteries in a well-ventilated, dry, and cool place in a non-conductive container in which the batteries cannot be short-circuited to each other or by other metal objects.
- Keep batteries clean and dry.
- Keep batteries away from water.
- If the device is not going to be used for a long period of time, remove the batteries from the battery compartment.

NOTICE

- Do not use rechargeable batteries.

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The device may not be fully usable in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment, such as antenna cables or external antennas) at least 30 cm away from all device parts.
- Failure to comply with the above could impair the performance of the device.

5. INFORMATION ON MEASURING BODY TEMPERATURE

5.1 Influences on body temperature

- A person's individual metabolism
- Age: Body temperature is higher in infants and toddlers than in adults. Greater temperature fluctuations occur faster and more often in children. Normal body temperature decreases with age.
- Clothing
- Outside temperature
- Time of day: Body temperature is lower in the morning and increases throughout the day towards evening.
- Activity: Physical and, to a lesser extent, mental activities increase body temperature.

5.2 Influences from handling

- Temperatures measured with different thermometers should never be compared with one another. The temperature measurement varies depending on the part of the body where the

measurement is taken. Take this into account when interpreting the measurement results.

- The device and the patient must remain in the room where the measurement is to be performed for at least 30 minutes before the measurement.
- After physical activity, after taking a bath, or after being outdoors, wait approximately 30 minutes before taking a measurement.
- Holding the thermometer in your hand for too long before taking a measurement can cause the device to warm up. This means the measurement could be incorrect.
- Repeat the measurement if you are unsure of the measured values.

5.3 Influences on ear measurement

- Some people produce different readings in their left and right ear. In order to record temperature changes, always measure a person's temperature in the same ear.
- This thermometer may only be used without single-use protective caps.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- Ensure that the ear is clean and dry. Clean it if necessary before taking the measurement.

5.4 Influences on forehead measurement

Remember:

- Physical activity, increased perspiration on the forehead, taking vasoconstrictive medication, wearing head coverings in the forehead area, staying outdoor and skin irritations can distort the measurement.
- The forehead must be uncovered and free from sweat, hair, cosmetics, and dirt. Expose or clean the forehead and wait 10 minutes to improve measurement accuracy.
- Keep the thermometer still during measurement. Moving the thermometer during this may cause inaccurate measurements.

NOTICE

If you are uncertain about interpreting the results or if the values are abnormal (e.g. fever), consult your doctor. This also applies in the case of slight temperature changes if there are other symptoms of illness such as agitation, heavy sweating, flushed skin, fast pulse rate, tendency to collapse, etc.

6. DEVICE DESCRIPTION

The associated drawings are shown on page 4+5.

Buttons		Button press	Functions
1		short	Setting the ear and forehead mode E
		long	Setting the surface and room mode E
2		short	Displaying stored measurements
		long	Acoustic signal on/off
3		short	Starts the temperature measurement/Switches the device on
		long	Switches the device off

4 LED ring

5 Battery compartment cover

6 Measuring tip

7 Measuring sensor

8 Forehead measurement cap

Display description

9 Forehead mode

10 Ear mode

11 Surface and room mode

12 Acoustic signal on / off

13 Battery level indicator

14 Measurement unit in degrees Celsius (°C)

15 Distance display in centimetres (cm)

16 Distance display in inches (in)

17 Measurement unit in degrees Fahrenheit (°F)

18 Memory space number

19 Cleaning indicator for measuring tip and sensor

20 Distance indication in forehead mode

21 Temperature value

7. INITIAL USE

Inserting/replacing batteries

1. Remove the battery compartment lid **5**.
2. Insert the batteries (see chapter “Technical specifications”) making sure the polarity is correct according to the label.
3. Close the battery compartment cover.

The device automatically switches on and the measurement unit flashes in the display. Select

the measurement unit (°C or °F) by briefly pressing / or / and confirm with .

When the batteries become flat, appears. It is still possible to measure temperature. The batteries must be replaced soon. When and appear, the batteries must be replaced. If the batteries are too flat, the device will switch off automatically.

Switching the thermometer on/off

To switch on the device, briefly press . All display elements are briefly displayed and a short beep sounds. The device is ready to measure the temperature.
The device automatically switches off after approx. 60 seconds or you can press and hold .

8. SETTINGS

Measurement unit

To set the measuring unit, see chapter “Initial use”.

Acoustic signal

To switch the acoustic signal on or off, press and hold /. When switched on, a short acoustic signal sounds and is displayed. When the acoustic signal is off, is displayed. In case of error messages, two long beeps will always sound.

9. USAGE

NOTICE

Before taking the measurement, observe the notices in the chapter “Information on measuring body temperatures”.

9.1 Measuring body temperature in the ear

Before starting the measurement, check that the positioning is correct A. If a measurement is taken with incorrect positioning, a valid measurement result cannot be guaranteed.

1. Remove the forehead measurement cap to automatically switch from forehead mode to ear mode E. When in surface and room mode, briefly press / to switch to ear mode. Make sure that the device is in ear mode. You can tell this from the symbol on the display.
 Tip: To make it easier to remove the forehead measurement cap, press it together slightly.
2. Make sure that the measuring tip and also the ear canal are clean. Insert the measuring tip carefully into the ear A and short-press . The charging symbol D is shown on the display during the measurement.
3. The end of the measurement is signalled by an acoustic signal, if the acoustic signal is switched on. The measured value appears on the display and the LED ring shows a coloured classification.
4. Press briefly to take another measurement.

9.2 Measuring body temperature on the forehead

Before starting the measurement, check that the measuring distance is correct as well as the measurement range B.
If a measurement is taken outside the specified range, a valid measurement result cannot be guaranteed.

1. Place the forehead measurement cap on the measuring tip to automatically switch from ear mode to forehead mode E. When in surface and room mode, briefly press / to switch to forehead mode. Make sure that the device is in forehead mode. You can tell this from the symbol on the display.
2. Hold the thermometer 1 cm (0.4 in) in front of the measurement point in the centre of the forehead B. Before the first forehead measurement, the required distance is shown on the display by . The distance is displayed in centimeters (°C) or inches (°F), depending on the selected unit. Press briefly to start the measurement. The charging symbol D is shown on the display during the measurement.
3. The end of the measurement is signalled by an acoustic signal, if the acoustic signal is switched on. The measured value appears on the display and the LED ring shows a coloured classification.
4. Press briefly to take another measurement.

9.3 Measuring the surface or room temperature

1. With the device switched on, press and hold / until appears on the display E.
 Measuring surface temperatures on objects:
 Hold the thermometer 1 cm (0.4 in) in front of the desired measuring point C.
 Measuring room temperatures:
 Keep the thermometer out in the room whose temperature you want to measure. Make sure that there are no objects in the vicinity that could falsify the measurement result.
2. Press briefly to start the measurement. The charging symbol D is shown on the display during the measurement.
3. The end of the measurement is signalled by an acoustic signal, if the acoustic signal is switched on. The measured value appears on the display. There is no classification of the measured values.

4. Press briefly to take another measurement.

9.4 Evaluating the results

Indicator		Temperature range	Temperature classification
LED ring	Acoustic signal*		
Red	3 short beeps	32.0 °C - 34.9 °C (89.6 °F - 94.8 °F)	Low temperature
Green	1 long beep	35.0 °C - 37.4 °C (94.9 °F - 99.3 °F)	Normal temperature
Orange	3 short beeps	37.5 °C - 37.9 °C (99.4 °F - 100.3 °F)	Elevated temperature
Red	3 short beeps	38.0 °C - 43.0 °C (100.4 °F - 109.4 °F)	Fever

* if acoustic signal is switched on

9.5 Stored measurements

Display measurements

The device automatically stores the values from the last 30 measurements. If the 30 memory spaces are full, the oldest value on memory location 01 is overwritten.

1. With the device switched on, briefly press /. The latest measurement is displayed. The measured value, the measurement mode, the coloured classification (if applicable) and the memory space number are displayed.
2. Press / briefly to scroll down through the memory spaces, and press / briefly to scroll up.
3. To exit the memory of measured values, press briefly.

Delete measurements

The memory of the measured values can be deleted. All measurements will be erased during the deletion process.
The measured values can be deleted as follows:

1. Call up the measured value memory (see chapter “Display measured values”).
2. Press and hold / and / at the same time until CLr appears on the display.
3. Press briefly to delete the values. To cancel the deletion process, briefly press / or /.

10. CLEANING AND MAINTENANCE

Cleaning the device

NOTICE

- The measuring sensor is the most sensitive part of the thermometer. Be very careful with the measuring sensor when cleaning the device.
- Do not use any harsh cleaning products. Antiseptics, warm water, thinner and benzene are not permitted.
- Always observe all safety notes for users and the device. (“4. Warnings and safety notes”)

Clean the measuring sensor regularly, at least after each application in the ear. After the ear measurement, CLEAN ME appears on the display to remind you to clean the measuring sensor before switching off the device. Use a cotton bud that is lightly moistened with Isopropyl alcohol (70%) for this. Use it to carefully remove dirt from the measuring sensor and measuring tip. Carefully dry the cleaned areas and allow the device to dry for at least 45 minutes after cleaning.

Inadequate or improper cleaning of the measuring sensor can influence the measurements. To clean the entire device, use a soft cloth slightly moistened with a mild soapy solution. Under no circumstances may liquid enter the device. Do not use the device again until it is completely dry.

Storing the device

NOTICE

The device must not be stored or used at an excessively high or low temperature or humidity („14. Technical specifications“), in sunlight, in association with an electrical current or in dusty locations. Measurement inaccuracies can otherwise occur. Always store the device in the storage box.

11. ACCESSORIES AND/OR REPLACEMENT PARTS

Accessories and/or replacement parts are available at www.beurer.com, under “Service”.

Designation	Item number
Forehead measurement cap white	110.346

12. TROUBLESHOOTING

In case of error messages, two long beeps will always sound and the following displays will appear:

Display	Cause	Solution
	Room temperature below 15 °C (59 °F) or over 40 °C (104 °F).	Room temperature must be between 15 °C - 40 °C (59 °F - 104 °F). Allow the thermometer to rest in a room for at least 30 minutes at room temperature.
	The thermometer is placed incorrectly or unsteady.	Follow the instructions in the chapter “Usage” thoroughly. Repeat the measurement.
	The temperature recorded is higher than: 1. Forehead/ear mode: 43 °C (109.4 °F), 2. Surface/room mode: 80 °C (176 °F)	Operate the thermometer only within the specified temperature ranges. Read instructions for use thoroughly and take new temperature measurement. In the event of a repeated error message in body temperature measurement, contact your doctor. In the event of a repeated error message in surface/room mode, contact your retailer or Customer Services.
	The temperature recorded is lower than: 1. Forehead/ear mode: 32 °C (89.6 °F), 2. Surface/room mode: 10 °C (50 °F).	
	The batteries are flat.	Replace the batteries.

13. DISPOSAL

For environmental reasons, do not dispose of the device in household waste at the end of its service life. Dispose of the device at a suitable local collection or recycling point in your country. Dispose of the device in accordance with EC Directive Waste Electrical and Electronic Equipment (WEEE). If you have any questions, please contact the local authorities responsible for waste disposal.

You can obtain the location of collection points for old devices from the local or municipal authorities, local waste disposal companies or your retailer, for example.

Empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

14. TECHNICAL SPECIFICATIONS

Model	MultiSense
Type	FT 66
Measurement method	Infrared measurement
Basic functions	Forehead temperature measurement Ear temperature measurement Surface and room temperature measurement
Measurement units	Celsius (°C) or Fahrenheit (°F)
Operating mode	Comparison mode
Operating conditions	15 °C - 40 °C (59 °F - 104 °F) at a relative humidity < 90%, 700 -1060 hPa ambient pressure
Storage and transport conditions	-25 °C - 55 °C (-13 °F to 131 °F) at a relative humidity < 90%
Measuring distance	Forehead and surface mode: 1 cm (0.4 in)
Minimum measurement duration	2 seconds
Measurment range	Ear- and forehead mode: 32 °C - 43 °C (89.6 °F - 109.4 °F) Surface and room mode: 10 °C - 80 °C (50 °F - 176 °F)
Memory function	Automatically stores the last 30 measurements.
Display	Inverse display
Display resolution	0.1 °C/°F
Energy-saving functions	The device switches off automatically after 60 seconds.
Dimensions	L 161 mm x W 41 mm x H 47 mm
Weight	76 g (without batteries)
Battery 2 x1.5V AAA (LR03)	Batteries last for approx. 1000 measurements
IP protection	IP22

Expected useful life	Information on the service life of the product can be found on the website.
Classification	Internally powered ME device. Continuous operation. BF applied part (measuring tip)

Clinical measurement accuracy

Mode	Clinical bias	Limit of agreement	Clinical repeatability
Ear	-0.14 °C (-0.25 °F)	0.80 °C (1.44 °F)	0.17 °C (0.31 °F)
Fore-head	-0.10 °C (-0.18 °F)	0.98 °C (1.76 °F)	0.19 °C (0.34 °F)

Clinical test according to ISO 80601-2-56. Thermometer with adapted mode (forehead/ear mode) with oral equivalence measurement.

Laboratory measurement accuracy

Mode	Measuring temperature range	Environmental conditions	Accuracy
Ear / Fore-head	34 °C - 42 °C (93.2 °F - 107.6 °F)	15 °C - 40 °C (59 °F - 104 °F)	± 0.2 °C (± 0.4 °F)
	For other measurement and ambient temperature ranges		± 0.3 °C (± 0.5 °F)
Surface / Room	34 °C - 42 °C (93.2 °F - 107.6 °F)	15 °C - 40 °C (59 °F - 104 °F)	± 0.2 °C (± 0.4 °F)
	For other measurement and ambient temperature ranges		± 4% or ± 2 °C (4 °F), whichever is greater

The serial number is located on the device or in the battery compartment.
We reserve the right to make technical changes to improve and develop the product.
The device conforms with the respective national regulations and European standard EN 60601-1-2 (Group 1, Class B, in accordance with CISPR-11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8, IEC 61000-4-39) and is subject to particular precautions with regard to electromagnetic compatibility.

15. WARRANTY/SERVICE

Further information on the warranty and warranty conditions can be found in the warranty leaflet supplied. Before submitting a claim, please check the batteries and replace them if necessary.

Notification of incidents

For users/patients in the European Union and identical regulation systems, the following applies: If a major incident occurs during or through use of the product, notify the manufacturer and/or their representative of this as well as the respective national authority of the member state in which the user/patient is located.

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CH REP Beurer Schweiz AG, Hauptstrasse 77, CH-6260 Reiden
UK: BEURER UK Ltd., Suite 16, Stonecross Place, Stonecross Lane North,
Beurer GmbH • Söflinger Str. 218 • 89077 Ulm, Germany • www.beurer.com