

GlucoMemory software in German and English: simple overview of results and trends in target range

Including "beurer HealthManager" app for download

"beurer HealthManager" software for download free of charge at: www.beurer.com > SERVICE > Download

Required volume of blood: 0.5 microlitres

Duration of measurement: approx. 5 seconds

Code-free

Measured value memory for 900 measurements with date, time and measurement marking  $\,$ 

Calculating the average value for 7, 14, 30 or 90 days for - all measurements - marked measurements

Blood sample: fresh whole blood (capillary) or venous whole blood

Plasma-calibrated

Measurement range: 20 - 630 mg/dL / 1,1 - 35,0 mmol/L

Fully fulfils EN ISO 15197:2015 (compulsory since May 2016)

Automatic switch-off

Product dimensions (L/H/D): 95 x 53 x 15 mm

3L49\_2020-04-27\_03\_DS\_EN\_BEU Subject to errors and chang




GL 49 mg/dL & mmol/L

## beurer medical



## CONNECTION




Innovative networking between smartphone and blood pressure monitor – keep track of your measurements at home and on the move. Simple data transfer to a smartphone via Bluetooth®









• iOS ≥ 12.0 Android™ ≥ 8.0 • Bluetooth® ≥ 4.0





List of supported devices

GL49\_2020-04-27\_03\_DS\_EN\_BEU