



GL 48 mg/dL
& mmol/L

beurer
medical



Coloured target scale and display to help you to interpret the measurement



Strip ejection function for hygienic disposal



Extra wide test strips for easy insertion of the strips into the meter



Electronic diary available via web version incl. USB-uploader free of charge

Large, easily readable display

Extra wide test strips, easy measurement marking for optimised diabetes management: before eating, after eating, general marking

Blood volume check: measurement only possible with sufficient sample of blood

Starter set: GL 48 measuring device, LD 04 lancing device, 10 lancet needles, 10 test strips, USB cable, case, instructions for use, brief instructions, blood glucose diary

Ketone warning on display

GlucoMemory software in German and English: simple overview of results and trends in target range

PC software for download free of charge at: www.beurer.com > SERVICE > Download Center

Required volume of blood: 0.5 microlitres

Duration of measurement: approx. 5 seconds

Code-free

Measured value memory for 900 measurements with date, time and measurement marking

Calculating the average value for 7, 14, 30 or 90 days for
- all measurements - marked measurements

Blood sample: fresh whole blood (capillary) or venous whole blood

Plasma-calibrated

Suitable for self-testing

Measurement range: 20 - 630 mg/dL / 1,1 - 35,0 mmol/L

Fully fulfils EN ISO 15197:2015 (compulsory since May 2016)

Automatic switch-off

Product dimensions (L/H/D): 94 x 52 x 14 mm

Product weight: 50 g

Two button cell batteries (CR2032)



GL 48 mg/dL
& mmol/L

beurer
medical

Medical device, 5-year guarantee

Pack: 5 / Shipping carton: 4 x 5

EAN no.:

4 211125 463663 (mg/dL)

4 211125 463670 (mmol/L)


Item no.:

463.66 (mg/dL)

463.67 (mmol/L)



CONNECTION



* tested online application "beurer HealthManager" Cloud

GL48_2022-02-16_05_EN_BEU Subject to errors and changes