

**Touchless –**  
a real affinity with hands



 **DÜRR  
DENTAL**  
THE BEST, BY DESIGN

# Touchless universal dispenser

The failsafe Touchless 1000 (for 800 ml bottles) and Touchless 400 (for 500 ml bottles) sensor dispensers have several advantages – functionality, reliability and durability. This is important as hands-free handling prevents cross-contamination. The dispensers provide everything required to ensure smooth procedures in a dental practice.

Touchless T 1000/400 dispensers have been developed from the tried-and-tested universal dispensers with a lever for manual handling. With its failsafe sensor technology, the touchless version retains the advantages of the manual dispenser but with the following added benefits:

- Optimum hand hygiene thanks to contactless operation
- Battery operation with battery power indicator by means of an audible battery replacement signal
- User-friendly dosing technology





### Autoclavable pump

Easy replacement from the front.

### Flexible sealing cap

Ensures that the bottle sits firmly in the dispenser.

### High standards of hygiene

The case is made from anodised aluminium - its bacteriostatic effect hinders the growth of bacteria.

#### Technical data

Battery voltage 7.5 V DC

Battery type (5x) mono cell ANSI D (LR 20/1.5 V)

Battery service life up to 2 years

Dosage (adjustable):

0.75 ml/1 ml/1.25 ml/1.5 ml

Range of detection of fail-safe sensor:

4 cm from top and bottom

Weight without batteries 1250 g

Dimensions W x H x D: 9.2 x 33.0 x 21.0 cm

# Accessories



## Drip tray

With stainless steel bracket and removable plastic dish.



## Table stand


Made of sturdy plastic for the flexible use of the Touchless universal dispenser.



## HD products from Dür System Hygiene

- HD 410 Hand disinfection (500 ml or 800 ml)
- HD 435 Washing lotion (500 ml or 1 l)
- HD 440 Care lotion (500 ml)

**Viruses destroyed? Yes!**  
**Skin damaged? No!**  
**Hand disinfection**  
**must be gentle.**



To the usage video!

DÜRR DENTAL AG  
Höpfigheimer Strasse 17  
74321 Bietigheim-Bissingen  
Germany  
[www.duerrdental.com](http://www.duerrdental.com)  
[info@duerr.de](mailto:info@duerr.de)

