



New matrix –
excellent look

Provicol® QM Aesthetic

Provicol® QM Aesthetic

TRANSLUCENT, TEMPORARY LUTING CEMENT
WITH CALCIUM HYDROXIDE, EUGENOL-FREE

Provicol® QM Aesthetic

TRANSLUCENT, TEMPORARY LUTING CEMENT WITH CALCIUM HYDROXIDE, EUGENOL-FREE

Indications

Temporary luting of provisional or definitive inlays, onlays, crowns, partial crowns / veneers and bridges


Temporary filling of small one-surface cavities

Advantages

- Highly aesthetic, translucent appearance
- Increased strength for more retention and secure adhesion
- Optimum flowability for easy cementation
- Pinpoint stability allows precise removal of excess material
- Accurate fit due to very low film thickness
- Eugenol-free: no inhibition of curing reaction of resin-based filling or luting materials
- Free from colophony: low allergy potential
- Contains calcium hydroxide and zinc oxide
 - Calcium hydroxide supports the formation of tertiary dentine und has bacteriostatic effect
 - Zinc oxide is known for its antibacterial effect
- Radiopaque

Presentations

- REF 1108 QuickMix syringe 5 ml, accessories
- REF 1109 Set II QuickMix syringe 3 x 5 ml
(1 x Provicol QM, 1 x Provicol QM Plus, 1 x Provicol QM Aesthetic), accessories
- REF 2189 Mixing tips type 10, 50 pcs.



Provicol QM Aesthetic reaches a translucency, which is usually only reached by resin-based products. At the same time, Provicol QM Aesthetic combines all advantages of a luting cement and is compatible with a wide range of temporary and permanent restoration materials.

Source: internal measurement

Comparison translucency



Visual comparison of translucency and shade – Provicol QM Aesthetic also meets the demands in the anterior region, which were previously reserved for resin-based products.

