

TECHNICAL SPECIFICATIONS

Materials

- Transparent shield: PET
- Headband: 50% Polyester and 50% Elastane
- Protector foam: Polyurethane
- Rivet: Acetal resin
- Individual packaging: PE zip pouch, re-sealable
- Secondary packaging: Cardboard box

Main dimensional properties

- Shield width: 325 mm
- Shield total length: 230 mm
- Visor length: 210 mm
- Visor thickness: 400 μm .
- Shield depth: 40 mm
- Protective foam length: 225 mm

Biocompatibility

- ISO 10993-5: Non cytotoxic
- ISO 10993-10: Non sensitizing
- ISO 10993-10: Negative intracutaneous irritation test

Made with hypoallergenic materials.

Sterilization

- Non sterile product.

Storage and maintenance

- Keep protected from direct sunlight / the outside.
- Keep in a dry place.
- Single use product. No maintenance required.
- Limitation of temperature: 5°C – 40°C

Packaging

- Unitary PE zip pouch.

