







### THE FASIEST WAY TO CREATE NATURAL GLOSS

The innovative shape with the flexible polishing lamellae creates a natural gloss on all zirconia oxide restorations. No matter if the surfaces are level, convex or concave, the diamond-impregnated lamellae polishes them to a high-gloss level.

### **BENEFITS**


- Flexible polishing twist for refining all types of tooth surfaces
- Highly suitable for occlusions
- Modellation of zirconia restoration is maintained
- Sterilizable and reusable
- Economic efficiency due to high durability and reusability


**SCAN TO VIEW PRODUCT IFUs** 

















#### **RECOMMENDED AND MAXIMUM ROTATION SPEEDS**

Recommended speed: 7.000 – 10.000 min for Coarse/Medium/Fine Maximum speed: 20.000 rpm

# TWIST™ DIA for Zirconia

Zirconia finishing and polishing system

#1302-IP COARSE 14mm (10pcs) #1305-IP FINE 14mm (10pcs)

#1303-IP MEDIUM 14mm (10pcs) #1301-IP Introductory Kit 14mm (3pcs)

# POLISH COMPOSITE WITH TWIST™ DIA

The innovative shape with flexible polishing spirals offers various application benefits to the dentist for excellent polishing results. Find out more: kuraraynoritake.eu/twist-dia

Before using this product, be sure to read the Instructions for Use of the product.

The specifications and appearance of the product are subject to change without notice. Printed color can be slightly different from actual color.

