

STATMATIC[™] smart



STATMATIC smart¹
S33001

STATMATIC smart²
S33002

STATMATIC smart³
S33003

THE SMART SOLUTION FOR HANDPIECE CLEANING

The cleaning and maintenance unit *STATMATIC* smart has been specifically designed to meet the challenges created by complex designs of modern handpieces.



Easy To Use

The display of the *STATMATIC* smart shows program details and allows for easy operation.



Efficient

STATMATIC smart takes care of a thorough internal cleaning and maintenance of dental instruments like contra angles and turbines. In the standard program "Cleaning and maintenance" the automated cleaning and maintenance process takes only about 3 minutes per instrument. During the gentle cleaning process a special, non fixing liquid cleaner efficiently removes bacteria, germs and other dirt from the water and air channels as well as the hollow cavities. The maintenance of the internal components helps to ensure the longevity of the instruments.



Flexibility

Four different cleaning and maintenance programs offer flexibility and correspond to the different requirements of instruments and offices as well as to the individual procedures.

In addition the unit is equipped with a separate tip for the specific maintenance of the chuck, in order to maintain its optimal retention force.



Capacity

STATMATIC smart is designed for three instruments and is available in three different versions, with 1, 2 or 3 fixed motorized stations for contra angles. Turbine adapters are available separately.



Compact Design

With it's compact design *STATMATIC* smart saves valuable space in the hygiene room. For convenience, the USB-Stick and test adapter can also be stored space-saving within the unit, under its cover.



Convenient

During the automated processes cleaner and oil are dosed exactly and economically so that dosing mistakes and costs can be reduced.



Safe

For safety and traceability purposes all cycles are documented automatically and can easily be transferred to a computer by using the USB key provided.

SciCan | Your Infection Control Specialist[™]

www.scican.com

SD-456 EU EN FR R1. STATMATIC and Your Infection Control Specialist are trademarks of SciCan Ltd. Manufactured by: SycoTec GmbH & Co. KG Wangener Strasse 78, 88299 Leutkirch Germany.