

THE PERMANENT CROWN FROM THE 3D PRINTER

saremco print - CROWNTEC

THE „GAME CHANGER“

Teeth and denture base
customized from the
3D-printer with
CROWNTEC and
DENTURETEC*

CROWNTEC is a qualitatively outstanding composite resin for the manufacturing of highly biocompatible permanent restorations.

- for denture teeth
- for permanent crowns, inlays, onlays, and veneers
- highly esthetic one-phase material with perfectly matched opacity
- outstanding biocompatibility: no elutable substances in an aqueous environment
- ready to use after minimal shaking time
- custom staining possible with saremco els paintart*

*see separate product information for DENTURETEC and els paintart (CE marking)



Technical Data

Consistency	medium viscosity, homogeneous
Color	according to VITA® standard
Depth of cure (DIN EN ISO 4049)	≥ 1.5 mm
Flexural strength (DIN EN ISO 4049)	> 135 MPa (average)
E-Modulus (DIN EN ISO 4049)	> 4000 MPa
Barcol hardness	> 40
Shelf life	3 years
Medical device	class IIa

Order information (500g bottle)

REF	Product Description
8063	CROWNTEC, A1
8052	CROWNTEC, A2
8051	CROWNTEC, A3
8065	CROWNTEC, B1
8066	CROWNTEC, SW (bleach-color)

Following shades per request: A3.5, A4, B2, B3, C3, C4, D3

