

PERFORMANCE. OPERATIONAL EASE. SAFETY.

Obtain high-quality images consistently, while improving safety

Compact, intuitive and easy to use, the CS 2200 offers the superior performance you need while providing improved safety for your patients and staff—year after year. The unit's high frequency DC X-ray generator ensures delivery of the image quality you're looking for, while the compact, lightweight design accommodates the space requirements of almost any practice.



We are the only dental company to engineer and manufacture X-ray generators from start to finish, enabling us to maintain complete control over the standards and quality of each system.

CS 2200: EXCEEDING YOUR EXPECTATIONS FOR 2D IMAGING



MAXIMUM USABILITY

The intuitive, handheld timer allows users to quickly select correct settings or choose from pre-programmed options



HIGH IMAGE QUALITY

With a focal spot of 0.7 mm, the system produces sharp, detailed images for diagnostic excellence



CONFIGURATION VERSATILITY

With multiple installation configurations available, the CS 2200 system easily adapts to your practice's needs**



MINIMAL RADIATION EXPOSURE

The CS 2200 system automatically determines the optimal dose, maximizes the amount of useful x-rays and generates up to 30% less radiation*



IMAGING COMPATIBILITY

Compatible with both film and digital imaging applications, the CS 2200 system works the way you do



PRECISE POSITIONING

The stable arm and lightweight tube head ensure correct, effortless alignment and reduce the risk of cone-cutting images



IMAGE OPTIMIZATION

Equipped with adjustable tension (60 or 70 kV), the CS 2200 system lets you accentuate contrast and levels of gray shading, according to your diagnostic needs



* When compared to conventional x-ray generators

**The CS 2200 system's wall mount was specifically designed to facilitate the replacement of the former Trophy Irix unit to minimize the impact and cost to your practice.

CS 2200 INTRAORAL X-RAY SYSTEM

© Carestream Health, Inc. 2017. 14881 DE 22 SS 1116

carestreamdental.com

800.944.6365