




avant

avant system: **the time you save is yours**

dentalart

what
will you do with the
time you save?



A more efficient way of working

The Avant system comprises a series of kits that enable precise and coordinated work, improving your performance and day-to-day organizational effectiveness.

Contamination reduction

Equipment and consumables remain separate from each other before, during and after surgical operations.

Add functionality to Dental Art furniture and discover a way of working and comfort level you've never experienced before.
The Avant system enhances your professionalism and the safety of your associates and patients.

Time saving in many ways

Everything you need, close to hand – Avant follows your movements, helping you work better and faster.

Better occupational health

Thanks to the ideal height of the worktop and operational area, you'll avoid bending and stretching, and will be able to maintain proper posture of shoulders, neck and back while you work.

Better organizational effectiveness of your dental practice

A predefined and coordinated workflow saves time and optimizes functionality for the entire workspace. Eliminating unforeseen circumstances and interruptions helps you maintain working concentration.




All the materials and equipment needed for procedures are prepared in advance and moved to the operational area.

Trays holding sterile instruments are arranged in the wall cabinet or directly on the tray.

avant_ obstacle-free organized workspace

Infected material and instruments are collected for disposal and/or cleaning.





Only the material necessary for the procedure is moved from the storage wall cabinet to the large bottom drawer, protected but still at hand.

The tray follows each movement and keeps the surrounding space free.

avant_use one or all of the three-module system:

avant is flexible

Everything at hand

Right- or left-handed, two- or four-handed. The mobile tray follows your movements and the instruments are always there, at hand, right where you look for them.

A

flexdesk®

Flexdesk® suction system

The ability to choose the optimal suction system position for every process allows you to maintain a better posture, facilitating your movements.

B

flexclip®

Managing the risk of contamination of consumables

The Tubs drawer, in its lowered position, is accessible while you are seated and keeps material protected from contamination during procedures.

C

flexwall®

Avant is an open system

Avant is a series of modules you can adapt to Dental Art furniture. Choose the size and modularity that suits you best.

+




A

flexdesk®

rotation 360° •
excursion 195 cm •
extension 70 cm •
payload 15 kg •

flexdesk[®]

full-length sliding desk



A

flexdesk®


The flexdesk® module includes a Corian® tray and movement mechanism integrated into a Corian® worktop.

Socket on the tray arm, to power small equipment (optional).



flexdesk®

fluid movement for flexible working positions



flexdesk[®] freedom at last



flexdesk[®]

Only where and when it's needed
without using mobile cabinets
or other worktops, flexdesk[®]
increases the available space
around you and your patient.

4 hands

Flexdesk® is easily repositionable, always attuned to the habits of the doctor-assistant team.



2 hands

Flexdesk® enhances the solo work of doctors and hygienists.




Bring the suction system into the cabinet,
maintain proper posture.

B

flexclip®

With flexclip® the suction system
is no longer located inside the
dental chair, but is integrated into
the cabinet.

A fixed-position suction system can
require assistants to make awkward
and uncomfortable movements.



The hoses of the suction
system are integrated
and grouped together
inside the cabinet.
Flexclip® is compatible
with the leading suction
systems on the market.



right-handed



left-handed

Flexclip® is easily
repositionable and is thus
adaptable to any working
method, 2-4-handed,
right- or left-handed.

flexclip[®]

suction on flexdesk[®]

Hooked to the tray,
flexclip[®] follows its movement,
improving ergonomics,
flexibility and versatility.




C

flexwall®


height and
ergonomic
position

safety against
contamination of
materials

The drawer is accessible even while sitting and thanks to the split front its contents are visible and at hand at all times. Open, take and close is now easy, comfortable and quick – a genuine guarantee of safety.



Centered on the dental chair, the bottom drawer is very wide (100 cm). It can hold up to three Tubs and is always at hand whatever the working position.




With flexwall® the consumables in use are stored in the large bottom drawer where they are protected at all times from environmental contamination.


Flexwall® can be used solo or in combination with different storage modules, ranging from conventional units to ones tailored for Tubs & Trays.

flexwall®


elements



flexwall-one® W17S*
(single distribution module)
middle drawer with "split" front for
housing up to 3 Tubs.



flexwall® W17*
wall cabinet with 2 height-
adjustable glass shelves.
Housing up to 6 Tubs.



flexwall® W17TUBS*
wall cabinet arranged for
housing 9 Tubs.

Vertical column
arranged for housing 8 Tubs
(V16TUBS) or 16 Trays with the
following dimensions:
W18.3 H3.9 D28.4cm (V16TRAYS).



flexwall® W17TT*
wall cabinet arranged for housing
6 Tubs and 6 Trays with the
following dimensions:
W18.3 H3.9 D28.4cm.

* also with holes for front screen (optional)

avant_modularity

Choose furniture
model and
distribution
module.

1




only flexdesk®



Customisation:

The Avant system can be integrated into furniture of any size and placed according to use or room space.


hanging installation with single bottom drawer



modules on plinths



slim modules with hanging installation




2

flexdesk® + flexwall®
single distribution module





3

flexdesk® + flexwall®
storage and distribution together



Avant can be
integrated into longer
compositions to
enhance functionality.





RIGHT HAND
AV. MAN

increase 14% when bent

profile of
heavy winter
gloves A.A.F.



dentalart

Dental Art SpA - I
tel. +39 0445 802000
info@dental-art.it
www.dental-art.it

