

THE „GAME CHANGER“ IN THE FIELD OF PROSTHESIS MANUFACTURING

saremco print - DENTURETEC

CUSTOMIZED FROM THE 3D-PRINTER

- easy, fast and precise production
- always the perfect fit

PERFORMANCE AT IT'S BEST

- very strong and elastic at the same time
- long durability and high wearing comfort

HIGHLY BIOCOMPATIBLE

- no residual monomers proven*
- MMA-free prosthesis material

THE „GAME CHANGER“

Denture base and teeth
customized from the
3D-printer with
DENTURETEC and
CROWNTEC™



*Univ.-Prof. Dr. Dr. Franz-Xaver Reichl, Walther-Straub-Institute for Pharmacology and Toxicology, Ludwig Maximilian University, Munich
**see separate product information

Technical Data

Consistency	medium viscosity, homogeneous
Color	gum colored light
Depth of cure (DIN EN ISO 4049)	$\geq 1,5$ mm
Flexural strength (DIN EN ISO 4049)	> 100 MPa
E-Modulus (DIN EN ISO 4049)	> 1000 MPa
Barcol hardness	> 20
Shelf life	2 years
Medical device	class IIa

Order information (500g bottle)

REF	Product Description
8062	saremco print - DENTURETEC

