



ZARC

ZDown & ZUp

OBTURATION SYSTEM



FROM ODONTOLOGISTS
TO ODONTOLOGISTS

ZDown

Vertical Condensation



TEMPERATURE

Temperature range from 90 °C - 250 °C with 5 memory programs to customize the temperature and cutting time based on the type of material used for obturation.



LONG-LASTING BATTERY

Its long-lasting lithium battery (2600 mAh) allows the device to be used even while charging.



ERGONOMIC

Device designed to start and stop at any angle 360°, and it features a comfortable and light grip.



TEMPERATURE CONTROL

Heats up and cools down rapidly, for more control and safety during treatment.



PLUGGERS

It has interchangeable pluggers of different sizes, for a perfect adaptation to the anatomy of canals.



360° ACTIVATION RING

The device can be activated at any angle, making it easy to use in all situations.

0,5
seconds

It only needs 0.5 seconds to reach 250 °C.

2
seconds

It cools down in 2 seconds to avoid unnecessary risks.

2
hours

It charges completely in 2 hours.

4
hours

It can be used continuously for 4 hours at 250 °C.



Available pluggers



Small (black)
40/0.025



Medium (yellow)
50/0.05



Large (blue)
60/0.06



FROM ODONTOLOGISTS
TO ODONTOLOGISTS



Three-dimensional thermoplastic obturation



WIDE TEMPERATURE RANGE

Wide temperature range from 90 °C to 250 °C. It has 5 memory programs to customize the temperature and cutting time based on the type of material used for obturation.



REUSEABLE TIPS

The tips can be used several times, only changing the gutta-percha pellets. It is the simplest and most economical way to extrude thermoplastic gutta-percha.



ADJUSTABLE TEMPERATURE AND SPEED

It features a screen that allows adjustments to temperature and exclusion speed settings (controlled by the operator).



POWERFUL BATTERY

It has a highly efficient lithium battery (2600 mAh). The gutta-percha is heated quickly thanks to the powerful battery.



CAN BE USED WHILE CHARGING



Constant flow of gutta-percha for a homogenous obturation.




Interchangeable, mouldable and ergonomic tips.



Eliminates hand fatigue thanks to its ergonomic design.

Available tips

Tip	Gauge
Silver Tip 	23G or 25G
Copper Tip 	23G or 25G



ZDown & **Z**Up

THE DEFINITIVE OBTURATION

The most effective and
complete system on the market



Z-Down Reference Codes

Z-Down Obturation System **6380025**

Z-Down + Z-Up Obturation Kit **63800012K**

Plugger small (40/0.025) **6351058**

Plugger medium (50/0.05) **6351059**

Plugger large (60/0.06) **6351060**



Z-Up Reference Codes

Z-Up Obturation System **6380117**

Z-Down + Z-Up Obturation Kit **63800012K**

23G Complete Silver Tip for Z-Up **6351207 + 6351209**

25G Complete Silver Tip for Z-Up **6351207 + 6351208**

23G Multi-torsion Copper Tip for Z-Up (2 units) **6351217**

25G Multi-torsion Copper Tip for Z-Up (2 units) **6351216**

Gutta-percha Pellets for Z-Up **6344003**



FROM ODONTOLOGISTS
TO ODONTOLOGISTS

hola@zarc4endo.com
www.zarc4endo.com

