




Dry-Back[®] Dental Bibs

Stylish
protection



LOSER & CO
öfter mal was Gutes...





Medicom Dry-Back® patient bibs are the perfect addition to the hygiene products in your practice. Dry-Back® bibs consist of two cellulose layers and a polyethylene backing. The three layers adhere securely to each other with no need for additional adhesives.

Thanks to their excellent absorption and specially designed surface and edges, Dry-Back® bibs provide optimal protection against moisture and splashes during treatment.


Dry-Back® bibs are suitable for use during any dental procedure.

Add colour to your practice

Choose from a variety of vibrant hues and create your own colour scheme to coordinate with accessories such as Medicom saliva ejectors, headrest covers, and paper cups.



Adding a splash of colour to your practice creates a cheerful and attractive environment for your patients.



Clothing protection with additional layer of paper and plastic backing

Quick and easy

A practical dispenser allows convenient and hygienic removal of a bib for each patient.

Dry-Back® dental bibs measure 33 x 45.5 cm.



How to order

+49 (0) 21 71 - 70 66 70
@ info@loser.eu



Dry-Back® Dental Bibs



Practical dispenser
Content: 100 pcs.



White	Dark blue	Aqua	Yellow	Lavender	Orange	Dusty Rose
8283-100	8290-100	8289-100	8284-100	8280-100	8291-100	8286-100

Green	Blue	Black	Peach	Lime	Turquoise
8283-100	8282-100	8292-100	8288-100	8293-100	8294-100

The illustrations may show deviations from the original product in terms of colour etc. Our general terms and conditions apply.

