

J-Temp™


TEMPORARY RESIN

U
ULTRADENT
PRODUCTS, INC.


4-IN-1



**BITE RAMPS &
TEMPORARY
OCCLUSAL BUILDUPS**



**STRUCTURE
FOR ISOLATION
CLAMPING**



**TEMPORARY
RESTORATIONS**



**SPLINTING BETWEEN
IMPLANT COPINGS**



MADE BY DENTISTS FOR DENTISTS

"J-Temp temporary resin started because there was always a need for something that didn't exist. The more I used J-Temp temporary resin for different procedures, the more I found that it is really nice to have one product that can cover a gamut of different needs versus needing five or six different materials that may expire before you use them next."

Dr. Jaleena Jessop
VICE PRESIDENT OF CLINICAL AFFAIRS



4897 - J-Temp Syringe Kit

4 x 1,2 ml syringes

20 x Black Mini Tips



J-Temp Temporary Resin

J-Temp temporary resin is premixed for ease of use. It is viscous enough that it will not run, but fluid enough to be self-leveling¹ and easy to manipulate. We recommend use with the Black Mini™ tip for simple, precise placement. When it needs to be removed, the light purple color of J-Temp temporary resin allows for easy identification.

1. Data on file.