

matrix strips and dispensers

Tear resistant, plastic universal strips for use with auto-curing and light-curing composites and cements.
To be used when shaping the restorations for class II, III, IV and V cavities

- Transparent and yellowish-transparent (not recommended for light polymerization of composites)
- Extra thin or standard thickness (retains its stability when being used for larger cavities)
- The yellowish-transparent band can easier be seen on the tooth or cabinet top

Note: the part of the transparent matrix strips with the colour mark must not enter the mouth of the patient.

Enjoy the one time purchase of the patented matrix band dispenser, refill rolls are available at an attractive price

universal strips, straight, length: 120 mm, width: 8.5 mm, package containing 50 pieces

- | | |
|-----------------|--|
| US-120KE | transparent
extra thin: 0.05 mm |
| US-120KN | transparent
standard thickness: 0.08 mm |
| US-120EE | yellowish-transparent
extra thin: 0.05 mm |



universal strips, curved, length: 75 mm, width: 8.5 mm, package containing 50 pieces

- | | |
|----------------|--|
| US-75KE | transparent
extra thin: 0.05 mm |
| US-75KN | transparent
standard thickness: 0.08 mm |
| US-75EE | yellowish-transparent
extra thin: 0.05 mm |



matrix band dispenser incl. matrix band roll, length: 15 m, width: 8.5 mm

- | | |
|-----------------|--|
| US-15PKE | transparent
extra thin: 0.05 mm |
| US-15PKN | transparent
standard thickness: 0.08 mm |
| US-15PEE | yellowish-transparent
extra thin: 0.05 mm |



matrix band refill roll, length: 15 m, width: 8.5 mm, package containing 3 rolls

- | | |
|--------------------|--|
| US-15PKEN | transparent, extra thin: 0.05 mm |
| US-15PKNN | transparent, standard thickness: 0.08 mm |
| US-15PEEN | yellowish-transparent, extra thin: 0.05 mm |
| US-15PKEN E | package with 1 refill roll, transparent, extra thin: 0.05 mm |
| US-15PKNN E | package with 1 refill roll, transparent, standard thickness: 0.08 mm |
| US-15PEEN E | package with 1 refill roll, yellowish-transparent, extra thin: 0.05 mm |

