

PRISMA

Minimalist design combines premium quality with a comfortable round seat. This stool is perfect for shorter procedures and features a height-adjust mechanism under the seat.

Description

Seat height	44-84 cm
Seat dimensions	30 cm
Gas spring options	3
Alternative castors	Yes
Base diameter	50 cm
Adjustable seat	No
Upholstery	Scandinavian leather, vinyl or fabric
Weight	8 kg
Max. load	150 kg
Warranty	5 years

Accessories

Foot ring	Foot rest while working in high chair positions
Foot control	Foot operated height adjustment
SWING arm rest	Arm rest that doubles as a back and abdominal support. Can swivel 360° for complete freedom and versatility



Seat Heights

User height	Gas spring	Seat height	Seat height with SWING
< 160 cm	Standard	44-58 cm	51-60 cm
160-175 cm	High	51-70 cm	58-72 cm
> 175 cm	Extra High	58-84 cm	65-86 cm

Key Features

Mechanism

Made of powder coated steel with a ring for height adjustment.

Seat

Frame is made with wood and moulded foam.

Padding

Upholstered with flexible polyurethane foam, recyclable when using a certain process.

Upholstery

Our standard vinyl contains 100% polyurethane with cotton and polyester backing. It is PVC and Phthalate free. It is also OEKO-TEX certified, which means that the upholstery is easy to clean, and resistant to disinfectants and common cleaning product.

Gas spring

Columns are available in three sizes. The spring is made of polished steel in Germany using nitrogen, a sustainable gas existing in 78% of normal air. Gas springs are maintenance-free and recyclable.

Base

Polished recyclable aluminium.

Castors

Cores are made from polyamide 6 and tread is made of recyclable polyurethane. (Note: larger castor base is required when using SWING support arm).