

ITALIAN QUALITY AND  
INNOVATION FOR  
ORAL HEALTH

**Orabloc®**

NEEDLESTICK SAFETY SYSTEM



**Orabloc®**

**Needlestick Safety  
System**

Quick reference

- 1) Remove protective film from back and place on solid surface. Ensure enough capsules for use. Insert needle vertically into device. After use, follow your standard protocol for needles disposal. Remove device after use. Replace device when empty.

**PIERREL**

Manufactured for:  
Pierrel Pharma srl  
Strada Statale Appia 46/48  
81043 Capua (Ce) - Italy

**Insert here**

Tel: 877-394-2465  
orabloc@pierrelgroup.com  
www.orabloc.com

Designed in the USA  
U.S. Pat. No. 8,083,098

Lot: US012020



Consult operating  
instructions



Recommended operating  
and Storage temperature  
+5°C to +40°C (41°F to 104°F)

Rx Only

Do not dispose of this product in the unsorted municipal  
waste stream. Dispose of this product according  
to local regulations.



# Orabloc®

## NEEDLESTICK SAFETY SYSTEM

### STOP NEEDLESTICK INJURIES

No worries, be Safe, be Smart.



**Over 3 million needlestick injuries happen every year, according to World Health Organization, making it essential to safeguard healthcare workers and meet regulatory requirements.**

Orabloc Needlestick Safety System is a Class I medical device that complies with key safety regulations, including:

#### USA 2000 - Needlestick Safety and Prevention Act

The first legislation specifically aimed at preventing needlestick injuries.  
(Referencing OSHA Bloodborne Pathogens Standards)

#### EU 2010 - Directive n. 2010/32/UE

Focused on the prevention of injuries caused by sharp medical instruments, including needles.  
It mandates the use of safety-engineered devices in all healthcare settings — both public and private — to prevent needlestick injuries.

#### TECHNICAL FEATURES

- › 100 rigid silicone caps for 100 safe injections
- › Prevents needle penetration
- › Compatible with all needle types — even bent or angled
- › Color-coded caps for clear visual tracking
- › Non-slip silicone base for surface stability
- › Dual-end protection: covers both tips of the needle

#### OPERATIONAL ADVANTAGES

- › Portable: ideal for use at multiple workstations
- › Single-handed recapping: simplifies handling and enhances safety
- › No sterilization required
- › Emission-free and energy-efficient: no power source needed



WWW.PIERRELGROUP.COM

U.S. Appia 7 bis 46/48 - 81043 Capua (CE) - Italy | ☎ +39 0823626111 | ✉ info@pierrelgroup.com

CE Material intended exclusively for use by medical and healthcare professionals. - Rev. 05/2025