

Perfect for drying your handpieces


Q-Dry provides an optimal internal drying solution for handpieces after manual and automatic cleaning processes. By effectively removing liquids from within the handpieces, Q-Dry ensures that no spots are left on the pouches after sterilization.

Ergonomic design
Enjoy an intuitive user experience and the ability to dry other hollow instruments.

Quick Connect by W&H
With the Quick Connect feature, you can switch effortlessly between different types of handpieces.



Compatibility
Benefit from using the same adaptors found in the Assistina portfolio.



Q-Dry
with Quick Connect by W&H

System overview adaptors Assistina Twin & Assistina One

Straight and contra-angle handpieces

All straight and contra-angle handpieces with ISO connection	
Contra-angles with removable heads	Adaptor REF 02679000
T1 Classic from Sirona Dental Systems for contra-angle handpieces	Adaptor REF 02691000
FG chuck system (W&H, Sirona Dental Systems)	Adaptor REF 02693000
W&H handpiece and contra-angle head	Adaptor REF 05204600

Adaptor Quick ISO
REF 06422600



Turbines

W&H turbines, air scaler and air motors with Roto Quick system	W&H Adaptor REF 02690400
T1, T2, T3 from Sirona Dental Systems	Adaptor REF 02692000
Turbines with Multiflex connection from Kaltenbach & Voigt	Adaptor REF 04713200
Turbines with Borden 2(3)-hole connection	Adaptor REF 02083500
Turbines or adaptors with standard 4-hole connection (with or without light)	

Adaptor Quick RM
REF 06395100

