

#whdentalwerk



wh.com



**Made in
AUSTRIA**
and used all
over the
world! ✓

PROMOTIONS

FOR BELGIUM & LUXEMBOURG

Prijzen geldig vanaf 16 september tot en met 31 december 2025

#whdentalwerk



wh.com



**Made in
AUSTRIA**
and used all
over the
world! ✓

PROMOTIONS

FOR BELGIUM & LUXEMBOURG

Prijzen geldig vanaf 16 september tot en met 31 december 2025

Synea Fusion

The Synea Fusion is the world's most widely used contra-angle handpiece – quiet, strong and durable.

The air turbines in this series are equipped with ceramic bearings, a quadruple spray system, and are available with various head sizes.

Safety first: All Synea turbine handpieces are equipped with a technology that ensures that aerosols do not enter the instrument.

RESTORATION & PROSTHETICS



The powerhouse of contra-angle handpieces!

Synea Fusion Power Edition

WG-900 LT

Red contra-angle handpiece 1:4 with light, Ø 10.7 mm triple spray

36*-month warranty

REF 30459000 A50376

Promo
749,00



synea FUSION



Package Deal PD0326 Tot 57% korting

Koop **3 x Synea Fusion instrumenten** naar keuze aan **promoprijs** (mengen mag) en koop de **WG-56 LT** aan **1 euro** (waarde: 825,00 euro).

Restoration & Prosthetics

WG-56 LT

Blue contra-angle handpiece 1:1 with light, Ø 9.5 mm single spray

36*-month warranty

REF 30009000 A37655

Promo
469,00



WG-66 LT

Green contra-angle handpiece 2:1 with light, Ø 9.5 mm single spray

36*-month warranty

REF 30011000 A37654

Promo
659,00



WG-99 LT

Red contra-angle handpiece 1:5 with light, Ø 9.5 mm Quattro spray

36*-month warranty

REF 30007000 A37653

Promo
635,00



Good spray and lighting

With a narrow head

TG-97 L* for Roto Quick

TG-97 LM** for Multiflex®

Air turbine *LED+

**with light, 18 W, Ø 10 mm Quattro spray

24*-month warranty

*REF 30004000 A37657-L

**REF 30006000 A37657-LM V37657

Promo
559,00



TG-98 L* for Roto Quick

TG-98 LM** for Multiflex®

Air turbine *LED+

**with light, 21 W, Ø 11.5 mm Quattro spray

24*-month warranty

*REF 30001000 A37658-L

**REF 30003000 A37658-LM V37658

Promo
559,00



HG-43 A

Blue handpiece 1:1 no light single spray

24*-month warranty

REF 30013000 A37656

Promo
405,00



Cast from a single block, no seams



Small head means excellent view of working area


Synea Vision

High-end contra-angle handpieces and air turbines; lightweight, quiet, with excellent grip, a slim head and a long warranty. The air turbines in this series are equipped with ceramic bearings, a quintuple spray system, a scratch-resistant surface and are available with four different head sizes.

synea VISION

Ring LED+

Work without shadows!



• WK-56 LT

Blue contra-angle handpiece 1:1 with light, Ø 9.5 mm single spray

36-month warranty

REF 30025000 A37580

Promo
679,00




• WK-66 LT

Green contra-angle handpiece 2:1 with light, Ø 9.5 mm single spray

36-month warranty

REF 30026000 A37596

Promo
739,00



Synea Vision Power Edition

• WK-900 LT

Red contra-angle handpiece 1:4 with light, Ø 10.7 mm triple spray

36-month warranty

REF 30418000 A50377

Promo
949,00



With a narrow head

• WK-93 LT

Red contra-angle handpiece 1:4.5 with light, Ø 9.1 mm quintuple spray ceramic bearings

36-month warranty

REF 30024000 A37549

Promo
879,00



• WK-99 LT

Red contra-angle handpiece 1:5 with light, Ø 9.5 mm quintuple spray ceramic bearings

36-month warranty

REF 30023000 A37652

Promo
899,00



- Extremely powerful and robust contra-angle handpiece
- For the removal of zirconium and attachments
- Optimal power to speed ratio thanks to 1:4 speed ratio
- Indispensable for difficult restorations
- *25% longer lifespan

Package Deal PD0419

Tot 54% korting

Koop **3 x Synea Vision instrumenten** naar keuze aan **promoprijs** (mengen mag) en koop de **WK-56 LT** aan **1 euro** (waarde: 1066,00 euro)

Very narrow head

TK-94 L* for Roto Quick
TK-94 LM** for Multiflex®

Air turbine *Ring LED+
**with light
14 W, Ø 9 mm
quintuple spray

24-month warranty

*REF 30016000 A38968
**REF 30031000 A38972

Promo
859,00



TK-97 L* for Roto Quick
TK-97 LM** for Multiflex®

Air turbine *Ring LED+
**with light
21 W, Ø 10 mm
quintuple spray

24-month warranty

*REF 30017000 A38969
**REF 30018000 A38973

Promo
859,00



TK-98 L* for Roto Quick
TK-98 LM** for Multiflex®

Air turbine *Ring LED+
**with light
24 W, Ø 11.5 mm
quintuple spray

24-month warranty

*REF 30019000 A38970
**REF 30020000 A38974

Promo
859,00



TK-100 L* for Roto Quick
TK-100 LM** for Multiflex®

Air turbine *Ring LED+
**with light
26 W, Ø 13 mm
quintuple spray

24-month warranty

*REF 30021000 A38971
**REF 30022000 A38975

Promo
859,00



Synea Vision Short series; 20% lighter and shorter, yet just as powerful as standard Vision instruments!

WK-56 LT S

Blue contra-angle handpiece 1:1 with light, Ø 9.5 mm single spray

36-month warranty

REF 30217000 A39134

Promo
679,00



WK-66 LT S

Green contra-angle handpiece 2:1 with light, Ø 9.5 mm single spray

36-month warranty

REF 30218000 A39102

Promo
739,00



WK-99 LT S

Red contra-angle handpiece 1:5 with light, Ø 9.5 mm quintuple spray

36-month warranty

REF 30215000 A50661

Promo
899,00



WK-93 LT S

Red contra-angle handpiece 1:4.5 with light, Ø 9.1 mm quintuple spray

36-month warranty

REF 30216000 A39175

Promo
879,00



Lightweight and compact!

Implantmed Plus SI-1023

The powerful **Implantmed Plus** from W&H guarantees a high level of safety during treatment for both you and your patients. With its intuitive user interface, colour touchscreen and glass surface, the Implantmed offers a truly user-friendly experience. The drive is ideal for surgical procedures with its automatic thread-cutting function, automatic torque control and adjustable parameters. For **piezoelectric surgery**, the **Piezomed Plus Module** can be fully integrated with the Implantmed Plus surgical motor, as shown in the photo.

Both devices can be operated with a single foot control (which can be wireless), and the instruments connected can be supplied with the same physiological saline solution.

implantmed



Work according to protocol



with automatic thread-cutting function



Piezomed Plus Module

Implantmed Plus

Piezomed Plus module

Piezomed brings all the benefits of modern ultrasonic technology within reach: High-frequency microvibrations facilitate incredibly precise bone surgery. What is more, thanks to the cavitation effect you can maintain a virtually blood-free surgical site.

- Sinus lift
- Bone harvesting
- Apical resection
- Retrograde endodontics
- Extraction
- Implant site preparation
- Periodontal (sterile!)
- Etc.

piezomed



Course on the piezoelectric technique in surgery

Renowned dental implantologist Irfan Abas will teach you the numerous possibilities of the piezoelectric technique in modern surgery. In collaboration with Megagen and W&H; see www.abas-academy.com for more information.

ABAS
implant academy



Visit video.wh.com for all applications

Piezomed Plus Module SA-435 M

The Piezomed Plus Module can be connected to the Implantmed Plus surgical motor via plug and play. The most common oral surgical treatments can be performed with this combination.

Features:

- › tip recognition
 - › LED light
 - › handpiece with a 1.8-metre cable
- REF 90000350 A47782

Promo
4375,00



Data:
21 November
2025
29 May
2026



PACKAGES

Package A PD0056 Implantmed Plus SI1023 + WS-75 L

- › Implantmed Plus **with light**
- › Wireless foot control
- › Surgical contra-angle handpiece WS-75 L (20:1, green)

REF 90000166 A40613
REF 30033000 A40683

Package Deal
4375,00

Package B PD0286 Implantmed Plus SI1023 + WS-75

- › Implantmed Plus **no light**
- › Wired foot control
- › Surgical contra-angle handpiece WS-75 (20:1, green)

REF 90000163 A49287
REF 30065000 A01881

Package Deal
3489,00

Package C PD0055 Implantmed Classic SI923 + WI-75 E/KM

- › Implantmed Classic **no light**
- › Wired foot control
- › Surgical contra-angle handpiece WI-75 E/KM (20:1, green)

REF 90000161 A44299
REF 10207550 A28895

Package Deal
3189,00

Package D PD0418 Implantmed Plus SI1023 + WS-75 L + S-11 L

- › Implantmed Plus **with light**
- › Wireless foot control
- › Surgical contra-angle handpiece WS-75 L (20:1 green)
- › Surgical handpiece S-11 L (1:1, blue)

REF 90000166 A40613
REF 30033000 A40683
REF 30058000 A41366

Package Deal
5093,50



WS-75 L | WS-75 | WI-75 E/KM | S-11 L



**Piezomed
Module Plus
SA-435 M**

Flight case

- › For the Implantmed Plus, Piezomed Module and accessories (in combination with the equipment)

A51720
REF 07962790

Promo
399,00



Implantmed Plus SI-1023



Implantmed Classic SI-923



W&H Surgical System

**Implantmed Plus & Piezomed
Plus module, a single system
for all your oral surgery**

**Compact
Complete
Compatible**

Package E PD0175

- › **Implantmed Plus SI-1023** (top model); with LED light, wireless foot control, documentation function and a 1.8-metre cable
- › **Piezomed Plus Module SA-435 M**; with LED light, tip recognition, high power and a 1.8-metre cable
- › **WS-75 L** surgical Green contra-angle handpiece (20:1)

REF 90000166 A40613
REF 90000350 A47783
REF 30033000 A40683

Package Deal
8500,00

Surgical trolley *

- › With drawer for instruments and ISQ equipment, etc.
- › Features large wheels for easy rolling
- › Adjustable height for the drawer

REF 8088730 A48820

Prijs
1399,00

*Only in combination with one of the packages



Easy integration



Simple replacement of irrigator elements



Wireless or wired foot control

Surgical handpieces & contra-angle handpieces

The instruments below boast lights when connected to an EM19 LC motor.
When connected to an EM19 motor there is no light source.

• WS-75 L

With light.
Implant site preparation
and implant placing

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 30033000 A40683



Promo
1049,00

• WS-56 L

With light.
Osteotomy, sequestrotomy,
germectomy, hemisection

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30040000 A40697



Promo
1049,00

• WS-91 L

With light.
Specially designed to
improve access to the
premolar and molar area

1:2.7, contra-angle handpiece
for surgical drills (red)
FG, with a 45° head tilt

REF 30049000 A40698



Promo
1049,00

• WS-75 LG

Boasts its own lighting
Light output depends on speed.
Applications such as WS-75 L

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 30032000 A34809



Promo
1090,00



Can be
disassembled
without tools
for optimal
hygiene

• WS-92 L

With light.
Osteotomy, sequestrotomy,
germectomy, hemisection

1:2.7, contra-angle handpiece
for surgical drills (red)
FG

REF 30043000 A41640



Promo
1049,00

• S-11 L

With light.
Apical surgery, osteotomy,
bone modelling and
wisdom tooth surgery

1:1, handpiece for
surgical drills (blue)
Ø 2.35 mm

REF 30058000 A41366



Promo
1049,00

• S-11 LG

Boasts its own lighting
Light output depends on speed.
Applications such as S-11 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30057000 A39292



Promo
1090,00

With mini-LED+
and optimal cooling



See wh.com for more models and versions

• WS-56 LG

Boasts its own lighting.
Light output depends on speed.
Application as WS-56 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30039000 A51712



Promo
999,00

• WS-91 LG

Boasts its own lighting.
Light output depends on speed.
Application as WS-91 L

1:2.7, contra-angle handpiece
for surgical drills (red)
FG, with a 45° head tilt

REF 30048000 A39226



Promo
1029,00

• WS-92 LG

Boasts its own lighting.
Light output depends on speed.
Application as WS-92 L

1:2.7, contra-angle handpiece
for surgical drills (red)
FG

REF 30042000 A48939



Promo
1049,00

• WS-56

No light.
Applications such as WS-56 L

1:1, contra-angle handpiece
for surgical drills (blue), Ø 2.35 mm

REF 30038000 A28896



Promo
799,00

• WS-91

No light.
Applications such as WS-91 L

1:2.7, contra-angle handpiece
for surgical drills (red)
FG, with a 45° head tilt

REF 30047000 A49284



Promo
899,00

• S-11

No light.
Applications such as
S-11 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30056000 A01882



Promo
749,00

• WS-75

No light.
Applications such as
WS-75 L

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 30065000 A01881



Promo
899,00

• WI-75 E/KM

No light.
Applications such as WS-75 L
Cannot be disassembled.

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 10207550 A28895



Promo
719,00

• WS-92

No light.
Applications such as WS-92 L

1:2.7, contra-angle handpiece
for surgical drills (red), FG

REF 30041000 A31341



Promo
899,00

See wh.com for more models and versions



Instrument Kit Bone

Surgery Tips included:
B1, B2R, B2L, B3, B4, B5

REF 07173100
A40330

Prijs
1062,00



Instrument Kit Sinus Special

For sinus lifts.
Tips included:
S1, S2, S3, S4, S5

REF 07173200
A40331

Prijs
746,00



Instrument Kit Endodontics

Tips included:
R1D, R2RD, R2LD, R3D,
R4RD, R4LD

REF 07173300
A51716

Prijs
833,00



Instrument Kit Periodontal

Tips included:
P1, P2RD, P2LD

REF 07173400
A51717

Prijs
468,00



Instrument Kit Extraction

Tips included:
EX1, EX2

REF 07173500
A40339

Prijs
328,00



Instrument Kit Implant Crestal P

Posterior. Tips included:
I1, I2P, I3P, I4P,
Z25P, Z35P

REF 07997110
A51718

Prijs
1456,00



Instrument Kit Implant Crestal A

Anterior. Tips included:
I1, I2A, I3A, I4A,
Z25P, Z35P

REF 07997090
A51719

Prijs
1456,00



See wh.com for the full range of piezoelectric tips



Fast
Precise
Powerful

piezomed

Instrument B6 & B7, surgical saws for, among other things, apical resection and the removal of root remnants or alveolar bone

B6 REF 07022500 A40340

Prijs
215,00

B7 REF 07022600 A40341

Prijs
215,00



Piezoelectric handpiece

featuring a light,
tip recognition
and a spray.

• SZ-75

Special handpiece
for placing zygoma
implants

20:1, contra-angle handpiece
for zygoma drills (green)
Ø 2.35 mm

REF 30095000 A49285

Promo
1499,00



• S-9

Slightly angled; same
applications as S-11 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm, max. length 45 mm

REF 30053000 A01883

Promo
1049,00

also available with
a generator, as



• S-9 LG

REF 30054000 A51711

Promo
1159,00

• S-15

1:1, contra-angle handpiece
for surgical drills (blue),
for hard-to-reach areas
Ø 2.35 mm, length 45 mm

REF 30060000 A51713

Promo
849,00



• S-10

1:1, contra-angle handpiece
for surgical drills (blue),
for hard-to-reach areas
Ø 2.35 mm, length 70 mm

REF 30059000 A37216

Promo
1069,00



• S-12

1:1, contra-angle handpiece
for surgical drills (red),
for hard-to-reach areas
Ø 2.35 mm, length 70 mm

REF 30061000 A51714

Promo
1069,00



• S-16

1:1, contra-angle handpiece
for surgical drills (red),
for hard-to-reach areas
Ø 2.35 mm, length 45 mm

REF 30411000 A51715

Promo
1194,00



Surgical saw handpieces

Fast, powerful and durable

The saw handpieces are specially designed for conversion osteotomies and bone harvesting. A consistently high cutting performance is guaranteed even after many uses at maximum speed.

The saw blades are available in different sizes and are supplied sterile.

S-8 R

Saw handpieces with reciprocating motion

REF 10100801
A33125

Promo
1299,00



S-8 O

Saw handpieces with oscillating motion

REF 10100802
A29430

Promo
1549,00



S-8 S

Saw handpieces with sagittal motion

REF 10100803
A49283

Promo
1499,00




Easy to disassemble
and clean really effectively

Osstell Beacon

Implant Stability Quotient (ISQ)

- › Monitoring of osseointegration
- › Makes an earlier prosthetic process possible
- › Improved predictability of treatment for patients with risk factors
- › Easy communication with patients through visualisation of the growing-in phase
- › Maximum flexibility thanks to wireless technology
- › The Osstell Beacon can be integrated with the latest generation of W&H Implantmed



Work in a validated and documented way; ideal for referral practices

Osstell Beacon

REF 103000 A43554

Promo
2129,00



NEW!

IMPLANTOLOGY



Osstell Classic

REF 08086148 A51729

Promo
1659,00

The Osstell Classic has the same functionalities as the Osstell Beacon, but without the option to use the Osstell Connect Portal, and it cannot be connected to W&H Implantmed.

Includes free use of the Osstell Connect portal:

- › build up the implant history
- › compare implants
- › compare patients
- › compare pathways
- › compare problems

SmartPegs are necessary in order to measure osseointegration. These are available from your dealer, see wh.com for an overview.



V46725

Proxeo Twist

Top quality handpieces and contra-angle handpieces for prophylaxis

For your treatment unit

proxeo ^{TWIST}

WP-66 W

Prophy 4:1 contra-angle handpiece
for W&H LatchShort prophyl cups

REF 30189000 A42771

Promo
325,00



WP-64 MU

Prophy 4:1 contra-angle handpiece
with snap-on adapter and key
for all standard screw-in
prophy cups and brushes

REF 10216443 A28731

Promo
325,00



More space
with the
LatchShort
attachment



WP-64 M

Prophy 4:1 contra-angle handpiece
for System Young Screw-in
prophy cups and brushes.
With key.

REF 10216403 A02406

Promo
325,00



WP-66 M

Prophy 4:1 contra-angle handpiece
for all prophy cups, brushes
and snap-on adaptors with a
2.35 mm shaft

REF 30190000 A39034

Promo
325,00



HP-44 M

Prophy 4:1 handpiece
for Doriot system
disposable prophy cups

REF 10114403 A02320

Promo
279,00



Proxeo Twist

Cordless Polishing System



Stand-alone polishing contra-angle handpiece

- Work flexibly without cables
- Easily regulate speed with wireless foot control
- Perfect control of speed for precision cleaning and efficient polishing
- Three different attachments
- Minimal use of aerosols thanks to rotary polishing

Handy for working on location!

PL-40 HW + PL-66 W

Proxeo Twist Cordless
+ LatchShort cups

REF 90000346 A46570

Promo
999,00



PL-40 HW + PL-64

Proxeo Twist Cordless
+ screw-in for
YOUNG cups

REF 90000347 A51709

Promo
999,00



PL-40 HW + PL-44

Proxeo Twist Cordless
+ Doriot attachment for
disposable contra-angle
handpieces with cup

REF 90000348
A51710

Promo
999,00



LatchShort
cups

Screw-in for
YOUNG cups

Doriot system
with cup



Wireless
foot control

NEW!

Lisa Mini

Compact, fast and easy to operate!

Just like the renowned Lisa large autoclave, this smaller version B autoclave is fast and features an easy-to-use touchscreen. Available with a 5- or 3-litre chamber.

Key features:

- › Fast contra-angle handpiece cycle in just 10 minutes (5-litre chamber)
- › Standard cycle time of 22 minutes (5-litre chamber)
- › Data can be logged via USB, LAN or the W&H Steri App
- › Can be connected to a demineralised water supply
- › A label printer can also be connected
- › Available in 2 or 3-tray format (5-litre chamber)

lisa MINI

**5 litres**

REF 19935100 A51721

Promo
4879,00**5 litres***

REF 19945100 A50654

Promo
4429,00**3 litres***

REF 19944100 A50653

Promo
3999,00

*requires compressed air connection



Handy as an extra autoclave, when you have limited space or for small volumes!

Other items can be placed on top of the unit to make the most of your sterilization space. Comes with 2 trays as standard.

Lisa

High-end type B autoclave with extremely short cycle times.

Standard cycle time of 30 minutes, fast contra-angle handpiece cycle in 13 minutes. Features Eco Dry+ technology that adjusts the drying time to the weight of the load. Boasts an LED display and a touchscreen. Available with a 17- or 22-litre chamber. Data can be logged via USB, LAN or the W&H Steri app. The sterilization process can be tracked and documented down to the individual instrument or instrument kit via EliTrace. No extra software or computer needed. Thanks to smart technology, the Lisa warns you when your load is overheated, even if the door is open – no more burnt fingers!

Fast and quiet!

17 litres

REF 19931180 A42621

Promo
5990,00**22 litres**

REF 19932180 A42687

Promo
6490,00

The Lisa from W&H is the most widely used autoclave in the world, and for good reason!



lisa

Lara, Lara XL and Lara XXL

The fastest 28-litre class B autoclave on the market!

- › Tried and trusted technology
- › Large capacity
- › Various upgrades available
- › Easy to operate

The Lara is also available in 17- and 22-litre versions:

17 litres

REF 19941120
A42622

Promo
5490,00

22 litres

REF 19942120
A42686

Promo
5990,00

Lara XXL 38 litres

REF 19936110
A51724

Promo
8960,00

Lara XL 28 litres

REF 19933110
A48819

Promo
6990,00

**38-litre capacity for
up to 12 kg per cycle**



**Fits
60 cm trays**

**Highly suitable
for surgical cassettes
thanks to its large size**



lara^{XL}
lara^{XXL}

Lina

**Type B autoclave with best-in-class
cycle times.**

- › Optimal sterilization programmes
- › The Eco B sterilization programme processes 0.5 kg in 25 minutes
- › Lina can reprocess loads of up to 6 kg to satisfy all the requirements of modern dental clinics
- › The B Universal 121°C cycle even allows you to reprocess delicate objects, including porous objects such as surgical clothing



17 litres

REF 19951131
A51722

Promo
4799,00

22 litres

REF 19952131
A51723

Promo
4990,00

**Excellent
value for
money!**



lina



The **Activation Code** system allows you to activate additional features to meet your practice's changing needs and future requirements.

Helix test kit

Complete set; consisting of the holder (the "snail shell") and 250 test strips.

Each test strip has a QR code that can be scanned and saved in our Steri app or in your network. Each code is linked to a sterilization cycle, which validates the process.

With QR-coded strips

REF T800207X

A34926

Prijs
159,00



Q-dry

Perfect solution for internal drying of contra-angle handpieces and air turbines.

Using the Assistina adaptors and air pressure, you can easily dry instruments prior to sterilization.

See wh.com for how-to videos.

q-dry

REF 19723118

A51689

Prijs
179,00



Multidem C27 CARTRIDGE

A36570

Promo
199,00



Multidem C27

Features a microfilter, capable of supplying two autoclaves with demineralised water simultaneously

Demineralised water supply for your autoclave

REF 19723112 A35531

Prijs
363,00



Water distiller

4-litre capacity, enough for 8 autoclave cycles. Ready in 5 hours.

REF 19723101 A32028

Prijs
359,00



NEW!

U-Son Ultrasonic cleaner

U-Son 6 litres

REF 19600005 A51687 + deksel + mandje

Promo
999,00

U-Son 3 litres

REF 19600004 A51686 + deksel + mandje

Promo
659,00



LISA Safe label printer

Can be connected directly to the W&H autoclave, no external computer required. The labels are only printed if the autoclave cycle has been completed successfully. Depending on your autoclave type, a connection kit may be required. One possible option is the barcode scanner.

Laser printer for steri bags.

REF 19721109 A46215

Promo
669,00



Seal²

Excellent basic sealer suitable for 7 roll sizes, from 5.5 up to 30 cm. No settings required, features transparent roller guide. The roll holder can optionally be wall-mounted.

REF 19722101 A33878

Promo
644,00



BePro soaking bath

including draining insert. 3-litre capacity.

REF 19500300 A51688

Prijs
56,00




SealVal²

This machine features a touchscreen through which settings can be easily saved or adjusted. The roll holder can also be wall-mounted to save space, and all data can be stored via a LAN connection. Suitable for 7 roll sizes, from 5.5 up to 30 cm.

Fast thermal sealer for rolls up to 30 cm wide

REF 19722103 A51685

Promo
1222,00



Assistina Twin

Reprocess up to 360 instruments per hour!

- › Automatic lubrication and maintenance device
- › Reprocess instruments in just 10 seconds
- › Short waiting times thanks to the two-chamber system
- › Simple and ergonomic operating concept
- › The Assistina Twin Care Set can handle the reprocessing of no fewer than 2,800 instruments
- › Add and remove instruments even more easily and faster thanks to the Quick Connect adaptor solution

Supplied with
TWIN Care Set:
HEPA filter, cleaning fluid and oil

REF 07484000 A41745
Prijs: 120,00

assistina^{TWIN}

NEW!
Now with optional
blow-through
function, for
bone-dry
instruments!

Assistina Twin

Assistina TWIN MB-302
+ 1 Quick ISO adaptor
+ 1 Quick RM adaptor*

REF 30310003 A49286

Promo
2274,00

Assistina TWIN MB-302
+ 2x Quick RM adaptors*

REF 30310002 A51728

Promo
2274,00

Assistina TWIN MB-302
+ 2x Quick ISO adaptors

REF 30310001 A51727

Promo
2274,00



*The RM adaptor requires an additional connecting piece for the air turbines:

for air turbines with
a W&H Roto Quick
connection –

REF 02690400

A01920

Promo
179,00



for air turbines with
a Multiflex (KAVO)
connection –

REF 04713200

A28198

Promo
99,00



for air turbines with
a Sirona T1, T2, T3
connection –

REF 02692000

A01844

Promo
189,00



Assistina ONE

The new 1-door
reprocessing device from W&H.

Superfast reprocessing time: The Assistina ONE can lubricate and clean a contra-angle handpiece or air turbine in just 18 seconds.

ONE Care Set:
HEPA filter, cleaning fluid and oil

REF 08064970 A48811
Prijs: 114,00

assistina^{ONE}

2x the speed
of the previous
model!

Assistina ONE

Assistina ONE MB-301
+ 1 Quick ISO adaptor
+ 1 Quick RM adaptor*

REF 30398003

A48810

Promo
1285,00

Assistina ONE MB-301
+ 1x Quick RM adaptor*

REF 30398002

A51726

Promo
1250,00

Assistina ONE MB-301
+ 1x Quick ISO adaptor

REF 30398001

A51725

Promo
1250,00



NEW!

The new Thed and Thed+

Top-class thermal disinfectors of the highest quality

- › With active drying, your precious instruments come out of the machine dry
- › Fast cycles
- › Compact, large-capacity machine
- › Offers single or double load level
- › The chemicals stay in the machine
- › Numerous accessories available for all instruments
- › Easy operation and connectivity
- › Available as a tabletop or undercounter model
- › Features quick cycle for contra-angle handpieces
- › Large, practical display


Thed and Thed+ tabletop model, with one rack



Thed and Thed+ tabletop model, with two racks



Space for 15 contra-angle handpieces



Thed+ also features:

- + light bar with status indication
- + programme advice
- + recognition of the racks and any injector rail in place
- + internal lighting
- + automatic adjustment of the following for each load:
 - cycle duration
 - water volume
 - chemical usage
 - power consumption

Thed and Thed+ undercounter model, with one rack



Thed and Thed+ undercounter model, with two racks



Prices*:

Thed tabletop model

REF 21050000 A51690

Promo
7495,00

Thed undercounter model

REF 21051000 A51691

Promo
8095,00

Thed+ tabletop model

REF 21055000 A51692

Promo
8995,00

Thed+ undercounter model

REF 21056001 A51693

Promo
9995,00

*excluding accessories and installation

Installation, servicing and maintenance by your dealer or certified W&H service partners. Request the brochure for specifications and all accessories!

thed

The tabletop models feature handy cassettes for cleaning agents



The undercounter models feature a drawer for cleaning agents

