

PROMOTIONS

Up to
-40%

Valid october - december 2025

ELSODENT in a few words



30 years of innovation and expertise in dental materials

A range of premium dental products, made in France.

Outstanding clinical qualities, validated by tests carried out in renowned universities.

Innovative formulations without endocrine disruptors or controversial degradable resins. BIO+ range without BPA, BisGMA, BisEMA, TEGDMA or added HEMA.

#innovation

#biocompatible

#madeinfrance

#performance



PURESHADE

Universal mimetic composite



V50474
Sprit 4 g

Prijs: 53,99

38% korting

**Promo
42,90**

**Per 3 / stuk
34,90**

- Simplified shade selection
3 shades = 16 VITA
- Perfect mimetics
- Antero-posterior



V50475
20 x 0,2 g

Prijs: 74,68


39% korting

**Promo
59,50**

**Per 3 / stuk
47,90**

PURESHADE	Shade 1	Shade 2	Shade 3
Equivalence VITA	A1, B1, B2, C1, D2	A2, A3, C2, D3, D4	A3.5, A4, B3, B4, C3, C4
Ref. 4g syringe	PURSH41	PURSH42	PURSH43
Ref. 20x0,2g unit dose	PURSH201	PURSH202	PURSH203

MD class IIa (0197)



Restauration classe I avec Pureshade 'Shade 2'
(M.SCHOOL -Paris)

PUREFILL Flow

Nanohybrid flowable composite



V41304
2 x 2 g

Prijs: 66,66

24% korting

**Promo
53,90**

**Per 2 / stuk
49,90**

Kit contains 2 x 2g syringes + 10 tips

PURF-4-* PUREFILL Flow - 2x 2g + 10 tips

PUREFILL Flow

high viscosity

*shade A1/B1 - A2 - A3 - A3,5 - B2 | MD class IIa (0197)

Price in € excl. VAT

BONDING

HEALBOND Universal



8th generation universal adhesive

A51405
5 ml

Prijs: 100,22

44% korting

Promo

59,90



PUBLI

MD class IIa (0197)



**New
design**



250 drops

- ▶ All types of etching
- ▶ All surfaces
- ▶ Contains 10-MDP
- ▶ Direct & Indirect restorations

CEMENT

INNOCEM



Adhesive, self-etching, dual-cure resin-cement

V42225
5ml syringe + 10 mixing tips + 10 regular intraoral tips

Prijs: 109,41

39% korting

Promo

82,90

Per 2 / stuk

69,90



Translucent A2	Opaque A2	Opaque A3
IC-10-TA2	IC-10-OA2	IC-10-OA3

Kit contains = 1 automix syringe 5ml + 10 tips
(mixing + intra-oral) MD class IIa (0197)

SOLO

Simple, reproducible protocol

- ▶ Adhesion to all surfaces
- ▶ Aesthetics | Total masking power for opaque shades.

DUO

For enhanced adhesion with HEALBOND Universal bonding

- ▶ **SMART CONCEPT**
HEALBOND Universal becomes dual on contact with INNOCEM, without activator
- ▶ **IPC technology**
"Touch cure" effect, for an instantly watertight joint

TEMPORIZATION

SEALTEMP & SEALTEMP S

100% resin self-cure temporary cements



Prijs: 81,37

39% korting

**Promo
65,00**

**Per 2 / stuk
57,90**

25 ans de
succès clinique !

Kit contains 2 syringes automix 5ml + 20 tips

STD-16	SEALTEMP	Strong retention	A38312
STDs-16	SEALTEMP S	Standard retention	A38268

MD class IIa



TEMPORARY CROWNS

PHOENIX

Resin for temporary crowns & bridges

V40301

Prijs: 99,76

32% korting

**Promo
82,50**

**Per 2 / stuk
69,90**

Universal, resistant &
aesthetic



Kit contains 1 x 50ml cartridge + 10 tips

SHADE	A1	A2	A3
PHOENIX	PH-50-A1	PH-50-A2	PH-50-A3

MD class IIa

CORE BUILD-UP

CORE D



Core build-up "dual" composite

A38249-D

Prijs: 179,39

25% korting

**Promo
135,90**

Kit contains 1 x 25ml cartridge + 20 tips

Ref. CDD-25 teinte dentine MD class IIa (0197)



▶ **SMART CONCEPT** The bonding becomes dual at contact with the resin, without activator.



without matrix

LIQUID DAM

NOTREDAM

Light-cure liquid dam



A50411

Prijs: 96,08

35% korting

Promo

64,90



Perfect for

- ▶ Endodontics
- ▶ Bleaching
- ▶ Prosthodontics
- ▶ Aesthetic & Restorative Dentistry

Ref. ND-10 2 syringes 3ml + 10 tips | MD class IIa (0197)

ENDODONTICS

ENDONEEDLE

Endocanal irrigation cannula

Side
opening



Box of 100,

V38275

G23 bleu	ED23-100	Prijs: 70,80	Promo: 59,90 Per 2: 52,90 stuk
G27 jaune	ED27-100	Prijs: 70,80	Promo: 59,90 Per 2: 52,90 stuk
G30 mauve	ED30-100	Prijs: 84,59	Promo: 69,90 Per 2: 59,90 stuk

MD class I



ETCHING

G-ETCH

Etching gel, phosphoric acid 38%

Prijs: 76,22

31% korting

Promo

52,00



- ▶ Ne coule pas
- ▶ Rinçage facile & rapide

Ref. JKB-50 JUMBO KIT 1 syringe 50ml + 5x empty syringes 3ml + 40 tips A07870
ETG-16 MAXI KIT 4 syringes 3ml + 20 tips | ~~43,69€~~ -30% 30,90€ A46371

MD class IIa

ULTIMA

MULTI-MODE & MULTI-INDICATION

A46331

Prijs: 1490,08

39% korting

**Promo
899,00**

LIGHT-CURING

SOFT
STANDARD
BOOST
SEQUENTIAL BOOST

DIAGNOSIS

UV DIAGNOSIS
TRANSILLUMINATION



**XL LENS !
12,3mm**



VIDEO

ULTI-G	Handpiece + charging base + Battery + 2 eye protectors + 100 sleeves.
ULTI-G500	Full protection sleeves x500
C-BAT	Universel charger (option)
BAT-3.7-L	additional battery (option)

MD class I



**4-YEARS
WARRANTY**

ENDY 7

EndoMotor with apex locator
Torque and speed control

Repair Index : **9/10**
Digital measurement: **100%**
Contra-angleweight : **25g**

A46619

Prijs: 1835,00

18% korting

**Promo
1490,00**



LA VALEUR SÛRE !
4-YEARS Warranty

EFFICIENT, RELIABLE, SIMPLE

- ▶ Programmable endomotor with integrated automatic apex locator. 9 memories.
- ▶ Speed control (200-800 rpm)
- ▶ Torque control (0.5-3.5 N.cm)
- ▶ Auto-disengage function & 2h quick charge.



ENDYC15 Control box, contra-angle handpiece, white cord (x2), labial electrode (x2), grip file for manual measurement, lubrication adapter, charger.

MD class IIa * 4-year warranty on control box | 2-year warranty on contra-angle and accessories

IONYX is part of
ELSODENT group