

# PROMOTIONAL OFFERS\*



Valid from 2<sup>nd</sup> of Feb till 30<sup>th</sup> of June 2026

Only for participating dealers Benelux.

All prices are in €, exclusive VAT.

# opus YOUR WAY

A51491

**Vraag de promoprijs !**

A51492

**Vraag de promoprijs !**

Opus cart

**Ontdek het OPUS-apparaat in  
uw praktijk en ervaar innovatie  
van dichtbij !**

**Vrijblijvende demonstratie via  
[marketing@denta.be](mailto:marketing@denta.be)**



# C20

## Easy to use. Easy to master.

More than just a tool, The C20 enhances the diagnostic process by creating an engaging and interactive experience for both the dental team and their patients.

By streamlining communication and boosting chairside efficiency, the C20 helps transform practices into high-performing, patient focused environments.



**Join the game  
OF MODERN  
DENTISTRY**



**Level up  
YOUR PRACTICE**



**Press start  
AND GO**




A51758

**Vraag de promoprijs !**

Ref. with AIS: S\_C20\_0002  
Ref. without AIS: S\_C20\_0001


# C50

## Daylight



Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

## Daylight+



**High contrasts** emphasizing anatomic fissures and details of the dentine structure / depth.

- Full HD: 1920x1080
- Ergonomic design with small tip (easy to use with small mouth/children)
- Autofocus
- Diagnostic feature: dental plaque and caries detection (C50 Fluo)

A49262

**Vraag de promoprijs !**

(daylight)

Ref. S\_C50\_0002 DAY LIGHT

A49261

**Vraag de promoprijs !**

(fluorescence technique)

Ref. S\_C50\_0001 FLUO



# newtron p5XS

1

Ref. F00246




1S

Ref. F00245



H3

Ref. F00369



P5XS b-LED incl. with the following tips:  
Ref. F62100

**47% korting**

A38077

Prijs: 2857,89

**Promo: 1490,00**



N°4

Ø22

Ref. MP3622PH

**48% korting**

**Package Deal PD0368**

Koop de Acteon **Newtron P5 XS B.LED** aan promoprijs

van **1490,00 euro** (waarde: 2857,89 euro) en

koop de **Mondspiegel Pure Reflect #4** of **#5**

aan **1 euro** (waarde: 56,53 euro)

## Newtron® handpiece



ref. F12281 - Autoclavable

A29690

Prijs: 651,37

**Promo: 268,50**

## Newtron® LED handpiece



ref. F12609

A32227

Prijs: 901,61

**Promo: 378,50**

## Newtron® Slim B. LED handpiece



• with blue LED ring: ref. F12900 A38204

• with white LED ring: ref. F12905 A39457

Prijs: 875,34 **Promo: 548,00**

## hygiene

### Hygiene Kit



ref. F00935


A31336

Prijs: 542,35

**Promo: 473,99**

## periodontics

### Perio Kit



ref. F00936

A33170

Prijs: 687,72

**Promo: 599,99**

## implantprotect

pure titanium

### Implant Protect Kit



ref. F02120

A39303

Prijs: 753,91

**Promo: 682,38**

## endosuccess

retreatment

### Endo Success Retreatment Kit



ref. F00737

A29752

Prijs: 719,91

**Promo: 632,40**

# piezotome cube

Dynamic power adaptation for  
safe bone surgeries

Ref. F50100

Delivered with a multifunction footswitch,  
a Piezotome® CUBE LED handpiece,  
1 handpiece holder, an Essential II kit  
(6 tips BS1S, BS4, SL1, SL2, SL3 et LC2),  
a Piezotome® dynamometric wrench,  
1 irrigation pole, 2 autoclavable irrigation  
lines, 5 clips and 15 perforators.



A41280

**Vraag de promoprijs !**



# implant center 2

## 3 modes:

- Piezotome® (bone surgery).
- Newtron® (perio, endo,...).
- I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.



**34% korting**

**PACKAGE DEAL PD0097**

Koop een **Implant Center 2** aan **promoprijs** van **9790,99 euro** (waarde: 14901,47 euro) en koop een **Newtron LED handstuk + Newtron LED snoer** aan **super promoprijs 1295,00 euro** (waarde: 2022,97 euro).

**34% korting**

A34090  
Prijs: 14901,47  
**Promo: 9790,99**

Ultrasonic control console for conventional dentistry and intraoral dentistry and motor control console for dental implantology  
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France



## Bone Surgery™ II Kit\*

Ref. F87509 - Tips BS1S II, BS2R II, BS2R II, BS4 II, BS5 II, BS6 II.

A34094 Prijs: 1468,12

**Promo: 1111,99**



## Sinus Lift II Kit\*

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

A34094 Prijs: 1104,84

**Promo: 1009,99**



## Intralift™ II Kit\*

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

A34089 Prijs: 1401,87

**Promo: 1111,99**



## Extraction™ II Kit\*

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

A34088 Prijs: 1114,76

**Promo: 1009,99**

\* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

## Accessoires

### Newtron® LED handpiece



A32227 Prijs: 901,61  
**Promo: 378,50**

Ref. F12609

For Piezotome® 2 and ImplantCenter™ 2

Class IIa medical devices (GMED) - CE 0459  
Manufacturer: SATELEC® a company of ACTEON® group - France

### Newtron® LED cord



A34077 Prijs: 922,81  
**Promo: 779,00**

Ref. F12830

For Piezotome® 2 and ImplantCenter™ 2

### Universal irrigation lines Sterile x10



A32384 Prijs: 188,67  
**Promo: 147,99**

Ref. F59905 - Single-use with integrated cassette and flow control system.

Class Is medical device (GMED) - CE 0459  
Manufacturer: SATELEC® a company of ACTEON® group - France

### PIEZOTOME® LED handpiece + cord



A46435 Prijs: 1676,00  
**Promo: 1324,99**

Ref. F12815

Class IIa medical devices (GMED) - CE 0459  
Manufacturer: SATELEC® a company of ACTEON® group - France



### Universal II Kit\*

Ref. F87539 - Tips BS6, LC1, LC2, CE1, CE3 et PZ1.

A46420 Prijs: 1428,27  
**Promo: 1111,99**



### Essential II Kit\*

Ref. F87528 - Tips BS1S II, BS4 II, SL1 II, SL2 II, SL3 II, LC2 II.

A41909 Prijs: 1468,12  
**Promo: 1111,99**



### Crest Splitting II Kit\*

Ref. F87567 - Tips CS1 II, CS2 II, CS3 II, CS4 II, CS5 II, CS6 II.

A39232 Prijs: 950,34  
**Promo: 938,70**



### Crown Extension II Kit\*

Ref. F87554 - Tips BS6 II, CE1 II, CE2 II, CE3 II.

A43497 Prijs: 719,03  
**Promo: 630,99**

\* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

# miniled ortho 2

32% korting

A34280

Prijs: 1486,25

Promo: 999,00

- **Customizable signal**

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

- **Power**

3000mW/cm<sup>2</sup>  
(Ø5,5mm light guide).

- **Long lasting battery**

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

- **Ortho friendly**

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

- **Efficiency & safety**

Cure any type of material used orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.



## Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from  
few brackets to a whole quadrant

# Simply Powerful

Outperforms in  
power, efficiency and speed

**33% korting**

A34219

Prijs: 1435,51

**Promo: 950,00**



# miniled supercharged

## • Visibility

Luminous screen with customizable modes and times, count curing time and load indicator.



## • Power

2000mW/cm<sup>2</sup>  
(Ø7,5mm light guide).

3000mW/cm<sup>2</sup>  
(Ø5,5mm light guide - optional).

## • Optimal result


3 modes to fit all clinical needs.



### FAST MODE

**The choice of rapidity**


Full power during 3, 4, 5 or 10 sec, according to the composite thickness.



### PULSE MODE

**Limits temperature rise**

5 or 10 flashes of 1 sec for a polymerization in security.



### RAMPING MODE

**Avoids composites retraction**

Progressive power increase during 9 or 15 sec to ensure perfect sealing of restoration.

# miniled standard

38% korting

A27384

Prijs: 1117,47

Promo: 690,00

- **Fast**

Only 6 to 12 sec to polymerize 2mm composite.

- **Power**

1250mW/cm<sup>2</sup>  
(Ø7,5mm light guide).

2000mW/cm<sup>2</sup>  
(Ø5,5mm light guide - optional).

- **Safety**


Base-charger includes a radiometer for light power intensity check-up.

- **Optimal result**

3 modes to fit all clinical needs.


## Highly Adaptable

The alliance of power & safety




### FAST MODE

Fast and easy  
Full power during 10 sec  
(bip at 5 sec) on regular depth.



### PULSE MODE

Ally rapidity & safety  
10 flashes of 1 sec  
(bip after 5 flashes)  
with times off to  
preserve dental pulp.



### RAMPING MODE

No stress cracking  
Progressive power  
increase during 20 sec  
(bip every 5 sec) for  
a complete and safe  
curing.



# servotome

Delivered with a single footswitch,  
a box with 10 electrodes,  
a conductive bracelet.

Ref. F57210



**28% korting**

A30972

Prijs: 2076,45

**Promo: 1490,00**

Tot 60% korting

Package Deal PD0335

# Discover our RISKONTROL® range

Koop **10 x Riskontrol classic** en/of **Riskontrol art** naar keuze aan **promoprijs** en koop een **Riskontrol adapter** aan **1 euro** naar keuze (waarde: 115,90 euro).

Adapt your tips to the taste of your patients



art by **riskontrol**

V38000  
Prijs: 62,76  
Promo: 32,99

-  Ref. 201 700 - Mint
-  Ref. 201 701 - Mandarine
-  Ref. 201 702 - Liquorice
-  Ref. 201 703 - Anise
-  Ref. 201 704 - Black currant

Non-sterile disposable tips (x 250)

Add color

to your patients' care experience



V10811  
Prijs: 57,78  
Promo: 27,50

classic by **riskontrol**

-  Ref. 201 664 - White
-  Ref. 201 665 - Blue
-  Ref. 201 666 - Yellow
-  Ref. 201 667 - Green
-  Ref. 201 668 - Pink

Non-sterile disposable tips (x 250)

Protect entirely your air/water syringe thanks to the plastic sheath



total protec by **riskontrol**

V33150  
Prijs: 33,59  
Promo: 29,99



-  Ref. 201 750
-  Ref. 201 751

Non-sterile disposable tips with plastic sheath (x 100)

Choose a totally sterile workflow



A33761  
Prijs: 38,89  
Promo: 30,50

sterile by **riskontrol**

-  Ref. 201 770

Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile by RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - FRANCE.



# Mouth mirrors

## PURE REFLECT®

- Very bright image
- Resistant

For all purposes

**RH**

- A greyish reflection



**TI**

- A blueish reflection

## FRONT SURFACE

Cross-sectional view





Image reflected directly on the surface of the glass, distortion-free.


**NO DOUBLE IMAGE**

## BACK SURFACE

Cross-sectional view of an aluminium mirror



Reflecting surface under the glass. **DOUBLE IMAGE**



## PURE REFLECT®

N°	Ref.	N.
4	MP3622PH	12 A42193
5	MP3624PH	12 A42194

Prijs: 56,53

**Promo: 39,99**



## A TRIPLE ACTION



Kaolin clay



Aluminum chloride

### 1 ASTRINGENT EFFECT

helps to completely dry the sulcus with no blood or fluids.

### 2 OPEN THE SULCUS

effectively with no pain or damage to tissues thanks to the viscosity carefully calculated with our clinical expertise.

### 3 PRESERVE EPITHELIAL ATTACHMENT\*

for a reduced risk of gingival recession and bone resorption.

*\*Canula must not be inserted into the sulcus*

**expasyl  
exact**



A44322  
Prijs: 147,22

**Promo: 99,85**

Pas leverbaar  
vanaf april 2026

**Paste for temporary sulcular opening**  
20 capsules (0,3g) Expasyl™ Exact - Ref. 261010



**Paste for temporary sulcular opening**  
20 capsules (1g) Expasyl™ - Ref. 261030

**expasyl**

A10829

Prijs: 248,87

**Promo: 127,50**



**Straight** - Ref. 261005 A26105  
**Bended** - Ref. 261015 A33371

Prijs: 111,16

**Promo: 79,99**

# x-mind dc



Bottom  
mounted

The X-MIND® generators can accommodate any operatory configuration.

3 arm lengths are available : 41 cm (16.2"), 82.5 cm (32.5") and 110 cm (43.5") and can be:

- Top wall mounted
- Bottom wall mounted

Dental unit, ceiling and mobile version are also available.

With its adaptable wall plate, it is easy to replace an older generator with the X-MIND® without having to do additional drilling on the wall.



**40% korting**

A37323-11B

Prijs: 4821,80

**Promo: 2849,00**

## Available with **Remote Wireless Control**

**Wireless Button Kit X-MIND DC**

Ref. 5.15.00135

All X-Mind DC devices delivered  
from April 2021 onwards

A50693



### No more wiring issues

Say goodbye to the complicated and expensive wiring work traditionally needed to connect the remote button in your practice.



### Cost effective

Save on installation expenses with the easy-to-integrate wireless solution.



### Enhanced convenience

Enjoy the freedom and flexibility of a wireless connection, making your practice more efficient and user-friendly.



Wall version	Top mounting	Top mounting square cone	Bottom mounting	Bottom mounting square cone
0,40 m extensor	Ref. 5.92.23137	Ref. 5.92.23167	Ref. 5.92.23130	Ref. 5.92.23160
0,80 m extensor	Ref. 5.92.23138	Ref. 5.92.23168	Ref. 5.92.23131	Ref. 5.92.23161
1,10 m extensor	Ref. 5.92.23139	Ref. 5.92.23169	Ref. 5.92.23132	Ref. 5.92.23162

All references for the X-MIND® DC


Choose from a variety of cones to better fit to your clinical needs



**Long cone** (31 cm - 12'')  
Recommended for the parallel technique



**Short cone** (20 cm - 8'')  
Recommended for bisector technique



**Rectangular Long cone** (31 cm - 12'')  
Reduces the X-rays applied to the patient's skin by 50%.

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy

# x-mind unity

A39620

Prijs: 5882,40

**Promo: 3419,00**

**41% korting**

## A sharp and contrasted image

- A contrasted and defined image
- Optimized protection for the patient
- The generator and sensor association

Fixing from above, arm 0.80 m, short cone.  
Other configurations: see public price 2024.

## Intraoral X-ray

- Very high quality
- Precise and flexible positioning
- Fast selection from traditional to digital mode
- Minimum exposure time: 0,02s
- Programmable timer
- Top or bottom mounting



Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy

# Meet the X-MIND<sup>®</sup> optima 3D Family.

A51759

**Vraag de promoprijs !**

## X-MIND<sup>®</sup> optima 2D/3D wall 100-240V

2D/3D PAN device 5x5 to 12x10 wall installation  
Ref. W1203005



A51760

**Vraag de promoprijs !**

## X-MIND<sup>®</sup> optima 2D/3D ceph 100-240V

2D/3D PAN CEPH device 5x5 to 12x10 floor installation  
Ref. W1203010



A51761

**Vraag de promoprijs !**

## X-MIND<sup>®</sup> optima 2D/3D floor 100-240V

2D/3D PAN device 5x5 to 12x10 floor installation  
Ref. W1203005.FL

New

AI



# Meet the 2D Prime Family

X-MIND® prime is available in three configurations: 2D wall mounted, 2D Ceph and 2D floor to suit all installation and clinical demands.

**Vraag de promoprijs !**



## X-MIND® prime WALL 110-240V

Order code: W1200001 A48189

The ACTEON®'s smallest extraoral 2D imaging solution ever for the maximum **accessibility in any context.\***

with one sensor: A48191

with two sensors: A51763



## X-MIND® prime CEPH 110-240V

Order code: W1200009

Cephalometric add-on for a full 2D system still encompassing the same philosophy of prime, **ideal for pediatric applications.**




## X-MIND® prime FLOOR 100-240V

Order code: W1200001.FL A51762


\*Belgian market only (W1200001/W120001.PB)

# Rely on a comprehensive set of 2D exams


X-MIND® prime provides a comprehensive range of 2D exams designed to meet the diverse needs of general dentistry, including dental panoramic imaging, temporomandibular joint assessments and maxillary sinus evaluations.




Dental panoramic



Maxillary sinuses



Temporomandibular joint



Child panoramic



Half-panoramic  
Left & right



Bitewing  
Left & right



Frontal

# PSPiX<sup>2</sup>

## is AUTONOMOUS

The fully automated process free yourself from constraints

Just drop the plate inside the slit, the PSPiX<sup>2</sup> will automatically detect the size, scan the plate, optimize the image, and send it instantly to your imaging software.

In the meantime, it will automatically erase and eject the plate for an immediate reuse or switch to standby mode.

**STUNTPRIJS**  
**4890,00 euro**  
*Bij inruil van eender welke fosforscanner !*

## Fully automated process

30% korting

A51747

Prijs: 8525,00

Promo: 5890,00

Ref. S\_702\_0013



1

### Auto-access:

The door only opens when the imaging plate is inserted.

2

### Auto-detect:

PSPiX<sup>2</sup> automatically detects the size of the inserted imaging plate.

3

### Auto-scan:

PSPiX<sup>2</sup> scans and optimizes the image before displaying it into the right workstation.

4

### Auto-eject:

The imaging plate is automatically erased and ejected for an immediate reuse.

5

### Eco-mode:

PSPiX<sup>2</sup> automatically switches to economic mode when not in use for an extended period of time.

## Protective bag for ERLM\*

For **PSPIX<sup>®2</sup>** and **PSPIX<sup>®</sup>**  
(previous generation)

A36130  
Prijs: 85,82  
**Promo: 65,99**



	<b>PSPIX<sup>®</sup></b> (previous generation)	
	Box of 250	
Size 0		700338L
Size 1		700339L
Size 2		700340L
Size 3		700341L

## ERLM plates\*

For **PSPIX<sup>®2</sup>**

V39947  
Prijs: 131,16  
**Promo: 96,99**



	<b>PSPIX<sup>®2</sup></b>	
	Box of 2	
Size 0 (22 x 35mm)		990215
Size 1 (24 x 40mm)		990216
Size 2 (31 x 41mm)		990217
Size 3 (27 x 54mm)		990218

## Protections for ERLM\*

Protective cover for  
**PSPIX<sup>®</sup>** and **PSPIX<sup>®2</sup>**  
(previous generation)



Compatible with	<b>PSPIX<sup>2</sup></b>	<b>PSPIX</b>
PSPIX <sup>®</sup>	✓	✓
SOREDEX <sup>®</sup> Digora	✓	✓
INSTRUMENTARIUM <sup>®</sup> Express	✓	✓
DÜRR DENTAL <sup>®</sup> Vistascan <sup>®</sup>	✓	✓
KAVO <sup>®</sup>	✓	✓

A36129  
Prijs: 48,06  
**Promo: 26,99**

	<b>PSPIX<sup>®</sup></b> (previous generation)	
	Box of 250	
Size 0		700342L
Size 1		700343L
Size 2		700344L
Size 3		700345L

## Protective bag for ERLM\*

For **PSPIX<sup>®2</sup>**

V39564  
Prijs: 100,86  
**Promo: 70,99**



	<b>PSPIX<sup>®2</sup></b>	
	Box of 300	
Size 0		700546
Size 1		700547
Size 2		700548
Size 3		700549

- Protection against scratches, light and cross contamination
- Ergonomic and easy to use
- Economic solution

\* Luminescent Radio Memory Display for digital phosphor plate systems

PSPIX<sup>®</sup> imaging plate - Class IIa Medical devices (GMED) - CE0459

Manufacturer: Sopro ACTEON<sup>®</sup> Group - France

Bag & cover: Class I medical devices - CE

Manufacturer: Sopro ACTEON<sup>®</sup> Group - France




**The thin and flexible imaging plates, available in 4 sizes, will allow you to easily scan all morphologies...**



Child, Adult Periapical	Adult, Bitewing Periapical	Child	Bitewing	Occlusal
				
<b>Size 0</b> 22x35mm	<b>Size 1</b> 24x40mm	<b>Size 2</b> 31x41mm	<b>Size 3</b> 27x54mm	<b>Size 4</b> 69x54mm

**... and all clinical cases**

<b>Periodontics</b> 	<b>Endodontics</b> 	<b>Orthodontics</b> 	<b>Pedodontics</b> 	<b>Cariology</b> 
<b>Bitewing</b> 	<b>Implantology</b> 	<b>Occlusal</b> 		

Recupel-bijdrage: 0,10. Prijzen zijn in euro exclusief btw en per stuk of verkoopseenheid. Druk- en zetfouten voorbehouden. Foto's zijn niet bindend. Bij bestellingen boven de 150 euro exclusief btw worden geen administratie - en verzendkosten aangerekend. Toestellen boven de 5000 euro (excl. btw) worden gratis geleverd en geplaatst. Plaatsing omvat steeds de levering, installatie (op de benodigde aansluitingen en leidingen, te voorzien door de klant zelf volgens onze instructies) de oplevering en een korte opleiding. Promoties geldig van **02/02/2026** tot **30/06/2026** en zolang de voorraad strekt. Acties niet cumuleerbaar met andere promoties.

Denta nv, Heiveldekens 2, 2550 Kontich, [www.denta.be](http://www.denta.be)