


#whdentalwerk



wh.com



PROMOTIONS

FOR DENTA

Valid from

01.02.2026 to
30.06.2026

NEW!



Implantmed Plus II surgical motor
with or without piezoelectric technology

www.denta.be
03 450 93 20
vente@denta.be



Synea Fusion

The Synea Fusion is the world's most widely used contra-angle handpiece – quiet, strong and durable.

The air turbines in this series are equipped with ceramic bearings, a quadruple spray system, and are available with various head sizes.

Safety first: All Synea turbine handpieces are equipped with a technology that ensures that aerosols do not enter the instrument.

synea FUSION



The powerhouse of contra-angle handpieces!

Synea Fusion Power Edition

• WG-900 LT

Red contra-angle handpiece 1:4
with light, Ø 10.7 mm
triple spray
36*-month warranty

REF 30459000

A50376

Promo: 799,00
Par 3 : 699,00 pièce



- Extremely powerful and robust contra-angle handpiece
- For the removal of zirconium and attachments
- Optimal power to speed ratio thanks to 1:4 speed ratio
- Indispensable for difficult restorations
25% longer lifespan

• WG-56 LT

Blue contra-angle handpiece 1:1
with light, Ø 9.5 mm
single spray
36*-month warranty

REF 30009000

A37655

Promo : 535,00

Par 3 : 455,00 pièce



• WG-66 LT

Green contra-angle handpiece 2:1
with light, Ø 9.5 mm
single spray
36*-month warranty

REF 30011000

A37654

Promo : 615,00

Par 3 : 525,00 pièce



• WG-99 LT

Red contra-angle handpiece 1:5
with light, Ø 9.5 mm
Quattro spray
36*-month warranty

REF 30007000

A37653

Promo : 795,00

Par 3 : 615,00 pièce



Good spray and lighting

• HG-43 A

Blue handpiece 1:1
no light
single spray

24*-month warranty

REF 30013000

A37656

Promo : 358,00,00

Par 3 : 315,00 pièce



TG-98 L* for Roto Quick

TG-98 LM** for Multiflex®

Air turbine *LED+, **with light,
29 W**, Ø 10 mm,
Quattro Spray

24*-month warranty

*REF 30001000

**REF 30003000

V37658

Promo : 635,00

Par 3 : 525,00 pièce



TG-97 L* for Roto Quick

TG-97 LM** for Multiflex®

Air turbine *LED+, **with light,
29 W**, Ø 11,5 mm,
Quattro Spray

24*-month warranty

*REF 30004000

**REF 30006000

V37657

Promo : 635,00

Par 3 : 525,00 pièce



Cast from a single block, no seams



NEW
More value for
no added
cost!

Now with an anti-scratch surface!



Synea Vision

High-end contra-angle handpieces and air turbines; lightweight, quiet, with excellent grip, a slim head and a long warranty. The air turbines in this series are equipped with ceramic bearings, a quintuple spray system, a scratch-resistant surface and are available with four different head sizes.

synea VISION

Ring LED+

Work without shadows!



• WK-56 LT

Blue contra-angle handpiece 1:1

with light, Ø 9.5 mm
single spray

36-month warranty

REF 30025000

A37580

Promo : 699,00



• WK-66 LT

Green contra-angle handpiece 2:1


with light, Ø 9.5 mm
single spray

36-month warranty

REF 30026000

A37596

Promo : 739,00



Synea Vision Power Edition

• WK-900 LT

Red contra-angle handpiece 1:4

with light, Ø 10.7 mm
triple spray

36-month warranty

REF 30418000

A50377

Promo : 949,00



With a narrow head

• WK-93 LT

Red contra-angle handpiece 1:4.5

with light, Ø 9.1 mm
quintuple spray
ceramic bearings

36-month warranty

REF 30024000

A37549

Promo : 878,00



• WK-99 LT

Red contra-angle handpiece 1:5

with light, Ø 9.5 mm
quintuple spray
ceramic bearings

36-month warranty

REF 30023000

A37652

Promo : 878,00



Extremely powerful and robust contra-angle handpiece

For the removal of zirconium and attachments

Optimal power to speed ratio thanks to 1:4 speed ratio

Indispensable for difficult restorations

***25% longer lifespan**



Optimal turbine design
for maximum power

Very narrow head

TK-94 L* for Roto Quick
TK-94 LM** for Multiflex®

**Air turbine *Ring LED+
 with light

14 W, Ø 9 mm
 quintuple spray

24-month warranty

*REF 30016000
 **REF 30031000



V38968
Promo : 899,00

TK-97 L* for Roto Quick
TK-97 LM** for Multiflex®

**Air turbine *Ring LED+
 with light

21 W, Ø 10 mm
 quintuple spray

24-month warranty

*REF 30017000
 **REF 30018000



V38966
Promo : 899,00

TK-98 L* for Roto Quick
TK-98 LM** for Multiflex®

**Air turbine *Ring LED+
 with light

36 W**, Ø 11.5 mm

24-month warranty

*REF 30019000
 **REF 30020000



V38970
Promo : 855,00

TK-100 L* for Roto Quick
TK-100 LM for Multiflex®

**Air turbine *Ring LED+
 with light

36 W**, Ø 13 mm

24-month warranty

*REF 30021000
 **REF 30022000



V38971
Promo : 899,00

Syne Vision Short series; 20% lighter and shorter, yet just as powerful as standard Vision instruments!

HK-43 LT

**Blue handpiece 1:1,
 with light,
 single spray**

36-month warranty

REF 30028000

A37641
Promo : 629,00



WK-56 LT S

**Blue contra-angle
 handpiece 1:1
 with light, Ø 9.5 mm
 single spray**

36-month warranty

REF 30217000

A39134
Promo : 699,00



WK-99 LT S

**Red contra-angle
 handpiece 1:5
 with light, Ø 9.5 mm
 quintuple spray**

36-month warranty

REF 30215000

A50661
Promo : 878,00



WK-93 LT S

**Red contra-angle
 handpiece 1:4.5
 with light, Ø 9.1 mm
 quintuple spray**

36-month warranty

REF 30216000

V39175
Promo : 878,00



**Lightweight
 and
 compact!**

Syne Vision Short series; 20% lighter and shorter, yet just as powerful as standard Vision instruments!

Implantmed Plus SI-1023


The powerful **Implantmed Plus** from W&H guarantees a high level of safety during treatment for both you and your patients. With its intuitive user interface, colour touchscreen and glass surface, the Implantmed offers a truly user-friendly experience. The drive is ideal for surgical procedures with its automatic thread-cutting function, automatic torque control and adjustable parameters. For **piezoelectric surgery**, the **Piezomed Plus Module** can be fully integrated with the Implantmed Plus surgical motor, as shown in the photo.

Both devices can be operated with a single foot control (which can be wireless), and the instruments connected can be supplied with the same physiological saline solution.

implantmed



Work according to protocol



with automatic thread-cutting function



**Piezomed Plus
Module**



Implantmed Plus

Package A **Implantmed Classic SI923 + 2x WI-75 E/KM**

- › Implantmed Classic no light
 - › Wired foot control
 - › Surgical contra-angle handpiece
2 x WI-75 E/KM (20:1, green)
- REF 90000161
REF 10207550
REF 10207550



Package Deal PD0445

Package B **Implantmed Plus SI1023 + WS-75**

- › Implantmed Plus no light
 - › Wired foot control
 - › Surgical contra-angle handpiece WI-75 E/KM (20:1, green)
- REF 90000163
REF 10207550



Package Deal PD0286

Piezomed Plus **Module SA-435 M**

The Piezomed Plus Module can be connected to the Implantmed Plus surgical motor via plug and play. The most common oral surgical treatments can be performed with this combination.

Features:

- › tip recognition
 - › LED light
 - › handpiece with a 1.8-metre cable
- REF 90000350

A47782
Promo : 4350,00



piezomed

Package C **Implantmed Plus SI1023 + WS-75 L**

- › Implantmed Plus with light
 - › Wired foot control
 - › Surgical contra-angle handpiece WS-75 L (20:1 green)
- REF 90000165
REF 30033000



Package Deal PD0056

Implantmed Plus II

The new Implantmed Plus II is the next evolution of the SI-1023, featuring the following key characteristics:

- › All motors are equipped with light
- › Large, 7 inch touchscreen
- › Wireless foot control at no extra cost
- › Easy connection to the Piezomed module Plus II:
 - No additional cables required
 - No additional cooling required
 - No additional foot control required
- › Easy connection to the W&H ioDent platform for software updates and patient documentation via built-in Wi-Fi
- › Option to scan all materials and generate a comprehensive implant report
- › Connection with Osstell ISQ measurement and reporting possible and straightforward

Package D Implantmed Plus II SI-2100 + WS-75 L

- › Implantmed Plus II with light
 - › Wired foot control
 - › Surgical contra-angle handpiece
WS-75 L (20:1, green)
- REF 90000418
REF 30033000



Package Deal PD0442

Package E Implantmed Plus II SI-2101 + WS-75 L

- › Implantmed Plus II with light
 - › Wireless foot control
 - › Surgical contra-angle handpiece
WS-75 L (20:1, green)
- REF 90000419
REF 30033000



Package Deal PD0443

Package F Implantmed Plus II SI-2102 + WS-75 L

- › Implantmed Plus II with light
 - › Scan edition
 - › Wireless foot control
 - › Surgical contra-angle handpiece
WS-75 L (20:1, green)
- REF 90000420
REF 30033000



Package Deal PD0444



NEW!

The new
Implantmed Plus II
surgical motor;
can be supplemented
with piezoelectric
technology if desired!



implantmed^{PLUS} II
piezomed^{PLUS} II
module

Flightcase

› For the Implantmed Plus,
Piezomed Module and
accessories (in combination
with the equipment)

REF 08237370

A52050

Promo : 499,00



Surgical motor only



Surgical motor with piezoelectric module



Surgical handpieces & contra-angle handpieces

The instruments below boast lights when connected to an EM19 LC motor.
When connected to an EM19 motor there is no light source.

• WS-75 L

With light.
Implant site preparation
and implant placing

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 30033000

A40683
Promo : 948,00



• WS-56 L

With light.
Osteotomy, sequestrotomy,
germectomy, hemisection

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30040000

A40698
Promo : 999,00



• WS-91 L

With light.
Specially designed to
improve access to the
premolar and molar area

1:2.7, contra-angle handpiece
for surgical drills (red)
FG, with a 45° head tilt

REF 30049000

A34809
Promo : 999,00



• WS-75 LG

Boasts its own lighting
Light output depends on speed.
Applications such as WS-75 L

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 30032000

A34809
Promo : 1089,00



Can be
disassembled
without tools
for optimal
hygiene

• WS-92 L

With light.
Osteotomy, sequestrotomy,
germectomy, hemisection

1:2.7, contra-angle handpiece
for surgical drills (red)
FG

REF 30043000

A41640
Promo : 999,00



• S-11 L

With light.
Apical surgery, osteotomy,
bone modelling and
wisdom tooth surgery

1:1, handpiece for
surgical drills (blue)
Ø 2.35 mm

REF 30058000

A41366
Promo : 955,00



• S-11 LG

Boasts its own lighting
Light output depends on speed.
Applications such as S-11 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30057000

A39292
Promo : 4350,00



With mini-LED+
and optimal cooling

See wh.com for more models and versions

• WS-56 LG

Boasts its own lighting.
Light output depends on speed.
Application as WS-56 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30039000

A51712

Promo : 1119,00



• WS-91 LG

Boasts its own lighting.
Light output depends on speed.
Application as WS-91 L

1:2.7, contra-angle handpiece
for surgical drills (red)
FG, with a 45° head tilt

REF 30048000

A39226

Promo : 1119,00



• WS-92 LG

Boasts its own lighting.
Light output depends on speed.
Application as WS-92 L

1:2.7, contra-angle handpiece
for surgical drills (red)
FG

REF 30042000

A48939

Promo : 1119,00



• WS-56

No light.
Applications such as WS-56 L

1:1, contra-angle handpiece for
surgical drills (blue), Ø 2.35 mm

REF 30038000

A28896

Promo : 798,00



• WS-91

No light.
Applications such as WS-91 L

1:2.7, contra-angle handpiece
for surgical drills (red)
FG, with a 45° head tilt

REF 30047000

A49284

Promo : 898,00



• S-11

No light.
Applications such as
S-11 L

1:1, contra-angle handpiece
for surgical drills (blue)
Ø 2.35 mm

REF 30056000

A01882

Promo : 698,00



• WS-75

No light.
Applications such as
WS-75 L

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 30065000

A01881

Promo : 898,00



• WI-75 E/KM

No light.
Applications such as WS-75 L
Cannot be disassembled.

20:1, contra-angle handpiece
for surgical drills (green)
Ø 2.35 mm

REF 10207550

A28895

Promo : 725,00



• WS-92

No light.
Applications such as WS-92 L

1:2.7, contra-angle handpiece
for surgical drills (red), FG

REF 30041000

A31341

Promo : 898,00



See wh.com for more models and versions

Instrument Kit Bone

Surgery Tips included:
B1, B2R, B2L, B3, B4, B5

REF 07173100

A40330

Promo : 999,00



Instrument Kit Sinus Special

For sinus lifts.
Tips included:
S1, S2, S3, S4, S5

REF 07173200

A40331

Promo : 705,00



Instrument Kit Endodontics

Tips included:
R1D, R2RD, R2LD, R3D,
R4RD, R4LD

REF 07173300

A51716

Promo : 790,00



Instrument Kit Periodontal

Tips included:
P1, P2RD, P2LD

REF 07173400

A51717

Promo : 444,00



Instrument Kit Extraction

Tips included:
EX1, EX2

REF 07173500

A40339

Promo : 311,00



Instrument Kit Implant Crestal P

Posterior. Tips included:
I1, I2P, I3P, I4P,
Z25P, Z35P

REF 07997110

A51718

Promo : 1375,00



Instrument Kit Implant Crestal A

Anterior. Tips included:
I1, I2A, I3A, I4A,
Z25P, Z35P

REF 07997090

A5171

Promo : 1375,00



Fast
Precise
Powerful

Instrument B6 & B7, surgical saws for,
among other things, apical resection and the
removal of root remnants or alveolar bone

B6 REF 07022500

A440330

Promo : 999,00

B7 REF 07022600

A40331

Promo : 705,00



Piezomed Plus module

W&H piezoelectric technology brings all the benefits of modern ultrasonic technology within reach: high-frequency vibrations facilitate extremely precise bone surgery. What is more, thanks to the cavitation effect you can maintain a virtually blood-free surgical site, while also minimising tissue damage.

- Sinus lift
- Bone harvesting
- Apical resection
- Retrograde endodontics
- Extraction
- Periodontal (sterile)
- Etc.

Visit video.wh.com for all applications

piezomed

Courses on the piezoelectric technique in dental surgery:


Internationally renowned dental implantologist Irfan Abas will teach you the numerous possibilities of the piezoelectric technique in modern dental surgery.

www.abas-academy.com

 **ABAS**
implant academy



The Dutch training institute Dental Best Practice also offers an intensive, hands-on piezoelectric course as part of their programme. For more information, visit dentalbestpractice.nl



Surgical saw handpieces

Fast, powerful and durable

The saw handpieces are specially designed for conversion osteotomies and bone harvesting. A consistently high cutting performance is guaranteed even after many uses at maximum speed.

The saw blades are available in different sizes and are supplied sterile.

S-8 R

Saw handpieces with reciprocating motion

REF 10100801

A33125

Promo : 1299,00



S-8 O

Saw handpieces with oscillating motion

REF 10100802

A29430

Promo : 1549,00



S-8 S

Saw handpieces with sagittal motion

REF 10100803

A49283

Promo : 1499,00




Easy to disassemble

and clean really effectively, ensuring a long service life and optimal patient safety

Osstell Beacon

Implant Stability Quotient (ISQ)

- › Monitoring of osseointegration
- › Makes an earlier prosthetic process possible
- › Improved predictability of treatment for patients with risk factors
- › Easy communication with patients through visualisation of the growing-in phase
- › Maximum flexibility thanks to wireless technology
- › The Osstell Beacon can be integrated with the latest generation of W&H Implantmed



Work in a validated and documented way; ideal for referral practices



Osstell Beacon

REF 103000

A43554

Promo : 2129,00



IMPLAN TOLOGY



Osstell Classic

A51729

Promo : 1699,00

REF 08086148

The Osstell Classic has the same functionalities as the Osstell Beacon, but without the option to use the Osstell Connect Portal, and it cannot be connected to W&H Implantmed.

Includes free use of the Osstell Connect portal:

- › build up the implant history
- › compare implants
- › compare patients
- › compare pathways
- › compare problems

SmartPegs are necessary in order to measure osseointegration. These are available from your dealer, see wh.com for an overview.



Proxeo Twist

Top quality handpieces and contra-angle handpieces for prophylaxis

For your treatment unit



More space
with the
**LatchShort
attachment**

WP-66 W

**Prophy 4:1 contra-angle
handpiece**

for W&H LatchShort
prophy cups

REF 30189000

A42771

Promo : 345,00

Par 3 : 259,00 pièce



WP-64 MU

Prophy 4:1 contra-angle handpiece

with snap-on adapter and key
for all standard screw-in
prophy cups and brushes

REF 10216443

A28731

Promo : 345,00

Par 3 : 259,00 pièce



WP-64 M

**Prophy 4:1 contra-angle
handpiece**

for System Young Screw-in
prophy cups and brushes.
With key.

REF 10216403

A02406

Promo : 345,00

Par 3 : 259,00 pièce



WP-66 M

**Prophy 4:1 contra-angle
handpiece**

for all prophy cups, brushes
and snap-on adapters with a
2.35 mm shaft

REF 30190000

A39034

Promo : 345,00

Par 3 : 259,00 pièce



HP-44 M

Prophy 4:1 handpiece
for Doriot system
disposable prophy cups

REF 1014403

A02320

Promo : 279,00

Par 3 : 219,00 pièce



Proxeo Twist

Cordless Polishing System



Stand-alone polishing contra-angle handpiece

- › Work flexibly without cables
- › Easily regulate speed with wireless foot control
- › Perfect control of speed for precision cleaning and efficient polishing
- › Three different attachments
- › Minimal use of aerosols thanks to rotary polishing

Handy for working on location!

V46570
Promo : 999,00

PL-40 HW + PL-66 W

Proxeo Twist Cordless
 + LatchShort cups

REF 90000346



PL-40 HW + PL-64

Proxeo Twist Cordless
 + screw-in for
YOUNG cups

REF 90000347



PL-40 HW + PL-44

Proxeo Twist Cordless
 + Doriot attachment for
 disposable contra-angle
 handpieces with cup

REF 90000348



Wireless foot control

Lisa Mini

Compact, fast and easy to operate!

Just like the renowned Lisa large autoclave, this smaller version B autoclave is fast and features an easy-to-use touchscreen. Available with a 5- or 3-litre chamber.

Key features:

- › Fast contra-angle handpiece cycle in just 10 minutes (5-litre chamber)
- › Standard cycle time of 22 minutes (5-litre chamber)
- › Data can be logged via USB, LAN or the W&H Steri App
- › Can be connected to a demineralised water supply
- › A label printer can also be connected
- › Available in 2 or 3-tray format (5-litre chamber)

lisa MINI

Compact,
stackable,
and fast!

Découvrez le Lisa Mini
dans notre showroom !



5 litres

REF 19935100

A51721

Promo : 4708,00

5 litres*

REF 19945100

A52653

Promo : 4247,00

3 litres*

REF 19944100

A50654

Promo : 3999,00

*requires compressed
air connection



Handy as an
extra autoclave,
when you have
limited space
or for small
volumes!



Other items can be placed on top of the unit to make the most of your sterilization space. Comes with 2 trays as standard.

Lisa

**Découvrez le Lisa
dans notre showroom !**

High-end type B autoclave with extremely short cycle times.

Standard cycle time of 30 minutes, fast contra-angle handpiece cycle in 13 minutes. Features Eco Dry+ technology that adjusts the drying time to the weight of the load. Boasts an LED display and a touchscreen. Available with a 17- or 22-litre chamber. Data can be logged via USB, LAN or the W&H Steri app. The sterilization process can be tracked and documented down to the individual instrument or instrument kit via EliTrace. No extra software or computer needed. Thanks to smart technology, the Lisa warns you when your load is overheated, even if the door is open – no more burnt fingers!

**Fast
and
quiet!**

17 litres

REF 19931180

A42621

Promo : 5990,00

22 litres

REF 19932180

A42687

Promo : 6390,00



**The Lisa from
W&H is the
most widely used
autoclave in the
world, and for
good reason!**

lisa



Lara, Lara XL and Lara XXL

The fastest 28-litre class B autoclave on the market!

- › Tried and trusted technology
- › Large capacity
- › Various upgrades available
- › Easy to operate

The Lara is also available in 17- and 22-litre versions:

17 litres

REF 19941120

A42622
Promo : 5446,00

22 litres

REF 19942120

A42686
Promo : 5935,00

Lara XXL 38 litres

REF 19936110

A51724
Promo : 4708,00

Lara XL 28 litres

REF 19933110

A48819
Promo : 6990,00

**38-litre capacity for
up to 12 kg per cycle**



**Fits
60 cm trays**



**Highly suitable
for surgical cassettes
thanks to its large size**

lara^{XL}
lara^{XXL}

Lina

Type B autoclave with best-in-class cycle times.

- › Optimal sterilization programmes
- › The Eco B sterilization programme processes 0.5 kg in 25 minutes
- › Lina can reprocess loads of up to 6 kg to satisfy all the requirements of modern dental clinics
- › The B Universal 121°C cycle even allows you to reprocess delicate objects, including porous objects such as surgical clothing



17 litres

REF 19951131

A51722

Promo : 5299,00

22 litres

REF 19952131

A51723

Promo : 5599,00

**Excellent
value for
money!**



The **Activation Code** system allows you to activate additional features to meet your practice's changing needs and future requirements.

lina

Helix test kit

Complete set; consisting of the holder (the "snail shell") and 250 test strips.

Each test strip has a QR code that can be scanned and saved in our Steri app or in your network. Each code is linked to a sterilization cycle, which validates the process.

With QR-coded strips

REF T800207X

A34926

Promo : 153,90



Q-dry

Perfect solution for internal drying of contra-angle handpieces and air turbines.

Using the Assistina adaptors and air pressure, you can easily dry instruments prior to sterilization.

See wh.com for how-to videos.

q-dry

REF 19723118

A51689

Promo : 173,85



Water distiller

4-litre capacity,
enough for
8 autoclave cycles.
Ready in 5 hours.

REF 19723101

A32028

Promo : 345,00



Multidem C27

Features a microfilter, capable of supplying two autoclaves with demineralised water simultaneously

Demineralised water supply
for your autoclave

REF 19723112

A35531

Promo : 363,00



NEW!

U-Son Ultrasonic cleaner

U-Son 6 litres
REF 19600005

51687
Promo : 999,00

U-Son 3 litres
REF 19600004

A51686
Promo : 659,00



EliSafe label printer

Can be connected directly to the W&H autoclave, no external computer required. The labels are only printed if the autoclave cycle has been completed successfully. Depending on your autoclave type, a connection kit may be required. One possible option is the barcode scanner, see photo on left-hand page.

REF 19721132, €189

Laser printer for
steri bags.

REF 19721109

A46215
Promo : 669,00



Seal²

Excellent basic sealer suitable for 7 roll sizes, from 5.5 up to 30 cm. No settings required, features transparent roller guide. The roll holder can optionally be wall-mounted.

REF 19722101

A33878
Promo : 4708,00



BePro soaking bath

including draining insert.
3-litre capacity.

REF 19500300

A51688
Promo : 54,25



SealVal²

This machine features a touchscreen through which settings can be easily saved or adjusted. The roll holder can also be wall-mounted to save space, and all data can be stored via a LAN connection. Suitable for 7 roll sizes, from 5.5 up to 30 cm.

Fast thermal sealer
for rolls up to 30 cm wide

REF 19722103

A51721
Promo : 4708,00



Assistina Twin

Découvrez
l'Assistina Twin
dans notre showroom !

Reprocess up to 360 instruments per hour!

- › Automatic lubrication and maintenance device
- › Reprocess instruments in just 10 seconds
- › Short waiting times thanks to the two-chamber system
- › Simple and ergonomic operating concept
- › The Assistina Twin Care Set can handle the reprocessing of no fewer than 2,800 instruments
- › Add and remove instruments even more easily and faster thanks to the Quick Connect adaptor solution

Supplied with
TWIN Care Set:
HEPA filter, cleaning fluid and oil

REF 07484000

A41745

Prijs : 122,00

NEW!
Now with optional
blow-through
function, for
bone-dry
instruments!

assistina  TWIN

Assistina Twin

Assistina TWIN MB-302

- + 1 Quick ISO adaptor
- + 1 Quick RM adaptor*

REF 30310003

A49286

Promo : 2185,00

Assistina TWIN MB-302

- + 2x Quick ISO adaptors

REF 30310001

A51721

Promo : 2185,00



Assistina ONE

**The new 1-door
reprocessing device from W&H.**

Superfast reprocessing time: The Assistina ONE can lubricate and clean a contra-angle handpiece or air turbine in just 18 seconds.

Includes:

ONE Care Set:

HEPA filter, cleaning fluid and oil

REF 08064970

A41745

Promo : 122,00



assistina^{ONE}

**2x the speed
of the previous
model!**

Assistina ONE

Assistina ONE MB-301

+ 1 Quick ISO adaptor

+ 1 Quick RM adaptor*

REF 30398003

A51721

Promo : 116,00

Assistina ONE MB-301

+ 1x Quick ISO adaptor

REF 30398001

A51721

Promo : 4708,00




NEW!

The new Thed and Thed+

Top-class thermal disinfectors of the highest quality

- › With active drying, your precious instruments come out of the machine dry
- › Fast cycles
- › Compact, large-capacity machine
- › Offers single or double load level
- › The chemicals stay in the machine
- › Numerous accessories available for all instruments
- › Easy operation and connectivity
- › Available as a tabletop or undercounter model
- › Features quick cycle for contra-angle handpieces
- › Large, practical display




Thed and Thed+ tabletop model, with one rack



Thed and Thed+ tabletop model, with two racks




Space for 15 contra-angle handpieces



Thed+ also features:


- + light bar with status indication
- + programme advice
- + recognition of the racks and any injector rail in place
- + internal lighting
- + automatic adjustment of the following for each load:
 - cycle duration
 - water volume
 - chemical usage
 - power consumption




**With
active
drying!**

Prix sur demande !

Thed and Thed+ undercounter model, with one rack



Thed and Thed+ undercounter model, with two racks



Thed tabletop model

REF 21050000

A51690

Thed undercounter model

REF 21051000

A51691

Thed+ tabletop model

REF 21055000

A51692

Thed+ undercounter model

REF 21056001

A51693

*excluding accessories and installation

Installation, servicing and maintenance by your dealer or certified W&H service partners. Request the brochure for specifications and all accessories!

thed

The tabletop models feature handy cassettes for cleaning agents



The undercounter models feature a drawer for cleaning agents






1. IMAGE ACQUISITION

SICAT
seethrough | studio


2. TREATMENT PLANNING



3. PROTOCOL PRESETTING



4. TREATMENT



6. DOCUMENTATION



5. MONITORING & DECISION SUPPORT

Original Osstell technology

ioDent
OsstellConnect

For further information about the products in this brochure, please contact your depot or dealer.

wh.com
+32 471268601
carine.bulte@wh.com



www.denta.be
03 450 93 20
vente@denta.be

These offers are valid for the Belgium & Luxembourg market only; prices are VAT excluded.
Valid from 01.02.2026 to 30.06.2026

The complete implantology workflow performed and documented using **W&H** equipment and software.