

#whdentalwerk



wh.com



# PROMOTIONS

FOR DENTA

**Valid from**  
01.02.2026 to  
30.06.2026

**NEW!**



**Implantmed Plus II** surgical motor  
with or without piezoelectric technology

www.denta.be  
03 450 93 20  
verkoop@denta.be



## Synea Fusion

The Synea Fusion is the world's most widely used contra-angle handpiece – quiet, strong and durable.

The air turbines in this series are equipped with ceramic bearings, a quadruple spray system, and are available with various head sizes.

**Safety first:** All Synea turbine handpieces are equipped with a technology that ensures that aerosols do not enter the instrument.

synea FUSION



The powerhouse of contra-angle handpieces!

# RESTORATION & PROSTHETICS

Synea Fusion Power Edition

• **WG-900 LT**

**Red contra-angle handpiece 1:4**  
with light, Ø 10.7 mm  
triple spray  
**36\*-month warranty**

REF 30459000

A50376

**Promo: 799,00**  
**Per 3 : 699,00 stuk**



- Extremely powerful and robust contra-angle handpiece
- For the removal of zirconium and attachments
- Optimal power to speed ratio thanks to 1:4 speed ratio
- Indispensable for difficult restorations  
25% longer lifespan



**WG-56 LT**

**Blue contra-angle handpiece 1:1**  
with light, Ø 9.5 mm  
single spray  
**36\*-month warranty**

REF 30009000  
A37655  
**Promo : 535,00**  
**Per 3 : 455,00 stuk**



**WG-66 LT**

**Green contra-angle handpiece 2:1**  
with light, Ø 9.5 mm  
single spray  
**36\*-month warranty**

REF 30011000  
A37654  
**Promo : 615,00**  
**Per 3 : 525,00 stuk**



**WG-99 LT**

**Red contra-angle handpiece 1:1**  
with light, Ø 9.5 mm  
Quattro spray  
**36\*-month warranty**

REF 30007000  
A37653  
**Promo : 795,00**  
**Per 3 : 615,00 stuk**



Good spray and lighting

**HG-43 A**

**Blue handpiece 1:1**  
no light  
**single spray**  
**24\*-month warranty**

REF 30013000  
A37656  
**Promo : 358,00,00**  
**Per 3 : 315,00 stuk**



**TG-98 L\*** for Roto Quick

**TG-98 LM\*\*** for Multiflex®  
**Air turbine** \*LED+, \*\*with light,  
29 W\*\*, Ø 10 mm,  
Quattro Spray

**24\*-month warranty**  
\*REF 30001000  
\*\*REF 30003000  
V37658  
**Promo : 635,00**  
**Per 3 : 525,00 stuk**



**TG-97 L\*** for Roto Quick

**TG-97 LM\*\*** for Multiflex®  
**Air turbine** \*LED+, \*\*with light,  
29 W\*\*, Ø 11,5 mm,  
Quattro Spray

**24\*-month warranty**  
\*REF 30004000  
\*\*REF 30006000  
V37657  
**Promo : 635,00**  
**Per 3 : 525,00 stuk**



Cast from a single block, no seams



**NEW**  
More value for  
no added  
cost!

Now with an anti-scratch surface!



\* when registered via wh.com

\*\* up to 29 Watt in accordance with ISO 14457

## Synea Vision

High-end contra-angle handpieces and air turbines; lightweight, quiet, with excellent grip, a slim head and a long warranty.

The air turbines in this series are equipped with ceramic bearings, a quintuple spray system, a scratch-resistant surface and are available with four different head sizes.

synea VISION

Ring LED+

Work without shadows!



### • WK-56 LT

**Blue contra-angle handpiece 1:1**  
with light, Ø 9.5 mm  
single spray

**36-month warranty**

REF 30025000

A37580

**Promo : 699,00**



### • WK-66 LT


**Green contra-angle handpiece 2:1**  
with light, Ø 9.5 mm  
single spray

**36-month warranty**

REF 30026000

A37596

**Promo : 739,00**



Synea Vision Power Edition

### • WK-900 LT

**Red contra-angle handpiece 1:4**  
with light, Ø 10.7 mm  
triple spray

**36-month warranty**

REF 30418000

A50377

**Promo : 949,00**



With a narrow head

### • WK-93 LT

**Red contra-angle handpiece 1:4.5**  
with light, Ø 9.1 mm  
quintuple spray  
ceramic bearings

**36-month warranty**

REF 30024000

A37549

**Promo : 878,00**



### • WK-99 LT

**Red contra-angle handpiece 1:5**  
with light, Ø 9.5 mm  
quintuple spray  
ceramic bearings

**36-month warranty**

REF 30023000

A37652

**Promo : 878,00**



**Extremely powerful and robust contra-angle handpiece**

**For the removal of zirconium and attachments**

**Optimal power to speed ratio thanks to 1:4 speed ratio**

**Indispensable for difficult restorations**

**\*25% longer lifespan**



Optimal turbine design for maximum power

**Very narrow head**

**TK-94 L\*** for Roto Quick  
**TK-94 LM\*\*** for Multiflex®

**Air turbine \*Ring LED+  
 \*\*with light**

14 W, Ø 9 mm  
 quintuple spray

**24-month warranty**

\*REF 30016000  
 \*\*REF 30031000



V38968  
**Promo : 899,00**

**TK-97 L\*** for Roto Quick  
**TK-97 LM\*\*** for Multiflex®

**Air turbine \*Ring LED+  
 \*\*with light**

21 W, Ø 10 mm  
 quintuple spray

**24-month warranty**

\*REF 30017000  
 \*\*REF 30018000



V38696  
**Promo : 899,00**

**TK-98 L\*** for Roto Quick  
**TK-98 LM\*\*** for Multiflex®

**Air turbine \*Ring LED+  
 \*\*with light**

36 W\*\*, Ø 11.5 mm

**24-month warranty**

\*REF 30019000  
 \*\*REF 30020000



V38970  
**Promo : 855,00**

**TK-100 L\*** for Roto Quick  
**TK-100 LM** for Multiflex®

**Air turbine \*Ring LED+  
 \*\*with light**

36 W\*\*, Ø 13 mm

**24-month warranty**

\*REF 30021000  
 \*\*REF 30022000



V38971  
**Promo : 899,00**

**Synea Vision Short series; 20% lighter and shorter, yet just as powerful as standard Vision instruments!**

**HK-43 LT**

**Blue handpiece 1:1,**  
 with light,  
 single spray

**36-month warranty**

REF 30028000

A37641  
**Promo : 629,00**



**WK-56 LT S**

**Blue contra-angle  
 handpiece 1:1**  
 with light, Ø 9.5 mm  
 single spray

**36-month warranty**

REF 30217000

A39134  
**Promo : 699,00**



**WK-99 LT S**

**Red contra-angle  
 handpiece 1:5**  
 with light, Ø 9.5 mm  
 quintuple spray

**36-month warranty**

REF 30215000

A50661  
**Promo : 878,00**



**WK-93 LT S**

**Red contra-angle  
 handpiece 1:4.5**  
 with light, Ø 9.1 mm  
 quintuple spray

**36-month warranty**

REF 30216000

V39175  
**Promo : 878,00**



**Lightweight  
 and  
 compact!**

**Synea Vision Short series; 20% lighter and shorter, yet just as powerful as standard Vision instruments!**

## Implantmed Plus SI-1023

The powerful **Implantmed Plus** from **W&H** guarantees a high level of safety during treatment for both you and your patients. With its intuitive user interface, colour touchscreen and glass surface, the Implantmed offers a truly user-friendly experience. The drive is ideal for surgical procedures with its automatic thread-cutting function, automatic torque control and adjustable parameters. For **piezoelectric surgery**, the **Piezomed Plus Module** can be fully integrated with the Implantmed Plus surgical motor, as shown in the photo.

Both devices can be operated with a single foot control (which can be wireless), and the instruments connected can be supplied with the same physiological saline solution.



Work according to protocol



with automatic thread-cutting function



**Piezomed Plus  
Module**



**Implantmed Plus**

**Package A**  
**Implantmed Classic SI923 + 2x WI-75 E/KM**

- › Implantmed Classic no light
  - › Wired foot control
  - › Surgical contra-angle handpiece  
 2 x WI-75 E/KM (20:1, green)
- REF 90000161  
 REF 10207550  
 REF 10207550



**Package Deal PD0445**

**Package B**  
**Implantmed Plus SI1023 + WS-75**

- › Implantmed Plus no light
  - › Wired foot control
  - › Surgical contra-angle handpiece WI-75 E/KM (20:1, green)
- REF 90000163  
 REF 10207550



**Package Deal PD0286**

**Piezomed Plus**  
**Module SA-435 M**

The Piezomed Plus Module can be connected to the Implantmed Plus surgical motor via plug and play. The most common oral surgical treatments can be performed with this combination.

**Features:**

- › tip recognition
  - › LED light
  - › handpiece with a 1.8-metre cable
- REF 90000350

A47782  
**Promo : 4350,00**



piezomed

**Package C**  
**Implantmed Plus SI1023 + WS-75 L**

- › Implantmed Plus with light
  - › Wired foot control
  - › Surgical contra-angle handpiece WS-75 L (20:1 green)
- REF 90000165  
 REF 30033000



**Package Deal 0056**

## Implantmed Plus II

The new Implantmed Plus II is the next evolution of the SI-1023, featuring the following key characteristics:

- › All motors are equipped with light
- › Large, 7 inch touchscreen
- › Wireless foot control at no extra cost
- › Easy connection to the Piezomed module Plus II:
  - No additional cables required
  - No additional cooling required
  - No additional foot control required
- › Easy connection to the W&H ioDent platform for software updates and patient documentation via built-in Wi-Fi
- › Option to scan all materials and generate a comprehensive implant report
- › Connection with Osstell ISQ measurement and reporting possible and straightforward

### Package D Implantmed Plus II SI-2100 + WS-75 L

- › Implantmed Plus II with light
  - › Wired foot control
  - › Surgical contra-angle handpiece  
WS-75 L (20:1, green)
- REF 90000418  
REF 30033000



Package Deal PD0442

### Package E Implantmed Plus II SI-2101 + WS-75 L

- › Implantmed Plus II with light
  - › Wireless foot control
  - › Surgical contra-angle handpiece  
WS-75 L (20:1, green)
- REF 90000419  
REF 30033000



Package Deal PD0443

### Package F Implantmed Plus II SI-2102 + WS-75 L

- › Implantmed Plus II with light
  - › Scan edition
  - › Wireless foot control
  - › Surgical contra-angle handpiece  
WS-75 L (20:1, green)
- REF 90000420  
REF 30033000



Package Deal PD0444



**NEW!**

The new  
**Implantmed Plus II**  
 surgical motor;  
 can be supplemented  
 with piezoelectric  
 technology if desired!



implantmed <sup>PLUS</sup> II  
 piezomed <sup>PLUS</sup> II  
 module

**Flightcase**

› For the Implantmed Plus,  
 Piezomed Module and  
 accessories (in combination  
 with the equipment)

REF 08237370

A52050  
**Promo : 499,00**



**Surgical  
 motor only**



**Surgical motor  
 with piezoelectric  
 module**



# Surgical handpieces & contra-angle handpieces

The instruments below boast lights when connected to an EM19 LC motor.  
When connected to an EM19 motor there is no light source.

## • WS-75 L

With light.  
Implant site preparation  
and implant placing

20:1, contra-angle handpiece  
for surgical drills (green)  
Ø 2.35 mm

REF 30033000

A40683

Promo : 948,00



## • WS-56 L

With light.  
Osteotomy, sequestromy,  
germectomy, hemisection

1:1, contra-angle handpiece  
for surgical drills (blue)  
Ø 2.35 mm

REF 30040000

A40698

Promo : 999,00



## • WS-91 L

With light.  
Specially designed to  
improve access to the  
premolar and molar area

1:2.7, contra-angle handpiece  
for surgical drills (red)  
FG, with a 45° head tilt

REF 30049000

A34809

Promo : 999,00



## • WS-75 LG

Boasts its own lighting  
Light output depends on speed.  
Applications such as WS-75 L

20:1, contra-angle handpiece  
for surgical drills (green)  
Ø 2.35 mm

REF 30032000

A34809

Promo : 1089,00



Can be  
disassembled  
without tools  
for optimal  
hygiene

## • WS-92 L

With light.  
Osteotomy, sequestromy,  
germectomy, hemisection

1:2.7, contra-angle handpiece  
for surgical drills (red)  
FG

REF 30043000

A41640

Promo : 999,00



## • S-11 L

With light.  
Apical surgery, osteotomy,  
bone modelling and  
wisdom tooth surgery

1:1, handpiece for  
surgical drills (blue)  
Ø 2.35 mm

REF 30058000

A41366

Promo : 955,00



## • S-11 LG

Boasts its own lighting  
Light output depends on speed.  
Applications such as S-11 L

1:1, contra-angle handpiece  
for surgical drills (blue)  
Ø 2.35 mm

REF 30057000

A39292

Promo : 4350,00



See wh.com for more models and versions

• **WS-56 LG**

Boasts its own lighting.  
Light output depends on speed.  
Application as WS-56 L

1:1, contra-angle handpiece  
for surgical drills (blue)  
Ø 2.35 mm

REF 30039000

A51712  
**Promo : 1119,00**



• **WS-91 LG**

Boasts its own lighting.  
Light output depends on speed.  
Application as WS-91 L

1:2.7, contra-angle handpiece  
for surgical drills (red)  
FG, with a 45° head tilt

REF 30048000

A39226  
**Promo : 1119,00**



• **WS-92 LG**

Boasts its own lighting.  
Light output depends on speed.  
Application as WS-92 L

1:2.7, contra-angle handpiece  
for surgical drills (red)  
FG

REF 30042000

A48939  
**Promo : 1119,00**



• **WS-56**

No light.  
Applications such as WS-56 L

1:1, contra-angle handpiece for  
surgical drills (blue), Ø 2.35 mm

REF 30038000

A28896  
**Promo : 798,00**



• **WS-91**

No light.  
Applications such as WS-91 L

1:2.7, contra-angle handpiece  
for surgical drills (red)  
FG, with a 45° head tilt

REF 30047000

A49284  
**Promo : 898,00**



• **S-11**

No light.  
Applications such as  
S-11 L

1:1, contra-angle handpiece  
for surgical drills (blue)  
Ø 2.35 mm

REF 30056000

A01882  
**Promo : 698,00**



• **WS-75**

No light.  
Applications such as  
WS-75 L

20:1, contra-angle handpiece  
for surgical drills (green)  
Ø 2.35 mm

REF 30065000

A01881  
**Promo : 898,00**



• **WI-75 E/KM**

No light.  
Applications such as WS-75 L  
Cannot be disassembled.

20:1, contra-angle handpiece  
for surgical drills (green)  
Ø 2.35 mm

REF 10207550

A28895  
**Promo : 725,00**



• **WS-92**

No light.  
Applications such as WS-92 L

1:2.7, contra-angle handpiece  
for surgical drills (red), FG

REF 30041000

A31341  
**Promo : 898,00**



See wh.com for more models and versions



**Instrument Kit Bone**

Surgery Tips included:  
B1, B2R, B2L, B3, B4, B5

REF 07173100

A40330  
**Promo : 999,00**



**Instrument Kit Sinus Special**

For sinus lifts.  
Tips included:  
S1, S2, S3, S4, S5

REF 07173200

A40331  
**Promo : 705,00**



**Instrument Kit Endodontics**

Tips included:  
R1D, R2RD, R2LD, R3D,  
R4RD, R4LD

REF 07173300

A51716  
**Promo : 790,00**



**Instrument Kit Periodontal**

Tips included:  
P1, P2RD, P2LD

REF 07173400

A51717  
**Promo : 444,00**



**Instrument Kit Extraction**

Tips included:  
EX1, EX2

REF 07173500

A40339  
**Promo : 311,00**



**Instrument Kit Implant Crestal P**

Posterior. Tips included:  
I1, I2P, I3P, I4P,  
Z25P, Z35P

REF 07997110

A51718  
**Promo : 1375,00**



**Instrument Kit Implant Crestal A**

Anterior. Tips included:  
I1, I2A, I3A, I4A,  
Z25P, Z35P

REF 07997090

A5171  
**Promo : 1375,00**

**Fast  
Precise  
Powerful**



**Instrument B6 & B7**, surgical saws for, among other things, apical resection and the removal of root remnants or alveolar bone

B6 REF 07022500

A440330  
**Promo : 999,00**

B7 REF 07022600

A40331  
**Promo : 705,00**

## Piezomed Plus module

W&H piezoelectric technology brings all the benefits of modern ultrasonic technology within reach: high-frequency vibrations facilitate extremely precise bone surgery. What is more, thanks to the cavitation effect you can maintain a virtually blood-free surgical site, while also minimising tissue damage.

- Sinus lift
- Bone harvesting
- Apical resection
- Retrograde endodontics
- Extraction
- Periodontal (sterile)
- Etc.

Visit [video.wh.com](http://video.wh.com) for all applications

piezomed

## Courses on the piezoelectric technique in dental surgery:


Internationally renowned dental implantologist Irfan Abas will teach you the numerous possibilities of the piezoelectric technique in modern dental surgery.

[www.abas-academy.com](http://www.abas-academy.com)

 **ABAS**  
implant academy



The Dutch training institute Dental Best Practice also offers an intensive, hands-on piezoelectric course as part of their programme. For more information, visit [dentalbestpractice.nl](http://dentalbestpractice.nl)



# Surgical saw handpieces

**Fast, powerful and durable**

The saw handpieces are specially designed for conversion osteotomies and bone harvesting. A consistently high cutting performance is guaranteed even after many uses at maximum speed.

The saw blades are available in different sizes and are supplied sterile.

## S-8 R

Saw handpieces with reciprocating motion

REF 10100801

A33125  
**Promo : 1299,00**



## S-8 O

Saw handpieces with oscillating motion

REF 10100802

A29430  
**Promo : 1549,00**



## S-8 S

Saw handpieces with sagittal motion

REF 10100803

A49283  
**Promo : 1499,00**



**Easy to disassemble** and clean really effectively, ensuring a long service life and optimal patient safety

## Osstell Beacon

### Implant Stability Quotient (ISQ)

- › Monitoring of osseointegration
- › Makes an earlier prosthetic process possible
- › Improved predictability of treatment for patients with risk factors
- › Easy communication with patients through visualisation of the growing-in phase
- › Maximum flexibility thanks to wireless technology
- › The Osstell Beacon can be integrated with the latest generation of W&H Implantmed



Work in a validated and documented way; ideal for referral practices



## Osstell Beacon

REF 103000

A43554

Promo : 2129,00



# IMPLANTOLOGY



## Osstell Classic

A51729

Promo : 1699,00

REF 08086148

The Osstell Classic has the same functionalities as the Osstell Beacon, but without the option to use the Osstell Connect Portal, and it cannot be connected to W&H Implantmed.

### Includes free use of the Osstell Connect portal:

- › build up the implant history
- › compare implants
- › compare patients
- › compare pathways
- › compare problems

SmartPegs are necessary in order to measure osseointegration. These are available from your dealer, see [wh.com](http://wh.com) for an overview.



# Proxeo Twist

Top quality handpieces and contra-angle handpieces for prophylaxis

For your treatment unit



More space with the **LatchShort attachment**

## WP-66 W

**Prophy 4:1 contra-angle handpiece**  
for W&H LatchShort prophy cups

REF 30189000

A42771

**Promo : 345,00**  
**Per 3 : 259,00 stuk**



## WP-64 MU

**Prophy 4:1 contra-angle handpiece**  
with snap-on adapter and key  
for all standard screw-in prophy cups and brushes

REF 10216443

A28731

**Promo : 345,00**  
**Per 3 : 259,00 stuk**



## WP-64 M

**Prophy 4:1 contra-angle handpiece**

for System Young Screw-in prophy cups and brushes.  
With key.

REF 10216403

A02406

**Promo : 345,00**  
**Per 3 : 259,00 stuk**



## WP-66 M

**Prophy 4:1 contra-angle handpiece**

for all prophy cups, brushes and snap-on adaptors with a 2.35 mm shaft

REF 30190000

A39034

**Promo : 345,00**  
**Per 3 : 259,00 stuk**



## HP-44 M

**Prophy 4:1 handpiece**  
for Doriot system disposable prophy cups

REF 1014403

A02320

**Promo : 279,00**  
**Per 3 : 219,00 stuk**





# Proxeo Twist

## Cordless Polishing System



### Stand-alone polishing contra-angle handpiece

- › Work flexibly without cables
- › Easily regulate speed with wireless foot control
- › Perfect control of speed for precision cleaning and efficient polishing
- › Three different attachments
- › Minimal use of aerosols thanks to rotary polishing

**Handy for working on location!**

V46570  
**Promo : 999,00**

### PL-40 HW + PL-66 W

Proxeo Twist Cordless  
 + LatchShort cups

REF 90000346



### PL-40 HW + PL-64

Proxeo Twist Cordless  
 + screw-in for  
**YOUNG cups**

REF 90000347



### PL-40 HW + PL-44

Proxeo Twist Cordless  
 + Doriot attachment for  
 disposable contra-angle  
 handpieces with cup

REF 90000348



LatchShort cups

Screw-in for YOUNG cups

Doriot system with cup



Wireless foot control



**Compact,  
stackable,  
and fast!**

**Ontdek de Lisa Mini  
in onze showroom !**



# Lisa Mini

## Compact, fast and easy to operate!

Just like the renowned Lisa large autoclave, this smaller version B autoclave is fast and features an easy-to-use touchscreen. Available with a 5- or 3-litre chamber.

Key features:

- › Fast contra-angle handpiece cycle in just 10 minutes (5-litre chamber)
- › Standard cycle time of 22 minutes (5-litre chamber)
- › Data can be logged via USB, LAN or the W&H Steri App
- › Can be connected to a demineralised water supply
- › A label printer can also be connected
- › Available in 2 or 3-tray format (5-litre chamber)

### 5 litres

REF 19935100

A51721

**Promo : 4708,00**

### 5 litres\*

REF 19945100

A52653

**Promo : 4247,00**

### 3 litres\*

REF 19944100

A50654

**Promo : 3999,00**

\*requires compressed air connection



**Handy as an  
extra autoclave,  
when you have  
limited space  
or for small  
volumes!**



Other items can be placed on top of the unit to make the most of your sterilization space. Comes with 2 trays as standard.



# Lisa

**Ontdek de Lisa  
in onze showroom !**

## High-end type B autoclave with extremely short cycle times.

Standard cycle time of 30 minutes, fast contra-angle handpiece cycle in 13 minutes. Features Eco Dry+ technology that adjusts the drying time to the weight of the load. Boasts an LED display and a touchscreen. Available with a 17- or 22-litre chamber. Data can be logged via USB, LAN or the W&H Steri app. The sterilization process can be tracked and documented down to the individual instrument or instrument kit via EliTrace. No extra software or computer needed. Thanks to smart technology, the Lisa warns you when your load is overheated, even if the door is open – no more burnt fingers!

**Fast  
and  
quiet!**

**17 litres**

REF 19931180

A42621

**Promo : 5990,00**

**22 litres**

REF 19932180

A42687

**Promo : 6390,00**



**The Lisa from  
W&H is the  
most widely used  
autoclave in the  
world, and for  
good reason!**





# Lara, Lara XL and Lara XXL

**The fastest 28-litre class B autoclave on the market!**

- › Tried and trusted technology
- › Large capacity
- › Various upgrades available
- › Easy to operate

The Lara is also available in 17- and 22-litre versions:

## 17 litres

REF 19941120

A42622  
**Promo : 5446,00**

## 22 litres

REF 19942120

A42686  
**Promo : 5935,00**

## Lara XXL 38 litres

REF 19936110

A51724  
**Promo : 4708,00**

## Lara XL 28 litres

REF 19933110

A48819  
**Promo : 6990,00**

**38-litre capacity for up to 12 kg per cycle**



**Fits 60 cm trays**



**Highly suitable for surgical cassettes thanks to its large size**



# Lina

**Type B autoclave with best-in-class cycle times.**

- › Optimal sterilization programmes
- › The Eco B sterilization programme processes 0.5 kg in 25 minutes
- › Lina can reprocess loads of up to 6 kg to satisfy all the requirements of modern dental clinics
- › The B Universal 121°C cycle even allows you to reprocess delicate objects, including porous objects such as surgical clothing



### 17 litres

REF 19951131

A51722

**Promo : 5299,00**

### 22 litres

REF 19952131

A51723

**Promo : 5599,00**

**Excellent value for money!**



The **Activation Code** system allows you to activate additional features to meet your practice's changing needs and future requirements.

lina

### Helix test kit

**Complete set; consisting of the holder (the “snail shell”) and 250 test strips.**

Each test strip has a QR code that can be scanned and saved in our Steri app or in your network. Each code is linked to a sterilization cycle, which validates the process.



**With QR-coded strips**

REF T800207X

A34926

**Promo : 153,90**



### Q-dry

**Perfect solution for internal drying of contra-angle handpieces and air turbines.**

Using the Assistina adaptors and air pressure, you can easily dry instruments prior to sterilization.

See wh.com for how-to videos.



q-dry

REF 19723118

A51689

**Promo : 173,85**

### Water distiller

4-litre capacity, enough for 8 autoclave cycles. Ready in 5 hours.

REF 19723101

A32028

**Promo : 345,00**



### Multidem C27

**Features a microfilter, capable of supplying two autoclaves with demineralised water simultaneously**

Demineralised water supply for your autoclave

REF 19723112

A35531

**Promo : 363,00**



**NEW!**

**U-Son Ultrasonic cleaner**

U-Son 6 litres  
REF 19600005

51687  
**Promo : 999,00**

U-Son 3 litres  
REF 19600004

A51686  
**Promo : 659,00**



**EliSafe label printer**

Can be connected directly to the W&H autoclave, no external computer required. The labels are only printed if the autoclave cycle has been completed successfully. Depending on your autoclave type, a connection kit may be required. One possible option is the barcode scanner, see photo on left-hand page.

REF 19721132, €189

Laser printer for  
steri bags.

REF 19721109

A46215  
**Promo : 669,00**



**Ontdek de U-Son & Seal<sup>2</sup>  
in onze schowroom !**



**Seal<sup>2</sup>**

Excellent basic sealer suitable for 7 roll sizes, from 5.5 up to 30 cm. No settings required, features transparent roller guide. The roll holder can optionally be wall-mounted.

REF 19722101

A33878  
**Promo : 4708,00**



**BePro soaking bath**

including draining insert.  
3-litre capacity.

REF 19500300

A51688  
**Promo : 54,25**



**SealVal<sup>2</sup>**

This machine features a touchscreen through which settings can be easily saved or adjusted. The roll holder can also be wall-mounted to save space, and all data can be stored via a LAN connection. Suitable for 7 roll sizes, from 5.5 up to 30 cm.

Fast thermal sealer  
for rolls up to 30 cm wide

REF 19722103

A51721  
**Promo : 4708,00**



# Assistina Twin

Ontdek de  
Assistina Twin  
in onze showroom !

Reprocess up to 360 instruments per hour!

- › Automatic lubrication and maintenance device
- › Reprocess instruments in just 10 seconds
- › Short waiting times thanks to the two-chamber system
- › Simple and ergonomic operating concept
- › The Assistina Twin Care Set can handle the reprocessing of no fewer than 2,800 instruments
- › Add and remove instruments even more easily and faster thanks to the Quick Connect adaptor solution

Supplied with  
TWIN Care Set:  
HEPA filter, cleaning fluid and oil

REF 07484000

A41745  
Prijs : 122,00

**NEW!**  
Now with optional  
blow-through  
function, for  
bone-dry  
instruments!

assistina  TWIN

## Assistina Twin

### Assistina TWIN MB-302

- + 1 Quick ISO adaptor
- + 1 Quick RM adaptor\*

REF 30310003

A49286  
Promo : 2185,00

### Assistina TWIN MB-302

- + 2x Quick ISO adaptors

REF 30310001

A51721  
Promo : 2185,00



# Assistina ONE

**The new 1-door reprocessing device from W&H.**

Superfast reprocessing time: The Assistina ONE can lubricate and clean a contra-angle handpiece or air turbine in just 18 seconds.

**Includes:**

ONE Care Set:

HEPA filter, cleaning fluid and oil

REF 08064970



A41745  
**Promo : 122,00**

**2x the speed  
of the previous  
model!**

assistina ONE

## Assistina ONE

### Assistina ONE MB-301

- + 1 Quick ISO adaptor
- + 1 Quick RM adaptor\*

REF 30398003

A51721  
**Promo : 116,00**

### Assistina ONE MB-301

- + 1x Quick ISO adaptor

REF 30398001

A51721  
**Promo : 4708,00**





# The new Thed and Thed+

## Top-class thermal disinfectors of the highest quality

- › With active drying, your precious instruments come out of the machine dry
- › Fast cycles
- › Compact, large-capacity machine
- › Offers single or double load level
- › The chemicals stay in the machine
- › Numerous accessories available for all instruments
- › Easy operation and connectivity
- › Available as a tabletop or undercounter model
- › Features quick cycle for contra-angle handpieces
- › Large, practical display




### Thed and Thed+ tabletop model, with one rack



### Thed and Thed+ tabletop model, with two racks




### Space for 15 contra-angle handpieces



### Thed+ also features:


- + light bar with status indication
- + programme advice
- + recognition of the racks and any injector rail in place
- + internal lighting
- + automatic adjustment of the following for each load:
  - cycle duration
  - water volume
  - chemical usage
  - power consumption




**With active drying!**

**Prijs op aanvraag !**

**Thed and Thed+ undercounter model, with one rack**



**Thed and Thed+ undercounter model, with two racks**



**Thed tabletop model**

REF 21050000

A51690

**Thed undercounter model**

REF 21051000

A51691

**Thed+ tabletop model**

REF 21055000

A51692

**Thed+ undercounter model**

REF 21056001

A51693

\*excluding accessories and installation

Installation, servicing and maintenance by your dealer or certified W&H service partners. Request the brochure for specifications and all accessories!

thed

**The tabletop models feature handy cassettes for cleaning agents**



**The undercounter models feature a drawer for cleaning agents**






### 1. IMAGE ACQUISITION



SICAT  
seethrough | studio

### 2. TREATMENT PLANNING



ioDent  
OsstellConnect

### 6. DOCUMENTATION



### 5. MONITORING & DECISION SUPPORT

Original Osstell technology



### 4. TREATMENT



### 3. PROTOCOL PRESETTING

The complete implantology workflow performed and documented using **W&H** equipment and software.

For further information about the products in this brochure, please contact your depot or dealer.

**wh.com**  
+32 471268601  
carine.bulte@wh.com



**www.denta.be**  
**03 450 93 20**  
**verkoop@denta.be**

These offers are valid for the Belgium & Luxembourg market only; prices are VAT excluded.  
**Valid from 01.02.2026 to 30.06.2026**