

PROMOTIONAL OFFERS*

valid from 31st August till 18th of December 2026

Only for participating dealers BeNeLux.

All prices are in €, exclusive VAT.



*BeNeLux area



Menu

OPUS	4
C20	12
C50	14
Newtron	17
Piezotome Cube	21
Implant Center 2	22
Miniled	24
Riskontrol	28
Mouth mirrors	30
Expasyl	31
X-mind Dc	32
X-mind Unity	33
Sopix2	34
X-mind Optima 3D	36
X-mind Prime2	40
PSPIX2	43
Contact	47

opus YOUR WAY

€ 7.968

€ 880

Opus cart



reddot winner 2026
medical devices and technology design

acteon | opus



OPUS®

Testimonials

Those who have tried it speak about it best

The device will completely change the way I treat my patients. The unit will streamline appointments and make it easier to switch between supra and subgingival air polishing. With everything available on one unit it will enhance the patients' hygienist experience.

DH Sam Mc Hale (UK)

The device is what I've been waiting for... It's a good working workflow!

**Dr. Frederik-Christian
Hollay (Germany)**

It's versatile, easy to use, easy to maintain, good looking, more efficient, less invasive, delicate, human-made with personal input.

**Dr. Tea Natsvlshvili
(Georgia)**

OPUS® has raised the bar, keeping the smooth reliable scaler and making a duo unit is a game-changer. OPUS® is not just a machine but a partner in delivering exceptional patient care.

**DH Charlotte Benson
(UK)**

It's a very well designed machine that helps giving patients the best quality treatment with high efficiency and pain free workflow.

**DH Kassem Sindian
(Belgium)**

The air polishing is a nice experience for a patient.

**DH Carlotta Becciani
(Italy)**

OPUS is my best friend and assistant during patient treatments. It is easy to maintain, very user-friendly, and a pleasure to work with. My patients are also very happy with the results, which confirms its effectiveness for me.

**DH Sylvi Boer
(Netherlands)**



From simple prophylaxis to complex treatments, OPUS® adapts to every procedure and patient condition. **Its versatility empowers you to work your way enhancing your skills, elevating standards, and meeting patient expectations.**

For every **indication**,

min

PROSTHETICS


min Intervention
MAX Precision

Restoration finishing
Minimal cavity preparation
& micro-abrasion
Prosthetic debonding

PERIODONTICS

min SRP
MAX Preservation

Subgingival debridement
Implant care
Peri-implantitis treatment



CHOOSE YOUR WAY,

We empower clinicians with evidence-based tools to create personalized treatments, aligned with EFP S3-level guidelines⁽¹⁾ for stage-based periodontal care.

1. Sanz M, et al. Treatment of stage I-III periodontitis-The EFP S3 level clinical practice guideline. *J Clin Periodontol*. 2020;47 Suppl 22(Suppl 22):4-60.



With OPUS®, you always find a minMAX solution: **minimal invasiveness and constraints, maximum performance and comfort** achieved your way.

one minMAX solution

MAX



PROPHYLAXIS

Stain & biofilm removal
PMPR
Supragingival scaling



min Effort
MAX Comfort

ENDODONTICS

Canal access & preparation
Irrigant activation
Retrograde preparation
Obturation
Retreatment



min Uncertainty
MAX Control

MAKE YOUR CARE UNIQUE

Our approach blends versatility and scientific rigor to support minimally invasive, effective therapies tailored to each patient.

NOW IT'S TIME

Smart Features, Effortless Control

Original Air polisher & Ultrasonic Scaling

Grab & Go Powder System

No handling – effortless bottle insertion and removal, with automatic pressurization triggered by pressing the foot pedal

No refilling – pre-filled bottles eliminate the need for manual powder transfer

No cleaning – sealed system prevents contamination and removes the need for cleaning

ColorDrive Control:

Intuitive color coding for instant treatment:




Ultrasonic mode:

Automatic power setting according to the selected indication



Airpolishing mode:

Automatic power and indication settings based on the inserted cleaning powder bottle




TO OPUS[®] YOUR WAY

Designed to intuitively make it yours



Developed with dental professionals and backed by decades of expertise, OPUS[®] adapts to you, not the other way around. From its ergonomic interface to intelligent automatic adjustments, every detail simplifies your workflow, making it smoother and truly your own.




Smart Handpiece Recognition

Automatically detects the handpiece and selects the appropriate mode.

Responsive Touchscreen Interface


Works perfectly even with wet gloves.

No hidden menu — every essential function is visible at a glance.



Smart Foot Pedal

Four intuitive controls for total freedom and smooth transitions during procedures.




Ergonomic Cart*

Compact, mobile, and designed to fit seamlessly into any workspace.

*Optional accessory

POWERED BY OPUS®

Exclusive OPUS® Ultrasonic features engineered to prevent vibration loss



Strict controlled amplitude of vibration

Using high-frequency technology and linear vibrations, OPUS® technology allows a **constant vibration amplitude** and a precise area of action, **regardless of the chosen power level** or the clinical conditions encountered.

It effectively removes calculus while preserving the integrity of dental tissues. Newtron® technology's controlled amplitude ensures:

- Efficient intervention
- Minimally invasive treatment
- Reduced risk of microtrauma




Controlled water flow for optimal cavitation

Precise irrigation control guarantees optimal cavitation at the end of the tip, maximizing biofilm degradation and improving disinfection of the treated area.^{1,2,3}

The synergy between ultrasonics and irrigation enhances **antibacterial action**, essential for effective patient care.


1. Vyas N, et al. High Speed Imaging of Cavitation around Dental Ultrasonic Scaler Tips. *PLoS One*. 2016;11(3):e0149804.
2. Felver B, et al. Cavitation occurrence around ultrasonic dental scalers. *Ultrason Sonochem*. 2009;16(5):692-7.
3. Gagnot G, Pobleto M-G. The proper use of ultrasonic devices: control of vibration. *Revue d'Odonto. Stomatologie* 2004;33:85-9.

Exclusive clogging control system engineered to minimize clogging from the bottle to the nozzle



Top-Aspiration Powder System

Unlike conventional systems that draw powder from the bottom, causing staging and agglomeration, OPUS® aspirates from the top, **creating a vortex** inside the powder bottle. This vortex ensures consistent powder flow into the pipe, minimizes waste, and prevents clogging.



2 Way delivery system

SILKY® air and powder pipes remain separate from the water channel until the top of the nozzle, **reducing the risk of clogging**.

C20

Easy to use. Easy to master.

More than just a tool, The C20 enhances the diagnostic process by creating an engaging and interactive experience for both the dental team and their patients.

By streamlining communication and boosting chairside efficiency, the C20 helps transform practices into high-performing, patient focused environments.



Join the game
OF MODERN
DENTISTRY



Level up
YOUR PRACTICE

START

Press start
AND GO

MADE IN
EUROPE



~~€1.499~~
€1.399

Ref. with AIS: S_C20_0002
Ref. without AIS: S_C20_0001

Press start, And Go.

The C20 offers a **plug-and-play experience**, making it **straightforward** to install in any dental office. Its **intuitive design** ensures that users can start capturing and managing images immediately, enhancing workflow without any disruption.



Plug & Play

Ready to use with
pre-installed driver.



Universal connectors

USB 2 & universal
socket.




Customize settings

A wide range of settings
to suit your needs.




Windows® & Mac® compatible




C50

Daylight



Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

Daylight*



High contrasts emphasizing anatomic fissures and details of the dentine structure / depth.

- Full HD: 1920x1080
- Ergonomic design with small tip (easy to use with small mouth/children)
- Autofocus
- Diagnostic feature: dental plaque and caries detection (C50 Fluo)

~~€3.241,41~~

€2.999

(daylight)

Ref. S_C50_0002 DAY LIGHT

~~€4.666,05~~

€3.999

(fluorescence technique)

Ref. S_C50_0001 FLUO




4 Modes

Daylight



Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

Perio



Enhancement of **dental plaque and gingival inflammation** thanks to chromatic amplification.

Daylight*



High contrasts emphasizing anatomic fissures and details of the dentine structure/ depth.

Cario



Caries detection by **fluorescence technique** by LEDs.

5 Visions


Macro



Intraoral



Portrait



Tooth



Smilew



Exclusive and patented Color Coding System™

Power selection is accurately determined with the Color Coding System™, which intuitively associates each tip color ring with one of the 4 power ranges



newtron

p5XS

1

Ref. F00246




1S

Ref. F00245



H3

Ref. F00369



P5XS b-LED incl. with the following tips:
Ref. F62100

~~€2.857,89~~

€ 2.140,99



N°4
Ø22


Ref. MP3622PH

NOW INCLUDED

PURE REFLECT®
with a value
of €39,99

Ref. MP3622PH

Simple and intuitive settings with the Color Coding System™



Newtron® handpiece



ref. F12281 - Autoclavable

~~€651,37~~

€345,78

Newtron® LED handpiece



ref. F12609

~~€901,61~~

€504,90

Newtron® Slim B. LED handpiece



- with blue LED ring: ref. F12900
- with white LED ring: ref. F12905

~~€875,34~~

€559,98

hygiene

Hygiene Kit



ref. F00935

~~€529,14~~

€473,99

periodontics

Perio Kit



ref. F00936

~~€670,97~~

€599,99

implantprotect

pure titanium

Implant Protect Kit



ref. F02120

~~€735,54~~

€682,38

endosuccess

retreatment

Endo Success Retreatment Kit





ref. F00737

~~€702,36~~

€632,40

- Excellent view, full LED ring
- Avoid necrosis up to 120ml irrigation with peristaltic pump
- Improve durability of the handpiece, external irrigation
- Maximum hygiene, sterile single use or autoclavable lines
- Work don't play, use of big bottles (1l) possible
- Large choice of clinical application, 38 tips
- Tips II are compatible with IC2, Piezotome Cube, Piezotome2, Piezotome solo




piezotome cube

**Dynamic power adaptation for
safe bone surgeries**

~~€7.799,04~~
€5.762,99

Ref. F50100

Delivered with a multifunction footswitch,
a Piezotome® CUBE LED handpiece,
1 handpiece holder, an Essential II kit
(6 tips BS1S, BS4, SL1, SL2, SL3 et LC2),
a Piezotome® dynamometric wrench,
1 irrigation pole, 2 autoclavable irrigation
lines, 5 clips and 15 perforators.



Intraoral surgery ultrasonic control console.
Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

implant center 2

3 modes:

- Piezotome® (bone surgery).
- Newtron® (perio, endo,...).
- I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.

**The alliance of technologies
for powerful and secure
bone surgeries**



~~€14.901,47~~
€9.790,99

Ultrasonic control console for conventional dentistry and intraoral dentistry and motor control console for dental implantology
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France



Bone Surgery™ II Kit*

Ref. F87509 - Tips BS1S II, BS2L II, BS2R II, BS4 II, BS5 II, BS6 II.

~~€1.432,35~~ **€1.111,99**



Sinus Lift II Kit*

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

~~€1.077,92~~ **€1.009,99**



Intralift™ II Kit*

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

~~€1.367,71~~ **€1.111,99**



Extraction™ II Kit*

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

~~€1.087,60~~ **€1.009,99**

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

Accessoires

Promotion valable only for ImplantCenter2

Newtron® LED handpiece



Ref. F12609

For Piezotome® 2 and ImplantCenter™ 2

Class IIa medical devices (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France

Newtron® LED cord



Ref. F12830

For Piezotome® 2 and ImplantCenter™ 2

€2.022,97

€1.489

Universal irrigation lines Sterile x10



€184,07

€147,99

as from 4 pieces

Ref. F59905 - Single-use with integrated cassette and flow control system.

Class Is medical device (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France

PIEZOTOME® LED handpiece + cord



Ref. F12815

Class IIa medical devices (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France

~~€1.635,16~~

€1.324,99



Universal II Kit*

Ref. F87539 - Tips BS6, LC1, LC2, CE1, CE3 et PZ1.

~~€1.393,47~~ **€1.111,99**



Essential II Kit*

Ref. F87528 - Tips BS1S II, BS4 II, SL1 II, SL2 II, SL3 II, LC2 II.

~~€1.432,35~~ **€1.111,99**



Crest Splitting II Kit*

Ref. F87567 - Tips CS1 II, CS2 II, CS3 II, CS4 II, CS5 II, CS6 II.

~~€927,18~~ **€938,70**



Crown Extension II Kit*

Ref. F87554 - Tips BS6 II, CE1 II, CE2 II, CE3 II.

~~€737,70~~ **€630,99**

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

miniled ortho 2

~~€1.486,25~~

€1.111,99

Ref. F05220

- **Customizable signal**

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

- **Power**

3000mW/cm²

(Ø5,5mm light guide).

- **Long lasting battery**

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

- **Ortho friendly**

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

- **Efficiency & safety**

Cure any type of material used orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.

Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from
few brackets to a whole quadrant



Simply Powerful

Outperforms in
power, efficiency and speed

~~€1.435,51~~
€1.111,99

Ref. F05214



miniled supercharged

• Visibility

Luminous screen with customizable modes and times, count curing time and load indicator.




• Power

2000mW/cm²
(Ø7,5mm light guide).

3000mW/cm²
(Ø5,5mm light guide - optional).

• Optimal result


3 modes to fit all clinical needs.



FAST MODE

The choice of rapidity


Full power during 3, 4, 5 or 10 sec, according to the composite thickness.



PULSE MODE

Limits temperature rise

5 or 10 flashes of 1 sec for a polymerization in security.



RAMPING MODE

Avoids composites retraction

Progressive power increase during 9 or 15 sec to ensure perfect sealing of restoration.

miniled standard

- Fast**

Only 6 to 12 sec to polymerize
2mm composite.

- Power**

1250mW/cm²
(Ø7,5mm light guide).


2000mW/cm²
(Ø5,5mm light guide - optional).

- Safety**

Base-charger includes a radiometer for light power
intensity check-up.

- Optimal result**


3 modes to fit all clinical needs.



FAST MODE

Fast and easy


Full power during 10 sec
(bip at 5 sec) on regular
depth.



PULSE MODE

Ally rapidity & safety

10 flashes of 1 sec
(bip after 5 flashes)
with times off to
preserve dental pulp.



RAMPING MODE

No stress cracking

Progressive power
increase during 20 sec
(bip every 5 sec) for
a complete and safe
curing.

~~€1.117,47~~

€907,99

Ref. F02507

Highly Adaptable

The alliance of power & safety





servotome

Delivered with a single footswitch,
a box with 10 electrodes,
a conductive bracelet.

Ref. F57210

~~€2.025,86~~
€1.732,99



Riskontrol® contributes to clinical protection




RISKONTROL® is also a **GREEN** solution

- Disposable tip
- Economic: simplify process to save time
- Recyclable plastic: environmentally friendly
- Cost-effective
- Flexible for a better use: respect of patient's anatomy
- Safety lock system: no risk of ejecting the tips



Sterilizable tips are not sufficient anymore

RISKONTROL® can be used for all procedures in any workflow



Easy to fit

1

It's installed.

2

It's fitted.

3

It's locked.

Discover our RISKONTROL® range

Adapt your tips to the taste
of your patients



art
by **riskontrol** ~~€62,76~~
€32,99



Ref. 201 700 - Mint



Ref. 201 701 - Mandarine



Ref. 201 702 - Liquorice



Ref. 201 703 - Anise



Ref. 201 704 - Black currant

Non-sterile disposable tips (x 250)

Add color

to your patients' care experience



~~€57,78~~
€32,90

classic
by **riskontrol**



Ref. 201 664 - White



Ref. 201 665 - Blue



Ref. 201 666 - Yellow



Ref. 201 667 - Green



Ref. 201 668 - Pink

Non-sterile disposable tips (x 250)

Protect entirely your air/water syringe

thanks to the plastic sheath



total protec
by **riskontrol**

~~€33,59~~
€29,99



Ref. 201 750



Ref. 201 751

Non-sterile disposable tips with plastic sheath (x 100)

Choose a totally sterile

workflow



~~€38,89~~
€30,50

sterile
by **riskontrol**



Ref. 201 770

Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile by RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - FRANCE.

Mouth mirrors

PURE REFLECT®

- Very bright image
- Resistant

For all purposes

RH

- A greyish reflection



TI

- A blueish reflection

FRONT SURFACE

Cross-sectional view





Image reflected directly on the surface of the glass, distortion-free.


NO DOUBLE IMAGE

BACK SURFACE

Cross-sectional view of an aluminium mirror



Reflecting surface under the glass. **DOUBLE IMAGE**



PURE REFLECT®


N°	Ref.	N.
4	MP3622PH	12
5	MP3624PH	12

~~€56,33~~

€39,99



A TRIPLE ACTION



Kaolin clay

&



Aluminum
chloride

1 ASTRINGENT EFFECT

helps to completely dry the sulcus with no blood or fluids.

2 OPEN THE SULCUS

effectively with no pain or damage to tissues thanks to the viscosity carefully calculated with our clinical expertise.

3 PRESERVE EPITHELIAL ATTACHMENT*

for a reduced risk of gingival recession and bone resorption.

*Canula must not be inserted into the sulcus

expasyl
exact



~~€144,33~~
€100,99

Paste for temporary sulcular opening
20 capsules (0,3g) Expasyl™ Exact - Ref. 261010



Paste for temporary sulcular opening
20 capsules (1g) Expasyl™ - Ref. 261030

expasyl

~~€250,09~~
€142,99



~~€111,70~~
€79,99

Straight - Ref. 261005
Bended - Ref. 261015

x-mind dc



The X-MIND® generators can accomodate any operatory configuration.

3 arms lengths are available : 41 cm (16.2"), 82.5 cm (32.5") and 110 cm (43.5") and can be:

- Top wall mounted
- Bottom wall mounted

Dental unit, ceiling and mobile version are also available.

With its adaptable wall plate, it is easy to replace an older generator with the X-MIND® without having to do additional drilling on the wall.



~~€4.821,80~~

€3.199

Available with **Remote Wireless Control**



No more wiring issues

Say goodbye to the complicated and expensive wiring work traditionally needed to connect the remote button in your practice.

Wireless Button Kit X-MIND DC

Ref. 5.15.00135

All X-Mind DC devices delivered from April 2021 onwards

€299



Cost effective

Save on installation expenses with the easy-to-integrate wireless solution.



Enhanced convenience

Enjoy the freedom and flexibility of a wireless connection, making your practice more efficient and user-friendly.



Wall version	Top mounting	Top mounting square cone	Bottom mounting	Bottom mounting square cone
0,40 m extensor	Ref. 5.92.23137	Ref. 5.92.23167	Ref. 5.92.23130	Ref. 5.92.23160
0,80 m extensor	Ref. 5.92.23138	Ref. 5.92.23168	Ref. 5.92.23131	Ref. 5.92.23161
1,10 m extensor	Ref. 5.92.23139	Ref. 5.92.23169	Ref. 5.92.23132	Ref. 5.92.23162

All references for the X-MIND® DC

Choose from a variety of cones to better fit to your clinical needs



Long cone (31 cm - 12")
Recommended for the parallel technique



Short cone (20 cm - 8")
Recommended for bisector technique



Rectangular Long cone (31 cm - 12")
Reduces the X-rays applied to the patient's skin by 50%.

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy

x-mind unity

~~€5.882,40~~
€3.999

A sharp and contrasted image

- A contrasted and defined image
- Optimized protection for the patient
- The generator and sensor association

Fixing from above, arm 0.80 m, short cone.

Intraoral X-ray

- Very high quality
 - Precise and flexible positioning
 - Fast selection from traditional to digital mode
- Minimum exposure time: 0,02s
 - Programmable timer
 - Top or bottom mounting

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy





sopix 2



A remarkable image quality, even for endodontics and the more complexes

SOPIX®2 provides a high resolution image by exploiting the latest technologies available on the market.



SOPIX®2 USB Sensor size 1
high definition with a sensor
holder, 10 sensor sheaths
Ref. S_802_0008

~~€7.058,90~~ **€3.999**

Also available in size 2
Ref. S_802_0009

~~€8.345,08~~ **€5.499**

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON® Group- France

Free yourself from constraints
with the embedded system and its
detachable sensor part.



SOPIX®2 Plug-in kit S1

Ref. SK60-1000

~~€6.655,03~~ **€5.999**


SOPIX®2 Plug-in kit S2

Ref. SK60-1001

~~€9.194,22~~ **€7.499**

Size 2
31x42mm

Size 1
25x39mm



sopix² plug-in

The latest technology with no
compromise on image quality with
more safety, freedom and comfort.

Resolution: >18lp/mm

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON® Group- France

Meet the X-MIND® optima 3D Family.

X-MIND® optima 2D/3D wall 100-240V

2D/3D PAN device 5x5 to 12x10 wall installation

VERSION CORE: W1203005.DRT · € 43.999,00

VERSION FLEX: W1203005 · € 46.999,00



X-MIND® optima 2D/3D ceph 100-240V

2D/3D PAN CEPH device 5x5 to 12x10 floor installation

VERSION CORE: W1203010.DRT · € 54.999,00

VERSION FLEX: W1203010 · € 57.999,00

X-MIND® optima 2D/3D floor 100-240V

2D/3D PAN device 5x5 to 12x10 floor installation

VERSION CORE: W1203005.FLDRT · € 45.999,00

VERSION FLEX: W1203005.FL · € 48.999,00



Order code	Option
W1000026	RealGUIDE® (lifetime license)
W1203107*	17x12 FOV *Available for FLEX version only
W0900153	DICOM option
W0900046	AIS standalone
W1200129	WeDiagnostix 1 year subscription



ACTEON CLUB
Welcome leaflet
+ Webkey



Workstation and monitor

x-mind optima3D



AIS



RealGUIDE®

Exceptional clarity.

At the heart of the X-MIND® optima 3D is **state-of-the-art detector technology**, which guarantees exceptional image clarity and outstanding anatomical coverage. **Its comprehensive range of Fields of View (FOV)**, from 5x5 to 17x12* cm (Ø x H), makes it suitable for a **wide range of dental procedures**, from general dentistry to orthodontics and maxillofacial surgery, while increasing tolerance for patient positioning and reducing the risk of marginal artifacts.



EXTRAORDINARY CLARITY

Get the most out of quality with the brand-new **ClearIQ** filter



OUTSTANDING VERSATILITY

From 5x5 to 17x12*, a tailored solution for every clinical need



TOP-LEVEL PERFORMANCE

Enhanced image depth with the new kV output range 60-90



MAXIMUM PATIENT COMFORT

Quick and easy acquisition workflow




AUTOMATIC POSITIONERS RECOGNITION

Maximum accuracy



ZERO FOOTPRINT


Maximum accessibility



* Option - The additional cost for the 17x12 version is €3,000


5x5 cm FOV

- Improves the **detection of root canals**, particularly in case of abnormal root canal shapes



8.5x5 cm FOV

- Useful for **implant surgery** and in **general dentistry** as a low dose exam as **limited to a single jaw**




8.5x9 cm FOV

- Useful for **complex implant surgery** as well as left or right TMJ diagnosis



12x10 cm FOV

- Suitable for most **oral surgery cases** including multiple implant surgery



17x12 cm FOV

New

- Represents the perfect solution for **complex orthognathic cases**, **zygomatic implant planning** and **orthodontic surgery**



Meet the 2D Prime Family

X-MIND® prime is available in three configurations: 2D wall mounted, 2D Ceph and 2D floor to suit all installation and clinical demands.



~~€27.786,47~~
€13.999

X-MIND® prime WALL 110-240V

Order code: W1200001

The ACTEON®'s smallest extraoral 2D imaging solution ever for the maximum **accessibility** in **any context**.*



~~€46.677,78~~
with one sensor: **€24.999**

~~€55.279,20~~
with two sensors: **€29.999**

X-MIND® prime CEPH 110-240V

Order code: W1200009

Cephalometric add-on for a full 2D system still encompassing the same philosophy of prime, **ideal for pediatric applications**.



~~€30.786,47~~
€14.999

X-MIND® prime FLOOR 100-240V

Order code: W1200001.FL


*Belgian market only (W1200001/W120001.PB)

Rely on ac omprehensive set of 2D exams


X-MIND® prime provides a comprehensive range of 2D exams designed to meet the diverse needs of general dentistry, including dental panoramic imaging, temporomandibular joint assessments and maxillary sinus evaluations.




Dental panoramic



Maxillary sinuses



Temporomandibular joint




Child panoramic



Half-panoramic
Left & right



Bitewing
Left & right



Frontal



PSPiX²

is AUTONOMOUS


The fully automated process free yourself from constraints

Just drop the plate inside the slit, the PSPiX² will automatically detect the size, scan the plate, optimize the image, and send it instantly to your imaging software.

In the meantime, it will automatically erase and eject the plate for an immediate reuse or switch to standby mode.

Scan
simply &
quickly

Fully automated process



—€8.525,20—

€6.299

Ref. S_702_0013

1

Auto-access:

The door only opens when the imaging plate is inserted.

2

Auto-detect:

PSPiX² automatically detects the size of the inserted imaging plate.

3

Auto-scan:

PSPiX² scans and optimizes the image before displaying it into the right workstation.

4

Auto-eject:

The imaging plate is automatically erased and ejected for an immediate reuse.

5

Eco-mode:

PSPiX² automatically switches to economic mode when not in use for an extended period of time.

Protective bag for ERLM*

For **PSPIX®2** and **PSPIX®**
(previous generation)

~~€85,82~~
€65,99



	PSPIX® (previous generation)	
	Box of 250	
Size 0		700338L
Size 1		700339L
Size 2		700340L
Size 3		700341L

ERLM plates*

For **PSPIX®2**

~~€131,16~~
€96,99



	PSPIX®2	
	Box of 2	
Size 0 (22 x 35mm)		990215
Size 1 (24 x 40mm)		990216
Size 2 (31 x 41mm)		990217
Size 3 (27 x 54mm)		990218

Protections for ERLM*

Protective cover for
PSPIX® and **PSPIX®2**
(previous generation)

~~€48,06~~
€26,99



Compatible with

PSPIX®

SOREDEX®
Digora

INSTRUMENTARIUM®
Express

DÜRR DENTAL®
Vistascan®

KAVO®

	PSPIX®2	PSPIX
PSPIX®	✓	✓
SOREDEX® Digora	✓	✓
INSTRUMENTARIUM® Express	✓	✓
DÜRR DENTAL® Vistascan®	✓	✓
KAVO®	✓	✓

	PSPIX® (previous generation)	
	Box of 250	
Size 0		700342L
Size 1		700343L
Size 2		700344L
Size 3		700345L

Protective bag for ERLM*

For **PSPIX®2**

~~€100,86~~
€70,99



	PSPIX®2	
	Box of 300	
Size 0		700546
Size 1		700547
Size 2		700548
Size 3		700549

- Protection against scratches, light and cross contamination
- Ergonomic and easy to use
- Economic solution

* Luminescent Radio Memory Display for digital phosphor plate systems

PSPIX® imaging plate - Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON® Group - France
Bag & cover: Class I medical devices - CE
Manufacturer: Sopro ACTEON® Group - France



The thin and flexible imaging plates, available in 4 sizes, will allow you to easily scan all morphologies...



Child, Adult Periapical	Adult, Bitewing Periapical	Child	Bitewing	Occlusal
Size 0 22x35mm	Size 1 24x40mm	Size 2 31x41mm	Size 3 27x54mm	Size 4 69x54mm

... and all clinical cases

Periodontics 	Endodontics 	Orthodontics 	Pedodontics 	Cariology
Bitewing 	Implantology 	Occlusal 		



CONTACTS



Area Sales Manager Equipment & Consumables BeLux

Sandra DELLA-FINA

Tel.: +32 (0) 494 116 499

sandra.della-fina@acteongroup.com

Area Sales manager Equipment & consumables the Netherlands

Roger VAN EIJSDEN

Tel.: +31 (0) 624 835 367

roger.vaneijsden@acteongroup.com

Area Sales Manager Imaging BeNeLux

Andy LUBERTI

Tel.: +31 (0) 613 948 984

andy.luberti@acteongroup.com

Imaging Specialist BeNeLux

Pierre DEKUSSCHE

Tel.: +32 494 111 253

pierre.dekussche@acteongroup.com

Administration Commerciale

Sabine PARIS

Belgium

Tel.: +33 (0) 556 349 332

Fax: +33 (0) 566 349 331

sabine.paris@acteongroup.com

Nirina INISAN

Luxembourg / The Netherlands

Tel.: +33 (0) 556 349 250

Fax: +33 (0) 566 349 331

nirina.inisan@acteongroup.com

ACTEON® CUSTOMER SERVICE

ACTEON's* customer service can be contacted at:

Imaging

supportbenelux@acteongroup.com

Equipment / Pharma

support@acteongroup.com

ACTEON®

17 avenue Gustave Eiffel

33700 MERIGNAC - FRANCE

Tel.: +33 (0) 556 340 607

Fax: +33 (0) 556 349 292

info@acteongroup.com

www.acteongroup.com

Non contractual document

Copyright© 2026 ACTEON®. All rights reserved. No information or part of this document may be reproduced or disclosed in any format without the prior authorisation of ACTEON®.

For professional dental use only.

Read the instructions in the notice carefully before use.

All information that is essential to ensure correct use of these devices is included in the user manual or the summary of the product characteristics available on www.acteongroup.com.

Creation date: August 2026

x-mind optima3D

EASY-JOIN QR Code

SCAN ME!



 **acteon club**