

PROMOTIONAL OFFERS*



valid from 31st August till 18th of December 2026

Only for participating dealers BeNeLux.
All prices are in €, exclusive VAT.



*BeNeLux area

opus YOUR WAY

A51491

Adviesprijs: 7968,00

A51492

Adviesprijs: 880,00

Opus cart

**Ontdek het OPUS-apparaat in
uw praktijk en ervaar innovatie
van dichtbij !**

**Vrijblijvende demonstratie via
marketing@denta.be**



Ontdek het OPUS-apparaat in onze toonzaal – ervaar innovatie van dichtbij !

opus® – Ultramoderne luchtpolijsteenheid en ultrasone scaler



Gecontroleerde reiniging en polijsten, alles in één

OPUS biedt een complete oplossing voor biofilmverwijdering, vlekverwijdering en tandsteencontrole.

Met precisiescaleren en polijsten wordt elke profylaxe- of parobehandeling gestroomlijnd, voor optimaal klinisch resultaat en meer werkcomfort.


Minimaal invasieve zorg, maximaal behoud

De nauwkeurig gecontroleerde ultrasone trillingen van OPUS® beschermen het harde en zachte weefsel.

Met meer dan 60 unieke tips biedt het systeem volledige ondersteuning voor al uw behandelbehoeften: van parodontologie tot endodontie.

Nieuwe air-polisher: ontworpen voor comfort

De SILKY-pen is een ultralicht air-polishing handstuk voor supra- en subgingivale toepassingen. Met vier verschillende mondstukken biedt het precisie, veelzijdigheid en gebruiksgemak, terwijl het slanke ontwerp een prettige werkervaring garandeert.



Opus-Cart
A51492

Prijs
880,00



Silky Supra

V51508 **Mint of lemon**

6 x 200 g

Prijs
109,00



9% korting

Package Deal PD0420

Koop de **Opus** aan **bundelprijs**
van **7570,00 euro** (waarde: 7968,00)
en koop **Opus-Cart** aan **50% korting**
of **440,00 euro** (waarde: 880,00).

OPUS[®]

Testimonials

Those who have tried it speak about it best

The device will completely change the way I treat my patients. The unit will streamline appointments and make it easier to switch between supra and subgingival air polishing. With everything available on one unit it will enhance the patients' hygienist experience.

DH Sam Mc Hale (UK)

The device is what I've been waiting for... It's a good working workflow!

**Dr. Frederik-Christian
Hollay (Germany)**

It's versatile, easy to use, easy to maintain, good looking, more efficient, less invasive, delicate, human-made with personal input.

**Dr. Tea Natsvlishvili
(Georgia)**

OPUS[®] has raised the bar, keeping the smooth reliable scaler and making a duo unit is a game-changer. OPUS[®] is not just a machine but a partner in delivering exceptional patient care.

**DH Charlotte Benson
(UK)**

It's a very well designed machine that helps giving patients the best quality treatment with high efficiency and pain free workflow.

**DH Kassem Sindian
(Belgium)**

The air polishing is a nice experience for a patient.


**DH Carlotta Becciani
(Italy)**

OPUS is my best friend and assistant during patient treatments. It is easy to maintain, very user-friendly, and a pleasure to work with. My patients are also very happy with the results, which confirms its effectiveness for me.

**DH Sylvi Boer
(Netherlands)**

POWERED BY OPUS®

Exclusive OPUS® Ultrasonic features engineered to prevent vibration loss




Strict controlled amplitude of vibration

Using high-frequency technology and linear vibrations, OPUS® technology allows a **constant vibration amplitude** and a precise area of action, **regardless of the chosen power level** or the clinical conditions encountered.

It effectively removes calculus while preserving the integrity of dental tissues. Newtron® technology's controlled amplitude ensures:

- Efficient intervention
- Minimally invasive treatment
- Reduced risk of microtrauma




Controlled water flow for optimal cavitation

Precise irrigation control guarantees optimal cavitation at the end of the tip, maximizing biofilm degradation and improving disinfection of the treated area.^{1,2,3}

The synergy between ultrasonics and irrigation enhances **antibacterial action**, essential for effective patient care.


1. Vyas N, et al. High Speed Imaging of Cavitation around Dental Ultrasonic Scaler Tips. *PLoS One*. 2016;11(3):e0149804.
2. Felver B, et al. Cavitation occurrence around ultrasonic dental scalers. *Ultrason Sonochem*. 2009;16(5):692-7.
3. Gagnot G, Pobleto M-G. The proper use of ultrasonic devices: control of vibration. *Revue d'Odonto. Stomatologie* 2004;33:85-9.

Exclusive clogging control system engineered to minimize clogging from the bottle to the nozzle



Top-Aspiration Powder System

Unlike conventional systems that draw powder from the bottom, causing staging and agglomeration, OPUS® aspirates from the top, **creating a vortex** inside the powder bottle. This vortex ensures consistent powder flow into the pipe, minimizes waste, and prevents clogging.



2 Way delivery system

SILKY® air and powder pipes remain separate from the water channel until the top of the nozzle, **reducing the risk of clogging**.

C20

Easy to use. Easy to master.

More than just a tool, The C20 enhances the diagnostic process by creating an engaging and interactive experience for both the dental team and their patients.

By streamlining communication and boosting chairside efficiency, the C20 helps transform practices into high-performing, patient focused environments.



Join the game
OF MODERN
DENTISTRY



Level up
YOUR PRACTICE



Press start
AND GO

MADE IN
EUROPE



A51758


~~€1.499~~

€1.399

Ref. with AIS: S_C20_0002
Ref. without AIS: S_C20_0001


C50

Daylight



Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

Daylight⁺



High contrasts emphasizing anatomic fissures and details of the dentine structure / depth.

- Full HD: 1920x1080
- Ergonomic design with small tip (easy to use with small mouth/children)
- Autofocus
- Diagnostic feature: dental plaque and caries detection (C50 Fluo)

A49262

~~€3.241,41~~

€2.999

(daylight)

Ref. S_C50_0002 DAY LIGHT

A49261

~~€4.666,05~~

€3.999


(fluorescence technique)

Ref. S_C50_0001 FLUO




4 Modes

Daylight




Transcript of the **natural colors** for a better interpretation of the dental tint and helps to identify pathologies and with post-treatment control.

Perio




Enhancement of **dental plaque and gingival inflammation** thanks to chromatic amplification.

Daylight⁺



High contrasts emphasizing anatomic fissures and details of the dentine structure/depth.


Cario



Caries **detection by fluorescence technique** by LEDs.

5 Visions


Macro



Intraoral



Portrait



Tooth



Smilew



newtron p5XS

1

Ref. F00246



1S

Ref. F00245



H3

Ref. F00369



P5XS b-LED incl. with the following tips:
Ref. F62100

47% korting

A38077

Prijs: 2857,89

Promo: 1490,00



N°4

Ø22

Ref. MP3622PH

48% korting

Package Deal PD0368

Koop de Acteon **Newtron P5 XS B.LED** aan promoprijs

van **1490,00 euro** (waarde: 2857,89 euro) en

koop de **Mondspiegel Pure Reflect #4** of **#5**

aan **1 euro** (waarde: 56,53 euro)

Newtron® handpiece



ref. F12281 - Autoclavable

A29690

Prijs: 651,37

Promo: 268,50

Newtron® LED handpiece



ref. F12609

A32227

Prijs: 901,61

Promo: 378,50

Newtron® Slim B. LED handpiece



• with blue LED ring: ref. F12900 A38204

• with white LED ring: ref. F12905 A39457

Prijs: 875,34 **Promo: 499,00**

hygiene

Hygiene Kit



ref. F00935


A31336

Prijs: 529,14

Promo: 473,99

periodontics

Perio Kit



ref. F00936

A33170

Prijs: 670,97

Promo: 599,99

implantprotect

pure titanium

Implant Protect Kit



ref. F02120

A39303

Prijs: 735,54

Promo: 682,38

endosuccess

retreatment

Endo Success Retreatment Kit





ref. F00737

A29752

Prijs: 702,36

Promo: 632,40

- Excellent view, full LED ring
- Avoid necrosis up to 120ml irrigation with peristaltic pump
- Improve durability of the handpiece, external irrigation
- Maximum hygiene, sterile single use or autoclavable lines
- Work don't play, use of big bottles (1l) possible
- Large choice of clinical application, 38 tips
- Tips II are compatible with IC2, Piezotome Cube, Piezotome2, Piezotome solo



piezotome cube

**Dynamic power adaptation for
safe bone surgeries**


A41280

~~€7.799,04~~

€5.490,00

Ref. F50100

Delivered with a multifunction footswitch,
a Piezotome® CUBE LED handpiece,
1 handpiece holder, an Essential II kit
(6 tips BS1S, BS4, SL1, SL2, SL3 et LC2),
a Piezotome® dynamometric wrench,
1 irrigation pole, 2 autoclavable irrigation
lines, 5 clips and 15 perforators.



Intraoral surgery ultrasonic control console.
Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

implant center 2

3 modes:

- Piezotome® (bone surgery).
- Newtron® (perio, endo,...).
- I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.

A34090

~~€14.901,47~~

€9590,00



43% korting

PACKAGE DEAL PD0097

Koop een **Implant Center 2** aan **promoprijs** van **9590,99 euro** (waarde: 14901,47 euro) en koop een **Newtron LED handstuk + Newtron LED snoer** aan **super promoprijs 1125,00 euro** (waarde: 2022,97 euro).

Ultrasonic control console for conventional dentistry and intraoral dentistry and motor control console for dental implantology
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France



Bone Surgery™ II Kit* A34094

Ref. F87509 - Tips BS1S II, BS2R II, BS2R II, BS4 II, BS5 II, BS6 II.

~~€1.432,35~~ **€1.111,99**



Sinus Lift II Kit* A34092

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

~~€1.077,92~~ **€1.009,99**



Intralift™ II Kit* A34089

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

~~€1.367,71~~ **€1.111,99**



Extraction™ II Kit* A34088

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

~~€1.087,60~~ **€1.009,99**

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

Accessoires

Promotion valable only for ImplantCenter2

Newtron® LED handpiece



Ref. F12609

For Piezotome® 2 and ImplantCenter™ 2

Class IIa medical devices (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France

Newtron® LED cord



Ref. F12830

For Piezotome® 2 and ImplantCenter™ 2

S32227

~~€2.022,97~~

€1125,00

Universal irrigation lines Sterile x10



A32384

~~€184,07~~

€125,79

Per 4: € 115,79

Ref. F59905 - Single-use with integrated cassette and flow control system.

Class Is medical device (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France

PIEZOTOME® LED handpiece + cord



A46435

~~€1.635,16~~

€1.324,99

Ref. F12815

Class IIa medical devices (GMED) - CE 0459

Manufacturer: SATELEC® a company of ACTEON® group - France



Universal II Kit*

A46420

Ref. F87539 - Tips BS6, LC1, LC2, CE1, CE3 et PZ1.

~~€1.393,47~~ **€1.111,99**



Essential II Kit*

A41909

Ref. F87528 - Tips BS1S II, BS4 II, SL1 II, SL2 II, SL3 II, LC2 II.

~~€1.432,35~~ **€1.111,99**



Crest Splitting II Kit*

A39232

Ref. F87567 - Tips CS1 II, CS2 II, CS3 II, CS4 II, CS5 II, CS6 II.

~~€950,34~~ **€938,70**



Crown Extension II Kit*

A43497

Ref. F87554 - Tips BS6 II, CE1 II, CE2 II, CE3 II.

~~€737,70~~ **€630,99**

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

miniled ortho 2

32% korting

A34280

Prijs: 1486,25

Promo: 999,00

- **Customizable signal**

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

- **Power**

3000mW/cm²
(Ø5,5mm light guide).

- **Long lasting battery**

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

- **Ortho friendly**

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

- **Efficiency & safety**

Cure any type of material used orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.

Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from
few brackets to a whole quadrant



Simply Powerful

Outperforms in
power, efficiency and speed

32% korting

A34219

Prijs: 1435,51

Promo: 975,00



miniled supercharged

• Visibility

Luminous screen with customizable modes and times, count curing time and load indicator.



• Power

2000mW/cm²
(Ø7,5mm light guide).

3000mW/cm²
(Ø5,5mm light guide - optional).

• Optimal result


3 modes to fit all clinical needs.



FAST MODE

The choice of rapidity


Full power during 3, 4, 5 or 10 sec, according to the composite thickness.



PULSE MODE

Limits temperature rise

5 or 10 flashes of 1 sec for a polymerization in security.



RAMPING MODE

Avoids composites retraction

Progressive power increase during 9 or 15 sec to ensure perfect sealing of restoration.

miniled standard

28% korting

A27384

Prijs: 1117,47

Promo: 799,00

- **Fast**

Only 6 to 12 sec to polymerize
2mm composite.

- **Power**

1250mW/cm²
(Ø7,5mm light guide).

2000mW/cm²
(Ø5,5mm light guide - optional).

- **Safety**


Base-charger includes a radiometer for light power
intensity check-up.

- **Optimal result**

3 modes to fit all clinical needs.

Highly Adaptable


The alliance of power & safety



FAST MODE

Fast and easy


Full power during 10 sec
(bip at 5 sec) on regular
depth.



PULSE MODE

Ally rapidity & safety

10 flashes of 1 sec
(bip after 5 flashes)
with times off to
preserve dental pulp.



RAMPING MODE

No stress cracking

Progressive power
increase during 20 sec
(bip every 5 sec) for
a complete and safe
curing.



servotome

Delivered with a single footswitch,
a box with 10 electrodes,
a conductive bracelet.

Ref. F57210



27% korting

A30972

Prijs: 2076,45

Promo: 1499,00

Riskontrol® contributes to clinical protection





RISKONTROL® is also a **GREEN** solution

- Disposable tip
- Economic: simplify process to save time
- Recyclable plastic: environmentally friendly
- Cost-effective
- Flexible for a better use: respect of patient's anatomy
- Safety lock system: no risk of ejecting the tips




Sterilizable tips are not sufficient anymore

RISKONTROL® can be used for all procedures in any workflow




**Easy
to fit**

1




It's installed.

2



It's fitted.

3



It's locked.

Tot 60% korting

Package Deal PD0335

Koop **10 x Riskontrol classic** en/of **Riskontrol art**
naar keuze aan **promoprijs**
en koop een **Riskontrol adapter** aan **1 euro** naar keuze
(waarde: 115,90 euro).

Discover our RISKONTROL® range

Adapt your tips to the taste
of your patients



art
by **riskontrol**

V38000
Prijs: 62,76
Promo: 27,90



Ref. 201 700 - Mint



Ref. 201 701 - Mandarine



Ref. 201 702 - Liquorice



Ref. 201 703 - Anise



Ref. 201 704 - Black currant

Non-sterile disposable tips (x 250)

Add color

to your patients' care experience



V10811
Prijs: 57,78
Promo: 27,90

classic
by **riskontrol**



Ref. 201 664 - White



Ref. 201 665 - Blue



Ref. 201 666 - Yellow



Ref. 201 667 - Green



Ref. 201 668 - Pink

Non-sterile disposable tips (x 250)

Protect entirely your air/water syringe
thanks to the plastic sheath



total protec
by **riskontrol**

V33150
Prijs: 33,59
Promo: 23,50



Ref. 201 750



Ref. 201 751

Non-sterile disposable tips with plastic sheath (x 100)

Choose a totally sterile
workflow



A33761
Prijs: 38,89
Promo: 23,45

sterile
by **riskontrol**



Ref. 201 770

Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile by RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: PRODUITS DENTAIRES PIERRE ROLLAND - ACTEON® Group - FRANCE.



Mouth mirrors

PURE REFLECT®

- Very bright image
- Resistant

For all purposes

RH

- A greyish reflection



TI

- A blueish reflection

FRONT SURFACE

Cross-sectional view





Image reflected directly on the surface of the glass, distortion-free.

NO DOUBLE IMAGE


BACK SURFACE

Cross-sectional view of an aluminium mirror



Reflecting surface under the glass. **DOUBLE IMAGE**

PURE REFLECT®



N°	Ref.	N.
4	MP3622PH	12 A42193
5	MP3624PH	12 A42194

Prijs: 56,53

Promo: 39,99



A TRIPLE ACTION



1 ASTRINGENT EFFECT

helps to completely dry the sulcus with no blood or fluids.



OPEN THE SULCUS

2 effectively with no pain or damage to tissues thanks to the viscosity carefully calculated with our clinical expertise.



3 PRESERVE EPITHELIAL ATTACHMENT*

for a reduced risk of gingival recession and bone resorption.

*Canula must not be inserted into the sulcus

**expasyl
exact**



A44322

Prijs: 144,33

Promo: 100,99

Paste for temporary sulcular opening
20 capsules (0,3g) Expasyl™ Exact - Ref. 261010



Paste for temporary sulcular opening
20 capsules (1g) Expasyl™ - Ref. 261030

expasyl

A10829

Prijs: 248,87

Promo: 109,00



Straight - Ref. 261005 A26105
Bended - Ref. 261015 A33371

Prijs: 111,16

Promo: 68,90

x-mind dc



Bottom
mounted

The X-MIND® generators can accomodate any operatory configuration.

3 arms lengths are available : 41 cm (16.2"), 82.5 cm (32.5") and 110 cm (43.5") and can be:

- Top wall mounted
- Bottom wall mounted

Dental unit, ceiling and mobile version are also available.

With its adaptable wall plate, it is easy to replace an older generator with the X-MIND® whitout having to do additional drilling on the wall.



A37323-11B

~~€4.821,80~~

€3.199

Available with **Remote Wireless Control**



No more wiring issues

Say goodbye to the complicated and expensive wiring work traditionnally needed to connect the remote button in your practice.



Cost effective

Save on installation expenses with the easy-to-integrate wireless solution.



Enhanced convenience

Enjoy the freedom and flexibility of a wireless connection, making your practice more efficient and user-friendly.

A50693

Wireless Button Kit X-MIND DC

Ref. 5.15.00135

All X-Mind DC devices delivered
from April 2021 onwards

Prijs: 299,00

Promo: 249,00



Wall version	Top mounting	Top mounting square cone	Bottom mounting	Bottom mounting square cone
0,40 m extensor	Ref. 5.92.23137	Ref. 5.92.23167	Ref. 5.92.23130	Ref. 5.92.23160
0,80 m extensor	Ref. 5.92.23138	Ref. 5.92.23168	Ref. 5.92.23131	Ref. 5.92.23161
1,10 m extensor	Ref. 5.92.23139	Ref. 5.92.23169	Ref. 5.92.23132	Ref. 5.92.23162

All references for the X-MIND® DC

Choose from a variety of cones to better fit to your clinical needs



Long cone (31 cm - 12'')
Recommended for the parallel technique



Short cone (20 cm - 8'')
Recommended for bisector technique



Rectangular Long cone (31 cm - 12'')
Reduces the X-rays applied to the patient's skin by 50%.

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy

x-mind unity

A39620

Prijs: 5882,40

Promo: 3290,00

44% korting

A sharp and contrasted image

- A contrasted and defined image
- Optimized protection for the patient
- The generator and sensor association

Fixing from above, arm 0.80 m, short cone.
Other configurations: see public price 2024.

Intraoral X-ray

- Very high quality
- Precise and flexible positioning
- Fast selection from traditional to digital mode
- Minimum exposure time: 0,02s
- Programmable timer
- Top or bottom mounting



Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy



sopix 2



A remarkable image quality, even for endodontics and the more complexes

SOPIX^{®2} provides a high resolution image by exploiting the latest technologies available on the market.



SOPIX^{®2} USB Sensor size 1
high definition with a sensor
holder, 10 sensor sheaths
Ref. S_802_0008

~~€7.058,90~~
A42119

**Prijs op
aanvraag**

Also available in size 2
Ref. S_802_0009

~~€8.345,08~~
A48185

**Prijs op
aanvraag**

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON[®] Group- France

Free yourself from constraints
with the embedded system and its
detachable sensor part.



A49822
SOPIX^{®2} Plug-in kit S1
Ref. SK60-1000

~~€6.655,03~~ **€5.999**

A49823
SOPIX^{®2} Plug-in kit S2
Ref. SK60-1001

~~€9.194,22~~ **€7.499**

Size 2
31x42mm

Size 1
25x39mm



sopix² plug-in

The latest technology with no
compromise on image quality with
more safety, freedom and comfort.

Resolution: >18lp/mm

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON[®] Group- France

Meet the X-MIND® optima 3D Family.

Prijs: 71603,20

44% korting

Promo

40049,00

* Option:

The additional cost for the 17x12 version is €3,000

X-MIND® optima 2D/3D wall 100-240V

2D/3D PAN device 5x5 to 12x10 wall installation

A51759 Ref. W1203005



A51760

Prijs: 90168,53

Promo

49499,00

45% korting

* Option:

The additional cost for the 17x12 version is €3,000

X-MIND® optima 2D/3D ceph 100-240V

2D/3D PAN CEPH device 5x5 to 12x10 floor installation

Ref. W1203010

* Option:

The additional cost for the 17x12 version is €3,000

New

X-MIND® optima 2D/3D floor 100-240V

2D/3D PAN device 5x5 to 12x10 floor installation

Ref. W1203005.FL A51761

Prijs: 73194,27

43% korting


Promo

41399,00



AI

Order code	Option
W1000026	RealGUIDE® (lifetime license)
W1203107*	17x12 FOV *Available for FLEX version only
W0900153	DICOM option
W0900046	AIS standalone
W1200129	WeDiagnostix 1 year subscription



ACTEON CLUB
Welcome leaflet
+ Webkey



Workstation and monitor

x-mind optima3D



AIS



RealGUIDE®


Exceptional clarity.

At the heart of the X-MIND® optima 3D is **state-of-the-art detector technology**, which guarantees exceptional image clarity and outstanding anatomical coverage. **Its comprehensive range of Fields of View (FOV)**, from 5x5 to 17x12* cm (Ø x H), makes it suitable for a **wide range of dental procedures**, from general dentistry to orthodontics and maxillofacial surgery, while increasing tolerance for patient positioning and reducing the risk of marginal artifacts.



EXTRAORDINARY CLARITY

Get the most out of quality with the brand-new **ClearIQ** filter



OUTSTANDING VERSATILITY

From 5x5 to 17x12*, a tailored solution for every clinical need



TOP-LEVEL PERFORMANCE

Enhanced image depth with the new kV output range 60-90




MAXIMUM PATIENT COMFORT

Quick and easy acquisition workflow




AUTOMATIC POSITIONERS RECOGNITION

Maximum accuracy



ZERO FOOTPRINT


Maximum accessibility



* Option - The additional cost for the 17x12 version is €3,000

5x5 cm FOV

- Improves the **detection of root canals**, particularly in case of abnormal root canal shapes



8.5x5 cm FOV

- Useful for **implant surgery** and in **general dentistry** as a low dose exam as limited to a single jaw



8.5x9 cm FOV

- Useful for **complex implant surgery** as well as left or right TMJ diagnosis



12x10 cm FOV


- Suitable for most **oral surgery cases** including multiple implant surgery



17x12 cm FOV


New

- Represents the perfect solution for **complex orthognathic cases** **zygomatic implant planning** and **orthodontic surgery**



Meet the 2D Prime Family

X-MIND® prime is available in three configurations: 2D wall mounted, 2D Ceph and 2D floor to suit all installation and clinical demands.



Promo: **13.999,00**

X-MIND® prime WALL 110-240V

Order code: W1200001 A48189

The ACTEON®s smallest extraoral 2D imaging solution ever for the maximum **accessibility in any context.***



with one sensor: Promo: **24.999,00**

with two sensors: Promo: **29.999,00**

X-MIND® prime CEPH 110-240V

Order code: W1200009 A48191

Order code: W1200009 - ds A51763

Cephalometric add-on for a full 2D system still encompassing the same philosophy of prime, **ideal for pediatric applications.**



Promo: **14.999,00**


X-MIND® prime FLOOR 100-240V

Order code: W1200001.FL A51762

*Belgian market only (W1200001/W120001.PB)

Rely on a comprehensive set of 2D exams


X-MIND® prime provides a comprehensive range of 2D exams designed to meet the diverse needs of general dentistry, including dental panoramic imaging, temporomandibular joint assessments and maxillary sinus evaluations.




Dental panoramic



Maxillary sinuses



Temporomandibular joint



Child panoramic



Half-panoramic
Left & right



Bitewing
Left & right



Frontal

PSPiX²

is AUTONOMOUS

The fully automated process free yourself from constraints

Just drop the plate inside the slit, the PSPiX² will automatically detect the size, scan the plate, optimize the image, and send it instantly to your imaging software.

In the meantime, it will automatically erase and eject the plate for an immediate reuse or switch to standby mode.

Scan
simply &
quickly

Fully automated process

42% korting

A51747

—€8.525,20—

€4890,00

Ref. S_702_0013



1

Auto-access:

The door only opens when the imaging plate is inserted.

2

Auto-detect:

PSPiX² automatically detects the size of the inserted imaging plate.

3

Auto-scan:

PSPiX² scans and optimizes the image before displaying it into the right workstation.

4

Auto-eject:

The imaging plate is automatically erased and ejected for an immediate reuse.

5

Eco-mode:

PSPiX² automatically switches to economic mode when not in use for an extended period of time.

Protective bag for ERLM*

For **PSPIX®** and **PSPIX®2**

(previous generation)

V36130

~~€85,82~~

€65,99



	PSPIX® (previous generation)	
	Box of 250	
Size 0		700338L
Size 1		700339L
Size 2		700340L
Size 3		700341L

ERLM plates*

For **PSPIX®2**

V39947

~~€131,16~~

€96,99



	PSPIX®2	
	Box of 2	
Size 0 (22 x 35mm)		990215
Size 1 (24 x 40mm)		990216
Size 2 (31 x 41mm)		990217
Size 3 (27 x 54mm)		990218

Protections for ERLM*

V36129

Protective cover for

PSPIX® and **PSPIX®2**

(previous generation)

~~€48,06~~

€26,99



	Compatible with	
	PSPIX®2	PSPIX®
PSPIX®	✓	✓
SOREDEX® Digora	✓	✓
INSTRUMENTARIUM® Express	✓	✓
DÜRR DENTAL® Vistascan®	✓	✓
KAVO®	✓	✓

	PSPIX® (previous generation)	
	Box of 250	
Size 0		700342L
Size 1		700343L
Size 2		700344L
Size 3		700345L

Protective bag for ERLM*

V39564

For **PSPIX®2**

~~€100,86~~

€70,99



	PSPIX®2	
	Box of 300	
Size 0		700546
Size 1		700547
Size 2		700548
Size 3		700549

- Protection against scratches, light and cross contamination
- Ergonomic and easy to use
- Economic solution

* Luminescent Radio Memory Display for digital phosphor plate systems

PSPIX® imaging plate - Class IIa Medical devices (GMED) - CE0459

Manufacturer: Sopro ACTEON® Group - France

Bag & cover: Class I medical devices - CE


Manufacturer: Sopro ACTEON® Group - France



The thin and flexible imaging plates, available in 4 sizes, will allow you to easily scan all morphologies...



... and all clinical cases



Recupel-bijdrage: 0,10. Prijzen zijn in euro exclusief btw en per stuk of verkoopseenheid. Druk- en zelffouten voorbehouden. Foto's zijn niet bindend. Bij bestellingen boven de 150 euro exclusief btw worden geen administratie - en verzendkosten aangerekend. Toestellen boven de 5000 euro (excl. btw) worden gratis geleverd en geplaatst. Plaatsing omvat steeds de levering, installatie (op de benodigde aansluitingen en leidingen, te voorzien door de klant zelf volgens onze instructies) de oplevering en een korte opleiding. Promoties geldig van **01/10/2026** tot **31/12/2026** en zolang de voorraad strekt. Acties niet cumuleerbaar met andere promoties.

Denta nv, Heiveldekens 2, 2550 Kontich, www.denta.be