

Laboratory

Spin Stand for 6 DP Series Micropipettes

PRECISION, COMFORT & RELIABILITY

 DragLab

The **DragLab Spin Stand** is designed for safe, organized, and convenient storage of up to **six DP Series micropipettes**. Its **360° rotating head** provides easy access from every working angle, while the stable base supports safe single-handed handling. Compatible with both **single-channel and multi-channel DP Series micropipettes**, its lightweight and compact design helps keep laboratory workspaces organized and efficient.

Precision Meets Organization

The **DragLab Spin Stand for 6 DP Series Micropipettes** provides a practical solution for safe, organized, and convenient micropipette storage. Its **six-pipette capacity** allows users to keep multiple DP Series micropipettes securely accessible, while the **360° rotating head** provides easy access from every working angle. The Spin Stand features a stable base for safe single-handed handling, supports both **single-channel and multi-channel DP Series micropipettes**, and combines a lightweight design with easy cleaning and a compact footprint for efficient laboratory organization.



Application:

- **Micropipette Storage** — Safe and organized storage of up to six DP Series micropipettes.
- **Routine Laboratory Work** — Keeps frequently used micropipettes readily accessible at the workstation.
- **Single-Channel Pipetting** — Provides convenient storage and access for single-channel DP Series micropipettes.
- **Multi-Channel Pipetting** — Supports organized storage of compatible multi-channel DP Series micropipettes.
- **Shared Laboratory Workstations** — Helps maintain an organized and space-efficient pipetting area.
- **Laboratory Bench Organization** — The compact footprint helps save valuable bench space while keeping pipettes accessible.

A vertical bar with the colors of the German flag (black, red, and gold) positioned to the left of the title.

Quality and Expertise



DragLab is certified to ISO 9001:2015, reflecting our unwavering commitment to quality and excellence. We continuously optimize our production processes, ensuring that every step is efficient and effective. Through ongoing development and innovation, we strive to meet and exceed industry standards. Our customer-centric approach prioritizes the needs and satisfaction of our clients, driving us to maintain rigorous quality assurance protocols. This dedication ensures that laboratory equipment manufactured by DragLab consistently meets the highest standards of quality and reliability.

DragLab Laboratory Devices are meticulously designed and manufactured from high-quality materials to meet the most rigorous requirements of our customers. Our commitment to excellence and innovation drives continuous development and improvement, informed by daily usage and ongoing technological advancements. 'Quality and Expertise' is not just about creating impressive product designs; it reflects our company's policy of prioritizing high preparedness and customer-centric service as our fundamental objective. This dedication extends to all DragLab products, which comply with current European standards and bear the CE mark.



The **DragLab Spin Stand** securely stores up to **six DP Series micropipettes**, with a **360° rotating head** for easy access. Its stable, compact design supports both **single- and multi-channel pipettes**, keeping laboratory workspaces organized and efficient.



Features:

6-Pipette Capacity

Stores up to **six DP Series micropipettes** safely and keeps them readily accessible.

360° Rotating Head

Provides easy access to each pipette from **every working angle**.

Single & Multi-Channel Compatibility

Designed to accommodate both **single-channel and multi-channel DP Series micropipettes**.

Stable Base

Provides secure support for convenient **single-handed pipette handling**.

Lightweight Design

Easy to position and move around the laboratory as needed.

Easy Cleaning

Can be cleaned with **water or laboratory disinfectants** for practical routine maintenance.

Compact Footprint

Helps save valuable **laboratory bench space** while keeping pipettes organized.

Specifications:

Product Overview	
Product Type	Micropipette Stand
Capacity	6 micropipettes
Compatibility	Single-channel and multi-channel DP Series micropipettes
Rotation	360° rotating head
Material	ABS (Acrylonitrile Butadiene Styrene)
Dimensions	Height: 323 mm × Base Diameter: 195 mm
Color	White & Blue
Autoclavability	Not autoclavable
Certification	CE
Quality Management	ISO 9001:2015
environmental Management	ISO 14001

Order information:

Description	Model	Article number
Spin Stand for 6- DP Series Micropipettes	DP 160	7200160