



# The universal choice for urinalysis.

Designed to work during manual urine microscopy, the Typenex® Medical Universal Tube | 12 mL and Universal Cap are made to fit your processes.

# Universal Tube | 12 mL



## Features and Benefits

**Universal** - The Universal Tube can be used in manual urine microscopy procedures and may also function as a centrifuge tube.

**Ease of Use** - Fill the tube with the desired amount of urine and perform urine microscopy per your facility's protocol.

**Safe** - The cap securely fits onto the tube and prevents spillage in the lab.

## Instructions

1. Fill the Universal Tube | 12 mL with urine to the desired graduation.
2. Place the cap onto the tube if needed.
3. Perform urinalysis per your facility's protocol.

**NOTE:** The Universal Tube | 12 mL is compatible with the KOVA® Urinalysis Kit.

- KOVA® is a registered trademark of KOVA International.

## Ordering Information

Order Number	Description	Unit/Case	Color
US0101QD	Universal Tube   12 mL	500	Clear
US0201QD	Universal Cap	500	Clear

For more information, visit us online at [www.typenexmedical.com](http://www.typenexmedical.com) or call us at 866.897.3639.