


Enhance control

Typenex® Surgical offers a variety of nitinol and stainless steel wires designed to gain access into the upper urinary tract and facilitate the positioning of devices during endourological procedures.

Wire Guides




RocaWire Nitinol

A 100% Nitinol core wire guide, the RocaWire Nitinol has a hydrophilic coating and is radiopaque.



RocaWire PTFE

Stainless steel wire guides are available with a straight or angled tip and coated with PTFE to facilitate smooth placement and exchange.



RioTracer®

The RioTracer® wire guide is available with a straight or angled tip, a stiff nitinol core, a hydrophilic floppy tip, and a barium-filled radiopaque tip to enhance fluoroscopic visualization. The nitinol core provides resistance to kinking for enhanced control. The RioTracer is the optimal choice when navigating tortuous anatomy.

Ordering Information

Product Code	Description
ROWG3500ST	RocaWire Steel Guide Wire with PTFE, 0.035 in x 150 cm, 3 cm straight flexible tip
ROWG3502ST	RocaWire Steel Guide Wire with PTFE, 0.035 in x 150 cm, removable stylet tip
ROWG3510ST	RocaWire Steel Guide Wire with PTFE, 0.035 in x 150 cm, 3 cm flexible "J" tip
ROWP3200ST	RocaWire Hydrophilic Nitinol Guide Wire, 0.032 in x 150 cm, standard, 3 cm straight flexible tip
ROWP3500ST	RocaWire Hydrophilic Nitinol Guide Wire, 0.035 in x 150 cm, standard, 3 cm straight flexible tip
ROWP3502ST	RocaWire Hydrophilic Nitinol Guide Wire, 0.035 in x 150 cm, stiff, 3 cm straight flexible tip
ROWP3503ST	RocaWire Hydrophilic Nitinol Guide Wire, 0.035 in x 180 cm, standard, 3 cm straight flexible tip
ROWP3510ST	RocaWire Hydrophilic Nitinol Guide Wire, 0.035 in x 150 cm, standard, 3 cm flexible "J" tip
ROWP3802ST	RocaWire Hydrophilic Nitinol Guide Wire, 0.038 in x 150 cm, stiff, 3 cm straight flexible tip
ROWR3500ST	RioTracer® Hydrophilic Dual Flexible, 0.035 in x 150 cm, 3 cm tip
ROWR3501ST	RioTracer® Hydrophilic Angled/Straight Dual Flexible, 0.035 in x 150 cm, 3 cm tip
ROWR3502ST	RioTracer® Hydrophilic and Floppy Tip, 0.035 in x 150 cm, 3 cm tip

For more information, visit www.typenexsurgical.com.
Products distributed by Typenex Surgical. Manufactured by Promepia.