

Barcode Plus™ Slim | Handwritten

For hospitals and patient care providers concerned with patient safety during blood transfusions, the Barcode Plus™ Slim Blood Band delivers a unique identifier to help ensure the right blood goes to the right patient.

Easy to Use

The large font size of the Blood Bank ID (BBID) makes it quick and easy to check, and the barcodes allow more efficiency in capturing the BBID in your Laboratory Information System (LIS).

Durable

A shield protects the carbon copy from deteriorating when in contact with water and/or fluid during normal wear for up to 5 days.

R3™ Compatible

Easily reattach, relocate or resize Barcode Plus Slim with the R3™ Reattachment System.

Order Number	Color	Band/Case	Durability	Barcode Symbology
4R5378	Red	150/case	Up to 5 days	GS1-128
4R5379	Green	150/case	Up to 5 days	GS1-128
4R5380	Yellow	150/case	Up to 5 days	GS1-128
4R5381	Blue	150/case	Up to 5 days	GS1-128



For more information, visit us online at www.typenex.com or call us at 866.897.3639.



303 E. Wacker Drive, Suite 1030
Chicago, IL 60601



866.897.3639
866. 897.3630 (fax)



sales@typenex.com



www.typenex.com