



The universal choice for urinalysis.

Designed to work during manual urine microscopy, the Typenex® Medical Universal Tube | 12 mL and Universal Cap are made to fit your processes.

Universal Tube | 12 mL



Features and Benefits

Universal - The Universal Tube can be used in manual urine microscopy procedures and may also function as a centrifuge tube.

Ease of Use - Fill the tube with the desired amount of urine and perform urine microscopy per your facility's protocol.

Safe - The cap securely fits onto the tube and prevents spillage in the lab.

Instructions

1. Fill the Universal Tube | 12 mL with urine to the desired graduation.
2. Place the cap onto the tube if needed.
3. Perform urinalysis per your facility's protocol.

NOTE: The Universal Tube | 12 mL is compatible with the KOVA® Urinalysis Kit.

- KOVA® is a registered trademark of KOVA International.

Ordering Information

Order Number	Description	Unit/Case	Color
US0101QD	Universal Tube 12 mL	500	Clear
US0201QD	Universal Cap	500	Clear

For more information, visit us online at www.typenexmedical.com or call us at 866.897.3639.



303 E. Wacker Drive, Suite 1030
Chicago, IL 60601



866.897.3639
866.897.3630 (fax)



sales@typenex.com



www.typenexmedical.com