

# VIEWSAFE™ MICROSCOPE SLIDES

---



TYPENEX   
MEDICAL

# INTRODUCTION


---

Typenex® Medical is proud to offer a full portfolio of exceptional microscope slides — all pre-cleaned, ready to use, and produced from high-quality glass. They may be used for a wide range of applications and can accommodate various marking methods such as conventional labeling, writing, or printing. With multiple colors, styles, and features available, our microscope slides are thoughtfully tailored to suit your laboratory needs.

## OVERVIEW


---

Our microscope slides are offered in various corner/edge types:




### TRADITIONAL STYLE

- 90° ground edge
- 90° corners



### SAFETY CORNERS

- 90° ground edge
- Slightly clipped corners



### CLIPPED CORNERS

- 45° beveled edge
- No sharp corners

## VIEWSAFE™ PREMIUM MICROSCOPE SLIDES

ViewSafe Premium Microscope Slides are made from chemically pure white sheets of glass for superior transparency. The slides are designed to help maintain specimen integrity and accurate results.


- Dimensions: 25 x 75 mm
- Thickness: 1.05 mm
- Packaging: 72/box, 20 boxes/case (10 gross)

## VIEWSAFE™ ECONOMY MICROSCOPE SLIDES

ViewSafe Economy Microscope Slides are produced from green float glass sheets for remarkable quality. This product line offers considerable cost savings and can be used for general laboratory applications.

- Dimensions: 25 x 75 mm
- Thickness: 1.10 mm
- Packaging: 72/box, 20 boxes/case (10 gross)


# VIEWSAFE™ PREMIUM MICROSCOPE SLIDES



## PREMIUM PLAIN

Premium Plain Microscope Slides provide a clear, classic surface produced from the finest quality of white glass. They can be used for all routine laboratory applications.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
PS0101	Premium Microscope Slides	Clear	Ground	90°	1440
PS0103	Premium Microscope Slides	Clear	Beveled	Clipped	1440




## PREMIUM FROSTED

Premium Frosted Microscope Slides feature a frosted coating that is resistant to stains and chemicals used in laboratories. This finish can be applied to one or both sides of one end.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
PS0201	Premium Microscope Slides	Frosted	Ground	90°	1440
PS0203	Premium Microscope Slides	Frosted	Beveled	Clipped	1440
PS0301	Premium Microscope Slides	Frosted	Ground	90°	1440
PS0303	Premium Microscope Slides	Frosted, Both Sides	Beveled	Clipped	1440


# VIEWSAFE™ PREMIUM MICROSCOPE SLIDES



## PREMIUM COLOR-CODED

Premium Color-Coded Microscope Slides are colored on one end to help ensure specimen identification and traceability. The colored ends are resistant to stains and chemicals commonly used in laboratories.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
PS0201A	Premium Microscope Slides	Aqua Frosted ●	Ground	90°	1440
PS0201B	Premium Microscope Slides	Blue Frosted ●	Ground	90°	1440
PS0201D	Premium Microscope Slides	Gold Frosted ●	Ground	90°	1440
PS0201G	Premium Microscope Slides	Green Frosted ●	Ground	90°	1440
PS0201L	Premium Microscope Slides	Lilac Frosted ●	Ground	90°	1440
PS0201N	Premium Microscope Slides	Orange Frosted ●	Ground	90°	1440
PS0201P	Premium Microscope Slides	Pink Frosted ●	Ground	90°	1440
PS0201T	Premium Microscope Slides	Tan Frosted ●	Ground	90°	1440
PS0201W	Premium Microscope Slides	White Frosted ○	Ground	90°	1440
PS0201Y	Premium Microscope Slides	Yellow Frosted ●	Ground	90°	1440
PS0203B	Premium Microscope Slides	Blue Frosted ●	Beveled	Clipped	1440
PS0203G	Premium Microscope Slides	Green Frosted ●	Beveled	Clipped	1440
PS0203N	Premium Microscope Slides	Orange Frosted ●	Beveled	Clipped	1440
PS0203P	Premium Microscope Slides	Pink Frosted ●	Beveled	Clipped	1440
PS0203W	Premium Microscope Slides	White Frosted ○	Beveled	Clipped	1440
PS0203Y	Premium Microscope Slides	Yellow Frosted ●	Beveled	Clipped	1440




## PREMIUM CHARGED

Premium Charged Microscope Slides feature a hydrophobic surface for optimal cell and tissue adhesion. The adhesive technology and color-coded ends make them ideal for use in histology, cytology, and pathology departments.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
PS0401A	Charged Microscope Slides	Aqua Frosted ●	Ground	90°	1440
PS0401B	Charged Microscope Slides	Blue Frosted ●	Ground	90°	1440
PS0401D	Charged Microscope Slides	Gold Frosted ●	Ground	90°	1440
PS0401G	Charged Microscope Slides	Green Frosted ●	Ground	90°	1440
PS0401L	Charged Microscope Slides	Lilac Frosted ●	Ground	90°	1440
PS0401N	Charged Microscope Slides	Orange Frosted ●	Ground	90°	1440
PS0401P	Charged Microscope Slides	Pink Frosted ●	Ground	90°	1440
PS0401T	Charged Microscope Slides	Tan Frosted ●	Ground	90°	1440
PS0401W	Charged Microscope Slides	White Frosted ○	Ground	90°	1440
PS0401Y	Charged Microscope Slides	Yellow Frosted ●	Ground	90°	1440
PS0402B	Charged Microscope Slides	Blue Frosted ●	Beveled	Clipped	1440
PS0402G	Charged Microscope Slides	Green Frosted ●	Beveled	Clipped	1440
PS0402N	Charged Microscope Slides	Orange Frosted ●	Beveled	Clipped	1440
PS0402P	Charged Microscope Slides	Pink Frosted ●	Beveled	Clipped	1440
PS0402W	Charged Microscope Slides	White Frosted ○	Beveled	Clipped	1440
PS0402Y	Charged Microscope Slides	Yellow Frosted ●	Beveled	Clipped	1440


# VIEWSAFE™ ECONOMY MICROSCOPE SLIDES



## ECONOMY PLAIN

Economy Plain Microscope Slides are ideal for use in all routine laboratory applications. These slides provide exceptional quality while offering considerable cost savings.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
ES0101	Economy Microscope Slides	Clear	Ground	90°	1440
ES0102	Economy Microscope Slides	Clear	Ground	Safety	1440
ES0103	Economy Microscope Slides	Clear	Beveled	Clipped	1440



## ECONOMY FROSTED






















Economy Frosted Microscope Slides feature a frosted coating that is resistant to stains and chemicals used in laboratories. This finish can be applied to one or both sides of one end.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
ES0201	Economy Microscope Slides	Frosted	Ground	90°	1440
ES0202	Economy Microscope Slides	Frosted	Ground	Safety	1440
ES0203	Economy Microscope Slides	Frosted	Beveled	Clipped	1440
ES0301	Economy Microscope Slides	Frosted, Both Sides	Ground	90°	1440
ES0302	Economy Microscope Slides	Frosted, Both Sides	Ground	Safety	1440
ES0303	Economy Microscope Slides	Frosted, Both Sides	Beveled	Clipped	1440



## ECONOMY COLOR-CODED

Economy Color-Coded Microscope Slides are colored on one end to help ensure specimen identification and traceability. They are available in a number of bright hues to accommodate various procedures.

PRODUCT CODE	DESCRIPTION	COLOR	EDGE	CORNER	UNITS/CASE
ES0201B	Economy Microscope Slides	Blue Frosted 	Ground	90°	1440
ES0201G	Economy Microscope Slides	Green Frosted 	Ground	90°	1440
ES0201N	Economy Microscope Slides	Orange Frosted 	Ground	90°	1440
ES0201P	Economy Microscope Slides	Pink Frosted 	Ground	90°	1440
ES0201T	Economy Microscope Slides	Tan Frosted 	Ground	90°	1440
ES0201W	Economy Microscope Slides	White Frosted 	Ground	90°	1440
ES0201Y	Economy Microscope Slides	Yellow Frosted 	Ground	90°	1440
ES0202A	Economy Microscope Slides	Aqua Frosted 	Ground	Safety	1440
ES0202B	Economy Microscope Slides	Blue Frosted 	Ground	Safety	1440
ES0202D	Economy Microscope Slides	Gold Frosted 	Ground	Safety	1440
ES0202G	Economy Microscope Slides	Green Frosted 	Ground	Safety	1440
ES0202L	Economy Microscope Slides	Lilac Frosted 	Ground	Safety	1440
ES0202N	Economy Microscope Slides	Orange Frosted 	Ground	Safety	1440
ES0202P	Economy Microscope Slides	Pink Frosted 	Ground	Safety	1440
ES0202T	Economy Microscope Slides	Tan Frosted 	Ground	Safety	1440
ES0202W	Economy Microscope Slides	White Frosted 	Ground	Safety	1440
ES0202Y	Economy Microscope Slides	Yellow Frosted 	Ground	Safety	1440
ES0203B	Economy Microscope Slides	Blue Frosted 	Beveled	Clipped	1440
ES0203G	Economy Microscope Slides	Green Frosted 	Beveled	Clipped	1440
ES0203N	Economy Microscope Slides	Orange Frosted 	Beveled	Clipped	1440
ES0203P	Economy Microscope Slides	Pink Frosted 	Beveled	Clipped	1440



# MICROSCOPE SLIDE COVER GLASS



## VIEWSAFE™ COVER GLASS

Our quality microscope slide cover glass is designed with exceptional optical properties, stable molecular structure, and uniform flatness. It is available in three different thickness formats for added convenience.

- No. 1: 0.13 – 0.16 mm
- No. 1.5: 0.16 – 0.19 mm
- No. 2: 0.19 – 0.22 mm
- Packaging: 1 oz/pack, 10 packs/box (10 oz)

PRODUCT CODE			SIZE	APPROXIMATE UNITS/PACK		
NO. 1	NO. 1.5	NO. 2		NO. 1	NO. 1.5	NO. 2
MC0101	—	MC0102	18 x 18 mm	248	—	168
MC0103	—	MC0104	20 x 20 mm	192	—	154
MC0105	MC0106	MC0107	22 x 22 mm	168	148	136
MC0110	—	MC0111	22 x 32 mm	114	—	84
MC0114	—	MC0115	22 x 40 mm	86	—	75
MC0118	—	MC0119	22 x 50 mm	70	—	63
MC0123	—	MC0124	22 x 60 mm	60	—	50
MC0108	—	MC0109	24 x 24 mm	141	—	102
MC0112	—	—	24 x 32 mm	104	—	—
MC0116	—	MC0117	24 x 40 mm	82	—	71
MC0120	MC0121	MC0122	24 x 50 mm	65	56	56
MC0125	—	MC0126	24 x 60 mm	55	—	46

Visit [www.typhenex.com](http://www.typhenex.com) or call 866-897-3639 to learn more.