




# Quality you can count on.

Typenex® Medical offers a quality selection of microscope slides designed to support routine applications in liquid-based cytology, such as TCT and LCT slide preparation. These microscope slides feature hydrophilic surfaces, ground edges, 45° clipped corners, and various configurations to suit your specific liquid cytology needs.



**ViewSafe™** | Liquid-Based Cytology Microscope Slides

# ViewSafe™ | Liquid-Based Cytology Microscope Slides



## Features and Benefits

**Various Configurations** — Printed brackets, circles, coated, and non-coated styles available to support a wide range of procedures.

**Marking Areas** — White frosted coating applied to one end.

**Standard Dimensions** — 25 x 75 mm, 1–1.2 mm thickness.

## Ordering Information

Part Number	Description	Circles	Edges	Corners	Charged	Units/Case
LC0101	Liquid-Based Cytology Microscope Slides	1	Ground	Clipped	No	1440
LC0102	Liquid-Based Cytology Microscope Slides	1	Ground	Clipped	Yes	1440
LC0201	Liquid-Based Cytology Microscope Slides	2, Stacked	Ground	Clipped	No	1440
LC0202	Liquid-Based Cytology Microscope Slides	2, Stacked	Ground	Clipped	Yes	1440
LC0203	Liquid-Based Cytology Microscope Slides	2, Offset	Ground	Clipped	No	1440
LC0204	Liquid-Based Cytology Microscope Slides	2, Offset	Ground	Clipped	Yes	1440
LC0301	Liquid-Based Cytology Microscope Slides	2 Brackets	Ground	Clipped	Yes	1440

Microscope slide mailers and storage sold separately.

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.