

— TAILORED —
CLOSED WOUND DRAINAGE
— SOLUTIONS —



TYPENEX 
MEDICAL

INDUSTRY CONTEXT

Typenex® Medical offers an extensive closed wound drainage portfolio of more than 40 products.

A closed wound drainage system can assist in the recovery process. The system includes a long, silicone tube that inserts into the wound site. This drain connects to an external reservoir that collects the drained exudate. The system is often specific to a particular procedure or patient.

PRODUCTS

KITS

- Kits efficiently package drains, evacuators, and an optional trocar.
 - Trocar: A round shaft with a sharp, three-sided point
 - Drain: Comes pre-attached to help the drain insertion process
 - Reservoir: Includes 100 cc bulb evacuator or 400 ml 3-spring reservoir



PRODUCT CODE	SIZE	DRAIN TYPE	DRAIN SHAPE	TROCAR	ADAPTER	BULB	MATERIAL	BARD®1 CROSS REFERENCE	CARDINAL®2 CROSS REFERENCE	MEDLINE®3 CROSS REFERENCE	J&J ETHICON®4 CROSS REFERENCE
WD0401K1	7 mm	Full Perf	Flat	No	No	100 cc	Silicone	0071430	SU130-1360	DYNJWE1360	N/A
WD0401PK1	7 mm	¾ Perf	Flat	No	No	100 cc	Silicone	0071360	SU130-1348	DYNJWE1348	N/A
WD0402K1	10 mm	Full Perf	Flat	No	No	100 cc	Silicone	0071440	SU130-1361	DYNJWE1361	N/A
WD0402PK1	10 mm	¾ Perf	Flat	No	No	100 cc	Silicone	0071370	SU130-1349	DYNJWE1349	N/A
WD0503TK1	15 FR	Full Fluted	Round	Yes	Yes	100 cc	Silicone	0073320	N/A	DYNJWE2264	N/A
WD0504TK1	19 FR	Full Fluted	Round	Yes	Yes	100 cc	Silicone	N/A	N/A	DYNJWE2265	N/A
WD0701TK5	7 FR	Full Perf	Round	Yes	No	400 ml	PVC	43600	SU130-401D	DYNJWE401	N/A
WD0702TK5	10 FR	Full Perf	Round	Yes	No	400 ml	PVC	43610	SU130-402D	DYNJWE402	N/A
WD0703TK5	15 FR	Full Perf	Round	Yes	No	400 ml	PVC	43630	SU130-403D	DYNJWE403	N/A
WD0704TK5	19 FR	Full Perf	Round	Yes	No	400 ml	PVC	43620	SU130-404D	DYNJWE404	N/A

BULBS AND 3-SPRING RESERVOIR

- Includes a clip for convenient garment and linen attachment
- Sterile




PRODUCT CODE	SIZE	MATERIAL	BARD®1 CROSS REFERENCE	CARDINAL®2 CROSS REFERENCE	MEDLINE®3 CROSS REFERENCE	J&J ETHICON®4 CROSS REFERENCE
WD0101	100 cc	Silicone	0070740	SU130-1305	DYNJWE1305	2160
WD0104	400 cc	Silicone	0070400	SU130-1000	DYNJWE1000	2162 2163
WD0204	400 ml	PVC	N/A	SU130-475	DYNJWE475	N/A

PRODUCTS

(CONTINUED)

FLUTED DRAINS (COMPARABLE TO BLAKE® DRAINS⁴)

- Channel design minimizes tissues in-growth, which may help reduce patient discomfort upon removal⁶
- Radiopaque
- Sterile
- ¾ and full-length options
- Channels make up the inserted portion of the drain to assist with drainage and to lower the risk of blockage



PRODUCT CODE	SIZE	DRAIN TYPE	DRAIN SHAPE	TROCAR	ADAPTER	MATERIAL	BARD®1 CROSS REFERENCE	CARDINAL®2 CROSS REFERENCE	MEDLINE®3 CROSS REFERENCE	J&J ETHICON®4 CROSS REFERENCE
WD0601	7 mm	Full Fluted	Flat	No	No	Silicone	072211	JP-2211	DYNJWE2211	2211
WD0601P	7 mm	¾ Fluted	Flat	No	No	Silicone	072210	JP-2210	DYNJWE2210	2210
WD0601T	7 mm	Full Fluted	Flat	Yes	No	Silicone	072212	JP-2212	DYNJWE2212	2212
WD0602	10 mm	Full Fluted	Flat	No	No	Silicone	072214	JP-2214	DYNJWE2214	2214
WD0602P	10 mm	¾ Fluted	Flat	No	No	Silicone	072213	JP-2213	DYNJWE2213	2213
WD0602T	10 mm	Full Fluted	Flat	Yes	No	Silicone	072215	JP-2215	DYNJWE2215	2215
WD0602PT	10 mm	¾ Fluted	Flat	Yes	No	Silicone	072217	JP-2217	N/A	2217
WD0502	10 FR	Full Fluted	Round	No	Yes	Silicone	072186	JP-2186 JP-2226	DYNJWE2186	2226
WD0502T	10 FR	Full Fluted	Round	Yes	Yes	Silicone	072187	JP-2187 JP-2227	DYNJWE2187	2227
WD0503	15 FR	Full Fluted	Round	No	Yes	Silicone	072188	JP-2188 JP-2228	DYNJWE2188	2228
WD0503P	15 FR	¾ Fluted	Round	No	No	Silicone	072228	N/A	N/A	N/A
WD0503T	15 FR	Full Fluted	Round	Yes	Yes	Silicone	072189	JP-2229 JP-2189	DYNJWE2189	2229
WD0503PT	15 FR	¾ Fluted	Round	Yes	No	Silicone	072223 072229	JP-2223	N/A	N/A
WD0503B	15 FR	Full Fluted	Round	Bendable	Yes	Silicone	N/A	JP-2233	N/A	2233
WD0504	19 FR	Full Fluted	Round	No	Yes	Silicone	072190	JP-2190 JP-2230	DYNJWE2190	2230
WD0504P	19 FR	¾ Fluted	Round	No	No	Silicone	072230	N/A	N/A	N/A
WD0504T	19 FR	Full Fluted	Round	Yes	Yes	Silicone	072191 072192	JP-2231 JP-2191	DYNJWE2191	2231
WD0504PT	19 FR	¾ Fluted	Round	Yes	No	Silicone	072225 072231	JP-2225	N/A	N/A
WD0504B	19 FR	Full Fluted	Round	Bendable	Yes	Silicone	N/A	JP-2232	N/A	2232
WD0505	24 FR	Full Fluted	Round	No	Yes	Silicone	072234 072235	JP-2234	DYNJWE2234	2234
WD0506	28 FR	Full Fluted	Round	No	Yes	Silicone	072290	JP-2290	N/A	N/A

PRODUCTS

(CONTINUED)

PERFORATED SILICONE DRAINS (COMPARABLE TO JP® DRAINS⁵)

- Perforations dispersed along the drain for effective drainage
- Radiopaque
- Not made with natural rubber latex
- Sterile
- ¾, mid-, and full-length perforated options



PRODUCT CODE	SIZE	DRAIN TYPE	DRAIN SHAPE	TROCAR	ADAPTER	MATERIAL	BARD® ¹ CROSS REFERENCE	CARDINAL® ² CROSS REFERENCE	MEDLINE® ³ CROSS REFERENCE	J&J ETHICON® ⁴ CROSS REFERENCE
WD0401	7 mm	Full Perf	Flat	No	No	Silicone	0070430	SU130-1310	DYNJWE1310	N/A
WD0401P	7 mm	¾ Perf	Flat	No	No	Silicone	0070360	SU130-1308	DYNJWE1308A	N/A
WD0401T	7 mm	Full Perf	Flat	Yes	No	Silicone	N/A	SU130-1410	DYNJWE1410	N/A
WD0402	10 mm	Full Perf	Flat	No	No	Silicone	0070440	SU130-1311	DYNJWE1311	N/A
WD0402P	10 mm	¾ Perf	Flat	No	No	Silicone	0070370	SU130-1309	DYNJWE1309	N/A
WD0402T	10 mm	Full Perf	Flat	Yes	No	Silicone	N/A	SU130-1411	DYNJWE1411	N/A
WD0301	7 FR	Full Perf	Round	No	Yes	Silicone	0070200	SU130-1320	DYNJWE1320	N/A
WD0302	10 FR	Full Perf	Round	No	Yes	Silicone	0070210	SU130-1321	DYNJWE1321A	N/A
WD0302T	10 FR	Full Perf	Round	Yes	Yes	Silicone	0070310	SU130-0321	DYNJWE0321	N/A
WD0303	15 FR	Full Perf	Round	No	No	Silicone	0070220	SU130-1323	DYNJWE1323A	N/A
WD0303T	15 FR	Full Perf	Round	Yes	No	Silicone	0070320	SU130-0323	DYNJWE0323	N/A
WD0303MT	15 FR	Mid Perf	Round	Yes	No	Silicone	0070720	SU130-1329	N/A	N/A
WD0304	19 FR	Full Perf	Round	No	No	Silicone	0070230	SU130-1325	DYNJWE1325A	N/A
WD0304T	19 FR	Full Perf	Round	Yes	No	Silicone	0070330	SU130-0325	DYNJWE0325	N/A
WD0701T	7 FR	Full Perf	Round	Yes	No	PVC	80700	SU130-0520	N/A	N/A
WD0702T	10 FR	Full Perf	Round	Yes	No	PVC	80710	SU130-0521	DYNJWE0521	N/A
WD0703T	15 FR	Full Perf	Round	Yes	No	PVC	80720	SU130-052	DYNJWE0523	N/A
WD0704T	19 FR	Full Perf	Round	Yes	No	PVC	80730	SU130-0525	DYNJWE0525	N/A

Contact your Typenex Medical sales representative at 866.897.3639 to learn how you can start improving outcomes and reducing costs.

For more information, visit www.typenexmedical.com.

¹ Bard is a registered trademark of C.R. Bard, Inc.

² Cardinal Health is a registered trademark of Cardinal Health Technologies, LLC.

³ Medline is a registered trademark of Medline Industries, Inc.

⁴ Blake and Ethicon are registered trademarks of Johnson + Johnson Corporation.

⁵ Jackson-Pratt is a registered trademark of Cardinal Health.

⁶ Scott-Conner, Carol E.H., editor. Chassin's Operative Strategy in General Surgery: An Expositive Atlas. 4th ed., Springer-Verlag New York, 2014, p. 64.