




HISTOLOGY CASSETTES








Made to match your process.

Typenex® Medical provides a complete portfolio of histology cassettes, base molds, and more to support a wide variety of laboratory applications. Our offerings are cost-effective, made with quality materials, and resistant to the chemical action of histological solvents. Discover our thoughtfully tailored selection of product features, packaging styles, colors, and supplementary accessories to choose the solution that best suits your testing process.



Biopsy Cassettes

- Features a 45° writing area angle
- Suitable for processing small specimens without a biopsy pad
- Maximize fluid exchange with .26 x .26 mm or .9 x .9 mm openings
- Offered in various packaging styles to accommodate most printers

BIOPSY CASSETTES



EC0104

EC0104 | 0.26 x 0.26 mm | 45° | Front Hinge | Lid Attached

Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0104WXB	EC0104WXS	EC0104WXT
Yellow	EC0104YXB	EC0104YXS	EC0104YXT
Green	EC0104GXB	EC0104GXS	EC0104GXT
Pink	EC0104KXB	EC0104KXS	EC0104KXT
Blue	EC0104BXB	EC0104BXS	EC0104BXT
Tan	EC0104TXB	EC0104TXS	EC0104TXT
Purple	EC0104PXB	EC0104PXS	EC0104PXT
Peach	EC0104HXB	EC0104HXS	EC0104HXT
Orange	EC0104NXB	EC0104NXS	EC0104NXT
Aqua	EC0104AXB	EC0104AXS	EC0104AXT
Gray	Available Upon Request	EC0104EXS	Available Upon Request



EC0203

EC0203 | 0.9 x 0.9 mm | 45° | Back Hinge | Lid Attached | One-Piece Set

Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0203WXB	-	EC0203WXT
Yellow	EC0203YXB	-	EC0203YXT
Green	EC0203GXB	-	EC0203GXT
Pink	EC0203KXB	-	EC0203KXT
Blue	EC0203BXB	-	EC0203BXT
Tan	EC0203TXB	-	EC0203TXT
Purple	EC0203PXB	-	EC0203PXT
Peach	EC0203HXB	-	EC0203HXT
Orange	EC0203NXB	-	EC0203NXT
Aqua	EC0203AXB	-	EC0203AXT


BIOPSY CASSETTES



EC0204

EC0204 | 0.9 x 0.9 mm | 45° | Front Hinge | Lid Attached

Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0204WXB	EC0204WXS	EC0204WXT
Yellow	EC0204YXB	EC0204YXS	EC0204YXT
Green	EC0204GXB	EC0204GXS	EC0204GXT
Pink	EC0204KXB	EC0204KXS	EC0204KXT
Blue	EC0204BXB	EC0204BXS	EC0204BXT
Tan	EC0204TXB	EC0204TXS	EC0204TXT
Purple	EC0204PXB	EC0204PXS	EC0204PXT
Peach	EC0204HXB	EC0204HXS	EC0204HXT
Orange	EC0204NXB	EC0204NXS	EC0204NXT
Aqua	EC0204AXB	EC0204AXS	EC0204AXT



EC0205

EC0205 | 0.9 x 0.9 mm | 45° | Back Hinge | Lid Detached

Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0205WXB	EC0205WXS	EC0205WXT
Yellow	EC0205YXB	EC0205YXS	EC0205YXT
Green	EC0205GXB	EC0205GXS	EC0205GXT
Pink	EC0205KXB	EC0205KXS	EC0205KXT
Blue	EC0205BXB	EC0205BXS	EC0205BXT
Tan	EC0205TXB	EC0205TXS	EC0205TXT
Purple	EC0205PXB	EC0205PXS	EC0205PXT
Peach	EC0205HXB	EC0205HXS	EC0205HXT
Orange	EC0205NXB	EC0205NXS	EC0205NXT
Aqua	EC0205AXB	EC0205AXS	EC0205AXT

BIOPSY CASSETTES



EC0212

EC0212 0.9 x 0.9 mm 45° Back Hinge Lid Attached			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0212WXB	-	-
Yellow	EC0212YXB	-	-
Green	EC0212GXB	-	-
Pink	EC0212KXB	-	-
Blue	EC0212BXB	-	-
Tan	EC0212TXB	-	-
Purple	EC0212PXB	-	-
Peach	EC0212HXB	-	-
Orange	EC0212NXB	-	-
Aqua	EC0212AXB	-	-



EC0404


EC0404 0.26 x 0.26 mm 45° Front Hinge Lid Attached Four Compartments			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0404WXB	EC0404WXS	EC0404WXT
Yellow	EC0404YXB	EC0404YXS	EC0404YXT
Green	EC0404GXB	EC0404GXS	EC0404GXT
Pink	EC0404KXB	EC0404KXS	EC0404KXT
Blue	EC0404BXB	EC0404BXS	EC0404BXT
Tan	EC0404TXB	EC0404TXS	EC0404TXT
Purple	EC0404PXB	EC0404PXS	EC0404PXT
Peach	EC0404HXB	EC0404HXS	EC0404HXT
Orange	EC0404NXB	EC0404NXS	EC0404NXT
Aqua	EC0404AXB	EC0404AXS	EC0404AXT



Embedding Cassettes

- Features a 45° writing area angle
- No lid included, see Stainless Steel Lids
- Designed with 2.4 mm round holes
- Offered in various packaging styles to accommodate most printers

EMBEDDING CASSETTES




EC0506AXB

EC0506 2.4 x 2.4 mm 45° Front Hinge No Lid			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0506WXB	EC0506WXS	EC0506WXT
Yellow	EC0506YXB	EC0506YXS	EC0506YXT
Green	EC0506GXB	EC0506GXS	EC0506GXT
Pink	EC0506KXB	EC0506KXS	EC0506KXT
Blue	EC0506BXB	EC0506BXS	EC0506BXT
Tan	EC0506TXB	EC0506TXS	EC0506TXT
Purple	EC0506PXB	EC0506PXS	EC0506PXT
Peach	EC0506HXB	EC0506HXS	EC0506HXT
Orange	EC0506NXB	EC0506NXS	EC0506NXT
Aqua	EC0506AXB	EC0506AXS	EC0506AXT
Gray	Available Upon Request	EC0506EXS	Available Upon Request



EC0814
Stainless Steel Lid on
Embedding Cassette

STAINLESS STEEL LIDS		
Product Code	Description	Units/Case
ECO814	Stainless Steel Lid for Embedding Cassette	10



Slotted Cassettes

- Available in 35° and 45° writing area angles
- Lid locks securely onto cassette to minimize risk of specimen loss
- Designed with 5 x 1 mm rectangular holes for optimal fluid exchange
- Offered in various packaging styles to accommodate most printers

SLOTTED CASSETTES



EC0602

EC0602 5 x 1 mm 35° Back Hinge Lid Attached One-Piece Set			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0602WXB	–	EC0602WXT
Yellow	EC0602YXB	–	EC0602YXT
Green	EC0602GXB	–	EC0602GXT
Pink	EC0602KXB	–	EC0602KXT
Blue	EC0602BXB	–	EC0602BXT
Tan	EC0602TXB	–	EC0602TXT
Purple	EC0602PXB	–	EC0602PXT
Peach	EC0602HXB	–	EC0602HXT
Orange	EC0602NXB	–	EC0602NXT
Aqua	EC0602AXB	–	EC0602AXT



EC0603

EC0603 5 x 1 mm 45° Back Hinge Lid Attached One-Piece Set			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0603WXB	–	EC0603WXT
Yellow	EC0603YXB	–	EC0603YXT
Green	EC0603GXB	–	EC0603GXT
Pink	EC0603KXB	–	EC0603KXT
Blue	EC0603BXB	–	EC0603BXT
Tan	EC0603TXB	–	EC0603TXT
Purple	EC0603PXB	–	EC0603PXT
Peach	EC0603HXB	–	EC0603HXT
Orange	EC0603NXB	–	EC0603NXT
Aqua	EC0603AXB	–	EC0603AXT
Gray	Available Upon Request	–	EC0603EXT

SLOTTED CASSETTES



EC0605

EC0605 5 x 1 mm 45° Back Hinge Lid Detached			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0605WXB	EC0605WXS	EC0605WXT
Yellow	EC0605YXB	EC0605YXS	EC0605YXT
Green	EC0605GXB	EC0605GXS	EC0605GXT
Pink	EC0605KXB	EC0605KXS	EC0605KXT
Blue	EC0605BXB	EC0605BXS	EC0605BXT
Tan	EC0605TXB	EC0605TXS	EC0605TXT
Purple	EC0605PXB	EC0605PXS	EC0605PXT
Peach	EC0605HXB	EC0605HXS	EC0605HXT
Orange	EC0605NXB	EC0605NXS	EC0605NXT
Aqua	EC0605AXB	EC0605AXS	EC0605AXT
Gray	Available Upon Request	EC0605EXS	Available Upon Request



EC0607

EC0607 5 x 1 mm 45° Front Hinge Lid Attached Rounded Tabs			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0607WXB	EC0607WXS	EC0607WXT
Yellow	EC0607YXB	EC0607YXS	EC0607YXT
Green	EC0607GXB	EC0607GXS	EC0607GXT
Pink	EC0607KXB	EC0607KXS	EC0607KXT
Blue	EC0607BXB	EC0607BXS	EC0607BXT
Tan	EC0607TXB	EC0607TXS	EC0607TXT
Purple	EC0607PXB	EC0607PXS	EC0607PXT
Peach	EC0607HXB	EC0607HXS	EC0607HXT
Orange	EC0607NXB	EC0607NXS	EC0607NXT
Aqua	EC0607AXB	EC0607AXS	EC0607AXT
Gray	Available Upon Request	EC0607EXS	Available Upon Request

SLOTTED CASSETTES



EC0612

EC0612 5 x 1 mm 45° Back Hinge Lid Attached			
Color	Bulk (500/Case)	Sleeved (750/Case)	Taped (2000/Case)
White	EC0612WXB	-	-
Yellow	EC0612YXB	-	-
Green	EC0612GXB	-	-
Pink	EC0612KXB	-	-
Blue	EC0612BXB	-	-
Tan	EC0612TXB	-	-
Purple	EC0612PXB	-	-
Peach	EC0612HXB	-	-
Orange	EC0612NXB	-	-
Aqua	EC0612AXB	-	-
Gray	EC0612EXB	-	-



Base Molds

- Designed to support routine paraffin embedding
- Features polished surfaces and well-ground corners
- Accommodates any tissue specimen size

BASE MOLDS



CM0102

STAINLESS STEEL BASE MOLDS





Product Code	Dimensions (mm)	Units/Case
CM0101	7 x 7 x 7	10
CM0102	15 x 15 x 7	10
CM0103	24 x 24 x 7	10
CM0104	30 x 24 x 7	10
CM0105	37 x 24 x 7	10
CM0106	7 x 7 x 9	10
CM0107	15 x 15 x 9	10
CM0108	24 x 24 x 9	10
CM0109	30 x 24 x 9	10
CM0110	37 x 24 x 9	10



CM0212

DISPOSABLE BASE MOLDS

Product Code	Dimensions (mm)	Units/Case
CM0211	7 x 7 x 5	500
CM0212	17 x 17 x 5	500
CM0213	25 x 25 x 5	500
CM0214	25 x 32 x 5	500
CM0215	25 x 40 x 5	500
CM0216	30 x 20 x 12	50



To learn more, visit www.typhenex.com or contact us at 866-897-3639.

 303 East Wacker Drive, Suite 1030
Chicago, Illinois 60601 USA

 866-897-3639
866-897-3630 (fax)

 lab@typhenex.com

 www.typhenex.com