





Configure your cleanup.

Typenex® Medical's Absorbent Floor Mat Rolls equip healthcare professionals with an easy, efficient approach to keep operating room floors clean. With varying sizes and styles, our Absorbent Floor Mat Rolls offer a tailored solution specific to your facility.

Absorbent Floor Mat Rolls



Choose from a variety of options to find the right floor fluid management solution for your facility. The Double-Sided Mat Roll effectively captures fluid during and after surgery with its modular design. The Waterproof Backsheet Mat Roll can be placed in the operating room before surgery, providing a tight grip on a dry surface. Both are available in standard thickness or a lightweight “E Roll” for low-fluid intensive cases. Floor Mat Roll Carts are also available for effortless transportation from room to room.

Features and Benefits

Safe — Keeps your operating room floors clean, lowering the risk of slipping.

Efficient — Can be configured to accommodate your needs and helps to reduce turnover time.

Optimal Coverage — Perforated across the middle and every 22” to best fit the shape and size of your operating room floor.

Mobile — Optional Floor Mat Roll Cart allows for easy transportation and convenient storage.

Ordering Information

Product Code	Description	Thickness	Color	Size	Units/Case
FM0501	Double-Sided	Standard	Gray	29" x 110'	1
FM0501P	Waterproof Backsheet	Standard	Gray	29" x 110'	1
FM0502	Double-Sided E Roll	Lightweight	Gray	29" x 110'	1
FM0502P	Waterproof Backsheet E Roll	Lightweight	Gray	29" x 110'	1
FM0601C	Floor Mat Roll Cart	—	Silver/Blue	40" x 20"	1

Visit www.typenexmedical.com or call 866.897.3639 to learn more.