




# Take a stand against floor fluid buildup.

The Fluid Fighter® Suction Mat by Typenex® Medical is thoughtfully engineered to capture fluid from the surface while suction is drawn through a system of channels within the mat. This low-profile floor fluid management solution is designed with high-quality materials and can be securely adhered to your operating room floor.

# Fluid Fighter® | Suction Mat



## Features and Benefits

**Dual Action** — Absorbent top layer works simultaneously with an internal suction system to better manage fluid buildup.

**Low-Profile** — Four adhesive patches grip to the floor, keeping the mat flat and bulk-free.

**Quality Design** — Size and features are optimized to provide the best coverage and suction.

**Efficient** — Easy cleanup and disposal allows for a quick turnover between procedures.

## Ordering Information

Part Number	Description	Size	Units/Case
FM0202QD	Fluid Fighter®   Suction Mat	47" x 35"	12

Visit [www.typenex.com](http://www.typenex.com) or call 866-897-3639 to learn more.

Document Approvals  
Approved Date: 8/23/2022

Approval Task Verdict: Approve	Dylan Ash, (dash@typenex.com) Functional Area Approval 23-Aug-2022 19:10:17 GMT+0000
-----------------------------------	---

QA Approval Task Verdict: Approve	Michael Choquette, (mchoquette@typenex.com) Quality Assurance Approval 23-Aug-2022 19:24:01 GMT+0000
--------------------------------------	---