


Universal Tube | 12 mL


Instructions

1. Fill the tube with urine to the desired graduation.
2. Place the cap onto the tube if needed.
3. Perform urinalysis per your facility's protocol.



NOTE:

The Universal Tube is compatible with the KOVA® Urinalysis System.*

 Do not reuse

 In vitro diagnostic medical device

*KOVA® is a registered trademark of KOVA International.

 US0101QD
 Typenex Medical LLC
303 E Wacker Drive, Suite 1030
Chicago, Illinois 60601
www.typenexmedical.com
866.879.3639

TYPENEX 
MEDICAL

AN-1506-004.1