



# Get a grip on your segment piercing needs

The Typenex® Medical Typesafe® 2.0 | Segment Piercing Device is made to make pre-transfusion procedures easier. With an ergonomic grip designed to make it easy to handle when wearing gloves, it eliminates the need for scissors when sampling a segment.

# Typesafe® 2.0 | Segment Piercing Device



We've listened to customer feedback and worked with lab technicians to evolve the Typesafe 2.0 into the best possible segment piercing device for your lab. Learn more below about how our diligent innovations can help your lab technicians.

## Features and Benefits

**Efficient** - The dividers inside the device help guide the blood segment to the needle, provide support, and allow for easy dispensing.

**Safe** - Eliminates the need for scissors, reducing the risk of cross contamination and mess in the lab.

**Stable** - Fits securely on the test tube, allowing for increased steadiness when inserting the segment.

**Convenient** - Ergonomic flared design helps lab technicians grip the device, even when wearing gloves.

## Ordering Information

Product Code	Description	Color	Unit/Case
SP0201QM	Typesafe® 2.0   Segment Piercing Device	Green	1000
SP0201QD	Typesafe® 2.0   Segment Piercing Device	Green	500

For more information, visit [www.typenexmedical.com](http://www.typenexmedical.com).