


FinalCheck™ | Recipient Verification System



Designed to help reduce blood transfusion errors associated with patient identification, FinalCheck™ is an easy-to-use mechanical barrier system that can only be opened using a code that is placed on the patient's admissions band.

Patient Safety

With The FinalCheck Recipient Verification System you can assure that pre-transfusion patient ID checks are done at the bedside by putting this extra layer of safety in place.

Efficient

This physical barrier to the blood component can only be opened by entering the code that is applied to the patient's admissions band, ensuring the right blood goes to the right patient.

Ease of Use

The FinalCheck lock is easy to set and open using the 4-digit code on the patient label set.

Product Code	Color	Units/Case	Detail
LOCK			
C3000	N/A	100	
LABELS			
C3100	Red	3000/ream	
C3101	Green		
C3102	Yellow		
C3103	Blue		
C3113*	Blue	150	2 specimen tube stickers
BAGS			
C3200	Clear	100	
C3201	Clear		Biohazard Symbol
BANDS			
C4678	Red	150	Additional stickers

For more information, visit us online at www.typenexmedical.com or call us at 866.897.3639.