



Sample a blood segment with ease.

Typesafe® | Segment Piercing Device allows lab technicians to sample a blood segment without using scissors.

Typesafe® | Segment Piercing Device



Quick and easy to use for cross matching procedures, Typesafe® allows for segment sampling without the use of scissors.

Efficient - Allows the lab technician to dispense blood from a blood segment into a collection tube without the use of scissors.

Convenient - Easily discard the blood segment and the Typesafe device in designated hazardous medical waste containers.

Features and Benefits

- Fits securely onto 10 and 12 mm tubes
- Easy to handle with gloved hands
- Fits most standard segment sizes
- Not made with natural rubber latex



Product Code	Description	Color	Unit/Case
TS500	Typesafe® Segment Piercing Device	Green	500
4R4750	Typesafe® Segment Piercing Device	Green	1000

For more information, visit us online at www.typenexmedical.com or call us at 866.897.3639.