






Rad-Control Irradiates the Competition

As a trusted name in the blood bank, Typenex® Medical continues to develop innovative solutions with Rad-Control Irradiation Tags. This customizable product portfolio serves to verify the irradiation of blood products while enhancing efficiency. Our tags are stored at room temperature, conserving storage space and reducing trips to and from the refrigerator.

With unique styles and quality designs, Rad-Control Irradiation Tags can help to streamline protocols, save time, and lower the risk of transcription errors.

	Rad-Control Standard	Rad-Control Scale	Rad-Control Barcode	Rad-Sure™	RadTag®
					
Temperature Storage	Room Temperature	Room Temperature	Room Temperature	Room Temperature	-20°C to +6 °C
Length of Time at Room Temperature					
Up to 1 Day	✓	✓	✓	✓	✓
Up to 1 Month	✓	✓	✓	✓	✗
Up to Expiration Date	✓	✓	✓	✓	✗
Customization Available	✓	✓	✓	✗	✓
Log Book Sticker	✓	✓	✗	✗	✗
Barcode Verification	✗	✗	✓	✗	✗