



Budget-friendly tumescent delivery.

Typenex® Medical's tumescent tubing is an affordable drop-in solution suitable for all peristaltic pumps and tumescent delivery systems. Our single spike infiltration tubing is competitively priced and designed to complement existing workflows, allowing you to maximize cost savings with ease.

Tumescent Tubing



Features and Benefits

Single Spike — Tumescent tubing featuring two side clamps and silicone pressure insert tube.

Cost-Effective — Affordable drop-in solution can be seamlessly integrated into existing tumescence delivery systems.

Compatible — Suitable for use with all peristaltic pumps.

Ordering Information

Part Number	Description	Length	Units/Case
AMSDST0001	Tumescent Tubing	16 ft.	25

Tumescent Tubing is manufactured by Aligned Medical Solutions and distributed by Typenex® Medical.

Visit www.typenex.com or call 866-897-3639 to learn more.