

INTRODUCER KITS



Stay on the right track.

Introducer Kits by Typenex® Medical are designed to assist with smooth vascular access. Featuring a needle, guidewire, sheath, and dilator, this tailored all-in-one system offers ease and efficiency. With a variety of sizes and styles available, you can choose the solution that best suits each procedure.



FEATURES AND BENEFITS

Efficient Design — Innovative all-in-one system consolidates inventory and streamlines workflows.

Multiple Sizes — Various options to accommodate different procedures, preferences, and patients.

Tailored Components — Echogenic needles and marked specialty sheaths available for added ease and convenience.

Maximum Flow — Hemostasis valve adapter options help minimize blood loss during vasculature introduction.

STANDARD INTRODUCER KITS

Product Code	Description	Needle	Guidewire	HV Valve	Units/Case
MS0203	4 Fr x 12 cm Elite HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	10
MS0303	4 Fr x 24 cm Elite HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	5
MS0403	4 Fr x 45 cm Elite HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	10
MS0511	4 Fr x 70 cm Elite HV Micro-Introducer Kit	Standard	Nitinol	Yes	10
MS0603	5 Fr x 10 cm Coax Micro-Introducer Kit	Standard	Stainless Steel	No	10

MARKED INTRODUCER KITS

Product Code	Description	Needle	Guidewire	HV Valve	Units/Case
MS0244	4 Fr x 12 cm HV Micro-Introducer Kit, Marked Sheath	Echogenic	Nitinol	Yes	10
MS0344	4 Fr x 24 cm HV Micro-Introducer Kit, Marked Sheath	Echogenic	Nitinol	Yes	5
MS0444	4 Fr x 45 cm HV Micro-Introducer Kit, Marked Sheath	Echogenic	Nitinol	Yes	10

ECHOGENIC INTRODUCER KITS

Product Code	Description	Needle	Guidewire	HV Valve	Units/Case
MS0236	4 Fr x 12 cm Elite HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS0336	4 Fr x 24 cm Elite HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	5
MS0436	4 Fr x 45 cm Elite HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS0536	4 Fr x 70 cm Elite HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS0636	5 Fr x 10 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10

STIFFEN INTRODUCER KITS


Product Code	Description	Needle	Guidewire	HV Valve	Units/Case
MS1235S	4 Fr x 30 cm Coax Stiffen Micro-Introducer Kit	Standard	Nitinol	No	10
MS1335S	4 Fr x 40 cm Coax Stiffen Micro-Introducer Kit	Standard	Nitinol	No	5
MS1835S	4 Fr x 60 cm Coax Stiffen Micro-Introducer Kit	Standard	Nitinol	No	10

RADIOFREQUENCY KITS

Product Code	Description	Needle	Guidewire	HV Valve	Units/Case
MS1512	6 Fr x 12 cm HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS1612	6.5 Fr x 12 cm HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	5
MS1712	7 Fr x 12 cm HV Micro-Introducer Kit	Echogenic	Nitinol	Yes	10
MS1912	6 Fr x 10 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10
MS2012	6 Fr x 7 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10
MS2112	7 Fr x 10 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10
MS2212	7 Fr x 7 cm Coax Micro-Introducer Kit	Echogenic	Nitinol	No	10
MS2203	7 Fr x 7 cm HV Micro-Introducer Kit	Standard	Stainless Steel	Yes	5

Manufactured for Typenex® Medical.

Visit www.typenex.com or call 866-897-3639 to learn more.

 303 East Wacker Drive, Suite 1030
Chicago, Illinois 60601 USA

 866-897-3639
866-897-3630 (fax)

 sales@typenex.com

 www.typenex.com