

Fluid Fighter[®] Suction Canisters and Liners

Clean. Compact. Compatible.



INTRODUCTION

Typenex® Medical is proud to offer a comprehensive portfolio of reusable suction canisters, disposable flexible liners, and auxiliary accessories for streamlined fluid waste retention and disposal in clinical settings. Designed to support standardization and minimize medical waste, this sturdy yet lightweight system is available in a variety of configurations to suit your specific fluid management needs.

FEATURES AND BENEFITS



Fluid Fighter Suction Canisters

Shatter-Resistant — Reusable canisters feature an integrated one-piece design; backed by a 20-year guarantee.

Compatible — Adapters and custom solutions offer seamless integration and system compatibility with Cardinal Health, Medline, and Bemis Health Care.¹

Increased Visibility — Measurement graduations are clearly marked in black text.



Fluid Fighter Suction Liners

Compact — Optimal packaging and minimal assembly required, improving turnover times and conserving valuable storage space.

Self-Seating — Suction liners seat themselves in about 8 seconds, 50% faster than competing brands.

Closed System — Single-use liners are easy to remove after use, preventing contact with fluid waste.

Durable — Liners have been drop tested² and stretch up to 270% of initial dimension.

Infection Prevention — Integrated hydrophobic filter to prevent backflow of fluids and bacterial/viral transfer to vacuum source.³

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE	QUANTITY
SC0101	Fluid Fighter Suction Liner	1000 cc	100/Case
SC0102	Fluid Fighter Suction Liner	1500 cc	100/Case
SC0103	Fluid Fighter Suction Liner	2000 cc	100/Case
SC0104	Fluid Fighter Suction Liner	3000 cc	70/Case
SC0201	Fluid Fighter Suction Canister	1000 cc	15/Case
SC0202	Fluid Fighter Suction Canister	1500 cc	15/Case
SC0203	Fluid Fighter Suction Canister	2000 cc	12/Case
SC0204	Fluid Fighter Suction Canister	3000 cc	9/Case
SC0405CVL	Fluid Fighter Suction Liner Vacuum Control, L-Valve	—	50/Bag
SC0405CVT	Fluid Fighter Suction Liner Vacuum Control, T-Valve	—	50/Bag
SC0405XVTC	Fluid Fighter Suction Liner Vacuum, T-Connector	—	50/Bag
SC0405ET	Fluid Fighter Suction Liner Extension Tube (Female Connector)	6 mm x 5.9 ft	1/Each
SC0405LT	Fluid Fighter Suction System Trolley	—	1/Case
SC0405RH	Fluid Fighter Suction System Rail Hanger	—	1/Each
SC0405SCT	Fluid Fighter Suction System Serial Connection Tube	—	100/Case
SC0405ST	Fluid Fighter Suction System Silicone Tube	100 m (6 mm I.D.)	1/Case
SC0405WH	Fluid Fighter Suction System Wall Hanger Set	—	1/Each
SC0405WHA	Fluid Fighter Suction System Wall Hanger Adapter	—	1/ Each
PER05TD1618	Solidifying Agent, Dissolvable Pouch	—	150/Bag

Solidifiers are compatible with 1000 cc, 1500 cc, 2000 cc, and 3000 cc Fluid Fighter® Suction Liners.

¹Compatibility verified on Cardinal Health roll stands and part number 65652-145.

Compatibility verified on Medline part numbers DYNDSCCLW, DYNDSCC, and ORBASE.

Compatibility verified on Bemis Health Care part number 530510.

Cardinal Health™ is a registered trademark of Cardinal Health Technologies, LLC.

Medline™ is a registered trademark of Medline Industries, LP.

Bemis Health Care™ is a registered trademark of Bemis Manufacturing Company.

²Tylenex® Medical's Fluid Fighter® Suction Liners are drop tested according to ASTM D5276-98.

³Bacterial filtration efficiency of >99.85% and viral filtration efficiency of >97.5% validated by Pergo Medikal.

FOR MORE INFORMATION

Visit www.typenex.com or call
866-897-3639 to learn more.



303 East Wacker Drive, Suite 1030
Chicago, Illinois 60601 USA