




Quality you can count on.

Typenex® Medical offers a quality selection of microscope slides designed to support routine applications in liquid-based cytology, such as TCT and LCT slide preparation. These microscope slides feature hydrophilic surfaces, ground edges, 45° clipped corners, and various configurations to suit your specific liquid cytology needs.



ViewSafe™ | Liquid-Based Cytology Microscope Slides

ViewSafe™ | Liquid-Based Cytology Microscope Slides



Features and Benefits

Various Configurations — Printed brackets, circles, coated, and non-coated styles available to support a wide range of procedures.

Marking Areas — White frosted coating applied to one end.

Standard Dimensions — 25 x 75 mm, 1–1.2 mm thickness.



Ordering Information

Style*	Part Number	Description	Circles	Edges	Corners	Charged	Units/Case
Ⓐ	LC0101	Liquid-Based Cytology Microscope Slides	1	Ground	Clipped	No	1440
Ⓐ	LC0102	Liquid-Based Cytology Microscope Slides	1	Ground	Clipped	Yes	1440
Ⓑ	LC0201	Liquid-Based Cytology Microscope Slides	2, Stacked	Ground	Clipped	No	1440
Ⓑ	LC0202	Liquid-Based Cytology Microscope Slides	2, Stacked	Ground	Clipped	Yes	1440
Ⓒ	LC0203	Liquid-Based Cytology Microscope Slides	2, Offset	Ground	Clipped	No	1440
Ⓒ	LC0204	Liquid-Based Cytology Microscope Slides	2, Offset	Ground	Clipped	Yes	1440
Ⓓ	LC0301	Liquid-Based Cytology Microscope Slides	2 Brackets	Ground	Clipped	Yes	1440

*See above photo for letter denotation
Microscope slide mailers and storage sold separately.

Visit www.typenex.com or call 866-897-3639 to learn more.