



Conveniently packaged with case needs in mind.

Typenex® Medical's Anesthesia Breathing Circuit Kit is designed to offer convenience and ease of use for clinicians and staff. Intended for use with anesthesia and respiratory machines, this budget-friendly kit includes the components needed to deliver medical gases to patients during procedures. Our Anesthesia Breathing Circuit Kit can help consolidate resources, simplify purchasing processes, and conserve valuable storage space.

Anesthesia Breathing Circuit Kit



Features and Benefits

Convenience — Includes a bacterial/viral filter, 90-in expandable dual-limb breathing circuit, 3L reservoir bag, anesthesia mask (size 5), and 10-ft male-male gas sampling line.

Compliance — Interface meets ISO-5356-1 and ISO 80369-7 standards; compatible with most other brands.

Cost Savings — Competitively priced to help secure substantial savings.

Ordering Information

Part Number	Description	Units/Case
KB0101A	Anesthesia Breathing Circuit Kit, Adult	20

Manufactured for Typenex® Medical.

Visit www.typenex.com or call 866-897-3639 to learn more.